



AliMed ONE™ Wireless Motion Sensor

Instructions for Use, Setup, and Care

Intended Use

The AliMed ONE Wireless Motion Sensor pairs with the Alarm Base to provide caregiver notification of patient falls or elopement as part of an overall fall management program. The Motion Sensor can be mounted to a wall for wandering protection anywhere in the room.

Key Features

- 1. Universal Status Indicator**
Illuminates at varying times:
when Motion Sensor is turned on,
pairing with Alarm Base is complete,
or infrared beam is broken
- 2. Infrared Motion Sensor**
Infrared beam detects motion
- 3. On/Off Switch** Turns Motion Sensor on or off
- 4. Pair Status Indicator**
Illuminates when Alarm Base is in range for pairing
- 5. Battery Compartment Button**
Press to open battery compartment
- 6. Universal Wall Mount and Bracket** Holds Motion Sensor and mounts to wall




Figure 1



Figure 2




Figure 3

Setup and Operating Instructions

The Motion Sensor is installed by Nurse or Biomedical staff.

1. Motion Sensor requires one (1) 9V battery to operate (sold separately, see reverse).
2. Press the Battery Compartment Button (Fig. 3) to open the battery compartment.
3. Align the male and female connectors on the top of the battery to the connectors attached to the unit. Place in battery compartment.
4. The Motion Sensor is now ready for pairing.
5. Close the battery compartment door until it snaps into place.

Pairing the Motion Sensor with the Alarm Base:

1. Ensure the On/Off Switch is set to the "OFF" position (Fig. 2).
2. Press Pair Button on the Alarm Base.
3. The Base will beep and the first open channel will blink green and red indicating channel is ready for pairing.
4. Place Motion Sensor next to the Base matching PAIR on the Motion Sensor to PAIR on the Base (Fig. 4). *Note: Pair Status Indicator on Motion Sensor (Fig. 3) will illuminate when in range of the Base.*
5. If pairing is successful, the Base will beep and the assigned channel will blink green once. The Universal Status Indicator (Fig. 1) on Motion Sensor will blink three (3) times.
6. The Motion Sensor is now paired with the Base. The assigned channel on the Base will continuously blink at approximately 5-second intervals indicating Motion Sensor is in use.

PAIR indicator located on right side of Motion Sensor



Match PAIR on Motion Sensor to PAIR on Base

Figure 4

Wall Mounting Instructions

1. Secure the bracket to the wall with the screws provided (Fig. 3).
2. Align mounting slot on Motion Sensor with mounted bracket.
3. Slide Motion Sensor from the left until secured in place on the mounted bracket (Fig. 5).
3. Adjust Motion Sensor angle as needed.

Note: Once mounted to wall, Motion Sensor can be easily removed following the same installation instructions.

4. The Motion Sensor is now ready to use.

Slide Motion Sensor from the left until secured in place on bracket




Figure 5

Care

1. Wipe with soft, non-abrasive cloth using a 1:10 bleach solution (contains 0.5% chlorine concentration) or follow your facility's cleaning protocol.
2. Do not immerse in fluids.
3. After cleaning, check unit is functioning properly.

Warning

Failure to comply with warning may result in injury or death. This device is not suitable for all individuals. Other devices may be required. This device is not a substitute for visual monitoring by caregiver. The manufacturer does not claim that this device will stop elopement and/or stop falls. This device is designed to augment caregiver's comprehensive resident mobility management program. This product is to be used by a licensed caregiver only. As such, it is the entire responsibility of the caregiver to ensure that the system is properly installed and tested. Further, the system is not designed to replace good caregiving practices including, but not limited to:

- a.) Direct patient supervision
- b.) Adequate training for personnel for fall prevention and elopement
- c.) Testing of the system before each use

Limited Warranty: Seller warrants that the goods sold will be free from defects of workmanship and manufacturing for a period of one year from the date of sale.

Reorder Information

- #713817 AliMed ONE Wireless Motion Sensor
 #71088 9V Battery, 12/pk

Additional AliMed ONE System Components

- #713815 Wireless Alarm Base
 #713816 Wireless Hall Alert
 #713818 Wireless Sensor Pad Transmitter (one required for each sensor pad)
 #713819 30-Day Wireless Bed Sensor Pad, ea #713866, 10/pk
 #713834 6-Month Wireless Bed Sensor Pad, ea #713867, 10/pk
 #713838 1-Year Wireless Bed Sensor Pad, ea #713868, 10/pk
 #713821 30-Day Wireless Chair Sensor Pad, ea #713822, 10/pk
 #713836 6-Month Wireless Chair Sensor Pad, ea #713837, 10/pk
 #713840 1-Year Wireless Chair Sensor Pad, ea #713841, 10/pk
 #713823 14-Day Wireless Toilet Sensor Pad, ea

| Environmental Conditions | | | |
|--------------------------|----------------|----------------|----------------|
| | Operational | Storage | Transport |
| Temperature | 10°-40°C | -25°-70°C | -25°-70°C |
| | 50°-104°F | -13°-158°F | -13°-158°F |
| Humidity | 30%-75% | 0%-95% | 0%-95% |
| | Non-condensing | Non-condensing | Non-condensing |
| Pressure | 700-1060 hPa | 700-1060 hPa | 700-1060 hPa |

CONFORMS TO AAMI STD ES60601-1, IEC STD 60601-1-6
 CERTIFIED TO CSA STD C22.2 NO. 60601-1, NO. 60601-1-6

Symbol Glossary

| | |
|--|---|
| | Manufacturer |
| | Batch Code/Lot Number |
| | Catalog Number/Item Number |
| | Serial Number |
| | European Conformity Mark |
| | Authorized representative in the European Community |
| | Tested/certification for compliance to published industry standards |
| | Consult instructions for use |
| | Caution: Read all warnings and precautions in instructions for use |
| | DO NOT THROW IN TRASH |
| | Indicates FCC has certified device |
| | Indicates compliance to RoHS and is free of substances banned by RoHS |



Manufactured for:

AliMed® Inc.
 297 High St., Dedham, MA 02026 USA
 PH 800.225.2610 • FAX 800.437.2966



MDSS GmbH
 Schiffgraben 41
 30175 Hannover,
 Germany

AliMed.com
 #7721REV0919 CCO982
 ©2019 AliMed® Inc. All rights reserved.
 Made in China