

Early warning without restricting motion

Alarm sounds before patient finishes unfastening belt



■ **Sounds before the patient gets up**

■ **Fewer false alarms**

■ **Patient can move freely in bed**

Roll-Control Belt

The most restless patients at high fall risk can also have the most false alarms. The Roll-Control Belt wraps around the patient and allows the patient to roll from side-to-side, even sit up without leaving the bed or triggering a false alarm. But as soon as they start to open the hook-and-loop front, the monitor sounds to alert the caregiver. After the monitor is triggered, the patient still has to open more hook-and-loop, giving caregivers more time to respond.


The Roll-Control Belt is not considered a restraint if the patient is aware they can open it in an emergency.

Setup is simple, the belt is comfortable, and when used with the AliMed TR2 or Worry-Free™ monitor, it is both early warning and tamper-resistant—the ideal combination for your most restless patients.

To order, see reverse side.



12" of hook-and-loop belt gives caregiver time to respond in a timely manner



When engaged, patient is able to roll freely from side-to-side, or even sit up without sounding alarm



Roll-Control Belt

Choose your system configuration:



#70088 Roll-Control Belt w/TR2 Alarm

Replacement Items

#70089 Roll-Control Belt



#71 3316 Roll-Control Belt w/Worry-Free™ Alarm

Replacement Items

#70089 Roll-Control Belt

Visit AliMed.com to learn more about the TR2 and Worry-Free™ Alarm.

Partners in the Continuum of Care

AliMed is a leading developer, manufacturer, and distributor of high-quality patient care, rehabilitation, orthopedic, operating room, medical imaging, and ergonomic products. For over 40 years, AliMed has served the hospital, long-term care, and rehabilitation markets with a broad selection of specialty products.

Through established partnerships with more than 500 suppliers, a reputation for exceptional customer service, and a comprehensive range of unique products, AliMed provides more than 70,000 different products to healthcare facilities and businesses worldwide.

AliMed provides the following:

- Service to all medical markets with an extensive array of products, including those that are difficult to find
- Competitive pricing through contracts
- Field sales and dedicated customer service support
- Same-day shipping for in-stock items*
- Reduced acquisition costs through EDI ordering
- E-commerce and specialized catalogs to easily find and order products

*Any order placed before 3 p.m. ET, if in stock, ships the same day!

Pub. #6363REV0215 Product #11145

©2015 AliMed® inc. All rights reserved. • Printed in U.S.A.



AliMed.com

Call: 800.225.2610

Fax: 800.437.2966