

Optimal positioning and user comfort



Solid Seat Insert™ with Pommel

Best for patients who:

- ✓ Have unique positioning needs
- ✓ Are prone to forward sliding or leaning
- ✓ Require post-operative positioning support or pressure relief

 **Choosing the right cushion — read our blog:**
[AliMed.com/choosing-a-cushion](https://www.ali-med.com/choosing-a-cushion)

■ Offer a variety of specialty positioning and support

■ Limit forward sliding and reduce pressure injuries

■ Promote upright posture and improved spinal alignment

AliMed Specialty Cushions effectively manage specific needs, including limiting forward sliding and promoting postural alignment, or giving extra pressure relief and post-operative support.

Postural Support

Users with muscle weakness are at risk of forward sliding, which leads to poor alignment. Our **Sit-Straight™ Cushions** promote upright posture by keeping the hips positioned slightly lower than the knees while **Pommel Cushions** limit forward sliding and reduce the risk of adduction. For additional support, **Anti-Thrust Cushions** maintain user's neutral pelvic position with a higher front and lower rear surface.

Pressure Relief

Contour Cushions offer advanced coccyx support with a sculpted surface for pressure distribution and deep cutout for sacral pressure relief.

Post-Op Support

AliHip™ Cushions feature a sloping cutout that accommodates left or right hip.

*See reverse for
ordering information* 

Cushions for unique needs



Convex Base — T-Gel™
version

Anti-Thrust Cushions

Maintain neutral pelvic position

Removable cover: water-resistant, polyester AliTex™ top with a vinyl nonskid base; Black • Rear ties • **Convex Base:** 5"H front, 3"H back • **Standard Flat Base:** 5½"H front; 3"H back • **Flat Base w/Support Board:** 7½"H front; 3"H back

CONVEX BASE

- #1804 Foam, 18"W x 16"D
- #1931 Foam, 16"W x 16"D
- #1975 T-Gel/Foam, 18"W x 16"D
- #1974 T-Gel/Foam, 16"W x 16"D

FLAT BASE, 18"W x 16"D

- #1532 Standard
- #1531 w/Solid Support Board




Solid Seat Insert™ with Pommel

Reduces risk of adduction and forward sliding

Removable cover (not shown): nonslip vinyl, incontinence-resistant, fully launderable, Black • Rear ties

DESCRIPTION	CUSHION HEIGHT		ITEM #	
	CENTER	EDGES	16"W x 16"D	18"W x 16"D
Standard Foam	1½"	¾"	#1311	#1312
w/T-Foam Layer	2½"	1¾"	#1313	#1314
w/T-Gel Layer	2½"	1"	#1315	#1316




Sit-Straight™ Basic Cushions

Wedge cushion effectively eliminates forward sliding and sacral sitting

¾"H tapers to 3"H rear • Firm Polyfoam • Removable cover (not shown): wipe-clean, water-resistant, nonskid bottom, fully launderable, Black • Rear ties

- #1250 16"W x 16"D
- #1251 18"W x 16"D



Contour Cushions

Cutout provides full sacral pressure relief
18"W x 16"D • Polyurethane foam • Meets CAL 117
• Solid Seat Insert™ base for sling-seat wheelchairs

- #1166 Flat Base
 - #1167 Solid Seat Insert Base
- Specify Firm (for users over 140-lbs.) or Soft.



AliHip™ Cushions, Home Model

Keep post-op hip properly positioned
18"W x 16"D • 4"H slopes to 1"H in front • Dense Polyfoam • Optional cover: rear ties for secure wheelchair attachment, Black polyester

- #11065 Right Hip
- #11066 Left Hip
- #11067 Optional Cover, Right
- #11068 Optional Cover, Left

Visit AliMed.com/cushions for our full line.

Any order placed before 3 p.m. EST, if in stock, ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #7671 Product #11509

©2019 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com / customerservice@AliMed.com / 800.225.2610