

Designed for versatility




Table Socket Family

Simple, durable clamps securely hold O.R. table accessories

From our **Original Clark Socket** to the **FREEDOM® 360 Socket** with multiple adjustment options, there's an AliMed Socket for every need. Each is designed for ease of use, durability, and versatility, and they're all backed by a 3-year limited warranty. Choosing the right Clark Socket depends on the following:

- **Accessory position:** whether rotating (angled) or fixed (perpendicular or parallel to the table) during procedures
- **Rotation range:** AliMed Clark Sockets adjust quickly and precisely in 12°, 9°, or 5° increments
- **Accessory compatibility:** choose between flat blades (1"W x ¼" thick) or round posts (1½" diameter), or save time by choosing a clamp that accepts both
- **Table accommodation:** fits all O.R. tables with U.S. side rails (¾"W x 1½"H)
- **Quality materials:** choose durable stainless steel or anodized aluminum, which is 50% lighter
- **Guaranteed** with a 3-year limited warranty

SHOP AliMed.com



EasyLock™ Blade Clamps

:: Complete line of Table Socket Products

Original Clark Socket



#936673 Original Clark Socket

- **Durable:** Stainless steel construction
- Fits round posts
- Slides on between standoffs
- **Table accommodation:** Fits all O.R. tables with U.S. side rails (¾"W x 1½"H)
- 3-year limited warranty

2¾"W x 3⅞"H x 1¾"D • Stainless steel • 2 lbs.
 • Accepts round posts up to 1½" diameter
 • Rotates 360° in 12° increments

Clark Socket Plus: Slide-On




#937130 Clark Socket Plus Slide-On

- **Durable:** Stainless steel construction
- **Versatile:** Fits both round posts and flat blades
- Slides on between standoffs
- **Table accommodation:** Fits all O.R. tables with U.S. side rails (¾"W x 1½"H)
- 3-year limited warranty

2¾"W x 3⅞"H x 1¾"D • Stainless steel • 2 lbs.
 • Accepts round posts up to 1½" diameter
 • Accepts flat blades up to 1"W x ¼" thick
 • Rotates 360° in 9° increments

Clark Socket Plus: EZ-On



#937344 Clark Socket Plus EZ-On

- **Durable:** Stainless steel construction
- **Versatile:** Fits both round posts and flat blades
- **Quick setup and breakdown:** EZ-On attachment connects anywhere along side rail
- **Table accommodation:** Fits all O.R. tables with U.S. side rails (¾"W x 1½"H)
- 3-year limited warranty

2¾"W x 3⅞"H x 1¾"D • Stainless steel • 2 lbs.
 • Accepts round posts up to 1½" diameter
 • Accepts flat blades up to 1"W x ¼" thick
 • Rotates 360° in 9° increments

FREEDOM® Clark Socket



#934725 FREEDOM® Clark Socket

- **Lightweight:** Anodized aluminum construction
- **Versatile:** Fits both round posts and flat blades
- **Quick setup and breakdown:** EZ-On attachment connects anywhere along side rail
- **Table accommodation:** Fits all O.R. tables with U.S. side rails (¾"W x 1½"H)
- 3-year limited warranty

3½"H x 2½" diameter • Anodized aluminum • 0.8 lb.
 • Accepts round posts up to 1½" diameter
 • Accepts flat blades up to 1"W x ¼" thick
 • Rotates 360° in 5° increments


FREEDOM® 360 Socket




#937932 FREEDOM® 360 Socket

- **Lightweight:** Anodized aluminum construction
- **Versatile:** Three independent adjustments; fits both round posts and flat blades
- **Quick setup and breakdown:** EZ-On attachment connects anywhere along side rail
- **Table accommodation:** Fits all O.R. tables with U.S. side rails (¾"W x 1½"H)
- 3-year limited warranty

5¼"W x 2½"H x 3⅞"D • Anodized aluminum • 2 lbs.
 • Accepts round posts up to 1½" diameter
 • Accepts flat blades up to 1"W x ¼" thick
 • Rotates 360° in 5° increments



Independently adjust **placement, height, and angle** for precise positioning.



Compare		Product #	Accepts Round Posts	Accepts Flat Blades	Independent Adjustments	Side Rail Attachment*	Calibration Marks	Rotational Increments	Material
A.	Original Clark Socket	#936673	●	—	No	Slide-On	No	12°	Stainless Steel
B.	Clark Socket Plus Slide-On	#937130	●	—	No	Slide-On	No	9°	Stainless Steel
C.	Clark Socket Plus EZ-On	#937344	●	—	No	EZ-On	No	9°	Stainless Steel
D.	FREEDOM Clark Socket	#934725	●	—	No	EZ-On	Yes	5°	Aluminum
E.	FREEDOM 360 Socket	#937932	●	—	Yes	EZ-On	Yes	5°	Aluminum

*Accepts up to an 1½" diameter round post and/or accepts up to a 1"W x ¼" thick flat blade.
 *EZ-On models connect anywhere along the rail; Slide-On models slide on between standoffs.



For AliMed Fixed Sockets and Clamps, see back cover.

AliMed Fixed Sockets and Clamps



#934706 Flat Blade Socket, 2/cs

- **Lightweight, durable** anodized aluminum

1³/₄"H x 2" diameter • Anodized aluminum construction • 5.6 oz. • Accepts 1"W x 1/4" thick flat blades • Fits all O.R. tables with standard U.S. side rails (3/8"W x 1¹/₈"H) • 3-year limited warranty



#934707 Round Post Socket, 2/cs

- **Lightweight, durable** anodized aluminum

1³/₄"H x 2" diameter • Anodized aluminum construction • 5.6 oz. • Accepts 5/8" diameter round posts • Fits all O.R. tables with standard U.S. side rails (3/8"W x 1¹/₈"H) • 3-year limited warranty



#934726 Combo Socket, 2/cs

- **Lightweight, durable** anodized aluminum
- **Versatile:** Fits both round posts and flat blades

1³/₄"H x 2" diameter • Anodized aluminum construction • 5.6 oz. • Accepts 5/8" diameter round posts • Accepts 1"W x 1/4" thick flat blades • Fits all O.R. tables with standard U.S. side rails (3/8"W x 1¹/₈"H) • 3-year limited warranty



#931899 Universal Rail Clamp

- **Fits anywhere along the rail**

2"W x 2"L x 2¹/₂"H • Anodized aluminum construction • 8 oz. • Accepts 5/8" diameter round posts • Accepts 1"W x 1/4" thick flat blades • Fits all O.R. tables with standard U.S. side rails (3/8"W x 1¹/₈"H) • 3-year limited warranty



#930752 Side Rail Socket

- **Lightweight, durable** anodized aluminum

1⁵/₈"W x 2⁷/₈"L x 2"H • Anodized aluminum construction • 5.6 oz. • Accepts 1"W x 1/4" thick flat blades • Fits all O.R. tables with standard U.S. side rails (3/8"W x 1¹/₈"H) • 3-year limited warranty



#95-196 Side Rail Socket

- **Extremely durable** stainless steel

1³/₄"W x 1³/₄"L x 1¹/₈"H • Stainless steel construction • 12 oz. • Accepts 1"W x 1/4" thick flat blades • Fits all O.R. tables with standard U.S. side rails (3/8"W x 1¹/₈"H) • 3-year limited warranty

Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

The AliMed Advantage:

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items*
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

*Any order placed before 3 p.m. EST, if in stock, ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #7223REV0518 Product #11206

©2018 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com

Call: 800.225.2610

Fax: 800.437.2966