

ULTRApadded™ BEND TO FIT Contracture Braces

Static progressive stretch helps gain or maintain ROM

Foam-covered, bendable cuffs and spine for quick, tool-free adjustments

Padded terrycloth cover absorbs moisture to help prevent skin maceration



ELBOW

Elbow Orthoses

KNEE

Knee Orthoses



WRIST

Grip Orthoses

AliMed ULTRApadded™ Contracture Braces

provide protection and support for weak or flaccid extremities, and are effective in achieving gradual extension of non-fixed contractures. In cases of neuromuscular impairment, these braces aid in positioning and joint protection, even immobilizing painful joints.

The bend-to-fit spine can be custom fit without the use of tools or heat, and can accommodate joint deformities.

One size fits left or right • Removable, machine-washable cover • Green

Elbow Orthosis

#52676

- Provides gradual extension of non-fixed elbow contracture
- Helps manage elbow flexor spasticity and deformity

Knee Orthosis

#66888

- Provides gradual extension of non-fixed knee contracture
- Helps maintain gains from surgery or other procedures

Grip Orthosis

#51 0327

- Comforts and positions weak or deformed hands
- Adjustable side wings prevent ulnar and radial deviation



Pub. #7525 ©2017 AliMed® Inc. All rights reserved.