



297 High Street • Dedham, MA 02026  
Ph: 888-625-4633 • Fax: 888-352-8391 • Web: [www.AliMed.com](http://www.AliMed.com)


---

## ***Custom Sandbag Information & AliMed Quote Process***

- Once the custom quote request form has been filled out, please fax to 888-352-8391.**
- A formal price quote will be generated and sent via fax or email to the customer contact as soon as possible. *(Please note that response times may vary depending on the complexity of the request. If you do not hear a response within 5 business days, please call 888-625-4633.)***
- If the quote request form is not filled out completely, a quote will *NOT* be generated. *(The form will be returned for completion and will need to be re-sent by the requestor for quote processing.)***
- All quotes are non-negotiable, FINAL, and expire after 60 days. *(Please fill out a new form for all requests after the expiration date.)***
- To place an order, please fax back a copy of the quote along with a purchase order to 888-352-8391. A fax or email confirmation will be sent once order has been processed. *(Please reference the quote number, item number, and product description on the PO when ordering.)***
- Please NOTE: The attached forms are for a *QUOTE REQUEST ONLY* and are not considered official purchase orders.**
- Custom orders are non-returnable. *Exceptions will be made to products made incorrectly due to an AliMed error.***

**AliMed** CUSTOM SANDBAG INFORMATION

**STYLE OPTIONS:**

<p><b><i>Envelope Style</i></b></p> 	<ul style="list-style-type: none"> <li>• Narrows at the ends and thicker in the middle</li> <li>• Specify length and width when ordering</li> </ul>
<p><b><i>Envelope Style w/grommets</i></b></p> 	<ul style="list-style-type: none"> <li>• Hanging grommets available for envelope style ONLY sandbags</li> <li>• Specify where grommet should be on sandbag (along width or length?)</li> </ul>
<p><b><i>Box Style</i></b></p> 	<ul style="list-style-type: none"> <li>• Has an even thickness throughout with squared edges</li> <li>• Specify length, width and thickness when ordering</li> </ul>

**COVER OPTIONS:**

<p><b><i>Black Conductive Vinyl</i></b></p>	<ul style="list-style-type: none"> <li>• Conductive vinyl, slip resistant, durable, smooth</li> <li>• Sewn seams are not radiolucent but are located on pad edges</li> </ul>
<p><b><i>Naugahyde</i></b></p>	<ul style="list-style-type: none"> <li>• Thicker alternative to black conductive vinyl</li> <li>• Slip resistant, textured</li> </ul>

**FILL MATERIAL:**

<p><b><i>Sand</i></b></p>	<ul style="list-style-type: none"> <li>• Silica sand, comparable to beach sand</li> <li>• MRI Safe</li> </ul>
<p><b><i>AliSand Lightweight Pellets</i></b></p>	<ul style="list-style-type: none"> <li>• Made of plastic beads that mold just like sand but weigh significantly less</li> <li>• Good for to use with a larger sandbag to obtain minimum weight</li> <li>• MRI Safe</li> </ul>
<p><b><i>Lead Shot</i></b></p>	<ul style="list-style-type: none"> <li>• Lead beads</li> <li>• Good for a smaller sandbag to obtain maximum weight.</li> </ul>



### Custom Sandbag Quote Request Form

*Please print this document, fill out the form and fax to 888-352-8391 to have a quote sent for the custom product requested.*

*Use the sandbag informational sheet as a guide while filling out this sheet and to learn the quote request & order process.*

**(ALL FIELDS MUST BE FILLED IN TO PROCESS QUOTE REQUEST)**

**Customer Contact Information:**

Contact Name (First & Last): \_\_\_\_\_

Customer #: \_\_\_\_\_ Company/Facility: \_\_\_\_\_


Address: \_\_\_\_\_ City/State: \_\_\_\_\_

Zip Code: \_\_\_\_\_ Best time to Contact: \_\_\_\_\_ (AM/PM)

Ph: \_\_\_\_\_ Fax: \_\_\_\_\_ Email: \_\_\_\_\_

*Sandbags are made with double wall construction for added strength, durability and to prevent leaking.*

**ENVELOPE STYLE SPECIFICATIONS: (measurements in inches)**



Sandbag Width A = \_\_\_\_\_ Sandbag Length B = \_\_\_\_\_

**Hanging Grommets (available with envelope style sandbags ONLY)**


Check One:

Yes

No

\* Hanging Grommets are available for envelope style **ONLY** sandbags. If you would like to place hanging grommets on your sandbag, please draw placement on the diagram.

**BOX STYLE SANDBAG SPECIFICATIONS: (measurements in inches)**



Sandbag Length A = \_\_\_\_\_ Sandbag Width B = \_\_\_\_\_

Sandbag Thickness C = \_\_\_\_\_

**FILL THIS SECTION OUT FOR BOTH BOX STYLE AND ENVELOPE STYLE CUSTOM SANDBAGS:**

**COVER OPTION: Check off one:** *(reference the sandbag informational sheet)*

Black Conductive Vinyl

Naugahyde

**DENSITY OPTION: Check off one**

Full Pack

Medium Pack

Loose Pack

**FILL OPTIONS: Check off one** *(reference the sandbag informational sheet)*

Sand

AliSand

Lead Shot

**Which component is more important?**

Weight

Size

Finished Weight Desired (in lbs): \_\_\_\_\_

**Desired requirements of product? (radiolucent, bacteriostatic, waterproof, etc):**

\_\_\_\_\_

\_\_\_\_\_

**Application of Product: (length of procedure, type of patient, type of procedure):**

\_\_\_\_\_

\_\_\_\_\_

**Questions/Comments:**

\_\_\_\_\_

\_\_\_\_\_