

AliMed®

ProStim™ Electrotherapy

High on features, priced as low as...

\$64^{ea}



Digital EMS 3500
#TY32411

\$42^{ea}



Digital TENS 3000
#TY32394

\$28^{ea}




Analog TENS 2500
#TY32396

**FREE CASE
OF ELECTRODES**

with purchase of
ProStim™ Digital TENS!

#TY32394A.....\$49.95 set*

*Set includes one #32394 TENS 3000
and one case of #31700 2" square
Prostim Electrodes (4/sh, 10 sh/cs).
Offer expires 7/31/09



Introducing ProStim™ Electrotherapy by AliMed®

Finally, high-quality electrotherapy units that allow you to outfit your facility without compromising features. All ProStim units feature dual-channels and come complete with timer, convenient storage case, four reusable electrodes, two lead wires, batteries, and instruction manual. **Free with every unit!**



Shown actual size

Digital Muscle Stimulator 3500

- Six preset protocols, plus one customizable protocol
- Two waveforms
- Patient-lock mode

As low as
\$64⁰⁰ea

AliMed® ProStim® Digital EMS 3500 is packed with features typically found in higher-priced models. Now you and your patients can have a full-featured muscle stimulator at an affordable price. Six pre-programmed modes allow Quick Start set-up and one custom mode for individual protocols. High-strength plastic case with slide cover protects controls. Easy-to-read LCD display. Lock mode prevents accidental change of treatment parameters. Includes carry case, four ProStim Electrodes, two AA alkaline batteries, belt clip and two lead wires.

EMS 3500 (#TY32411)

Channel	Dual, isolated
Wave Form	Biphasic Symmetrical (galvanic-free), Symmetrical alternated monophasic
Pulse Rate	10 Hz-60 Hz
Pulse Width	50 µs-400 µs in steps of 50 µs
Timer	Continuous, 10, 20, 30, 45, 60 or 90 minutes
Intensity	5 to 100mA
Frequency	1 to 110 Hz
Ramp Time	0-5 sec. in steps of 0.25 sec
On/Off Time	1-60 sec. in steps of 1 sec.
Power Source	(2) AA alkaline batteries
Dimensions	2½"W x 5¼"L x ¾"D
Weight	0.25 lbs

#TY32411 EMS 3500\$ 79.95 ea
#TY32 41106 6/cs383.76 cs (**\$64.00^{ea!}**)

Caution: Federal law restricts devices on this page to sale by or on the order of a licensed professional.

Analog TENS 2500



Timer

- Three modes
- Single waveform
- Plastic shield protects dials from accidental changes in intensity

As low as
\$28⁰⁰ ea

AliMed® ProStim® 2500 TENS w/Timer for easy-to-use pain management therapy. Controls are protected behind front panel cover. Separate intensity controls are protected from accidental changes by the clear, flip-top cover. High-visibility LED indicates power status. Includes carry case, four ProStim Electrodes, one 9V alkaline battery, belt clip and two lead wires.

TENS 2500 (#TY32396)

Channel	Dual, isolated
Modes	Continuous, Burst, Modulated
Wave Form	Asymmetrical Bi-phasic square
Amplitude	0-80mA peak into 500 load each channel, constant current
Pulse Rate	2 Hz-150 Hz
Pulse Width	30 µs-260 µs
Timer	Continuous, 30, or 60 minutes
Power Source	9V alkaline battery
Dimensions	2½"W x 3¾"L x 1" D
Weight	0.25 lbs

#TY32396 TENS 2500\$ 34.95 ea
 #TY32 39610 10/cs280.00 cs (**\$28.00 ea!**)

Shown actual size

Digital TENS 3000



Quick Start-up

- Five modes
- Two waveforms
- Patient-lock mode
- High-visibility digital display

As low as
\$42⁰⁰ ea

AliMed® ProStim® 3000 Digital TENS No complicated controls to learn, so your standard protocol can start with just one press of the button. Easy-to-read LCD with large-button control panel. Choose from two waveforms and five operation modes. Lock mode prevents accidental change of treatment parameters. Unit has low-battery indicator, countdown timer, and automatic shut off. Includes belt clip, carry case, four ProStim Electrodes, one 9V alkaline battery, and two lead wires.

TENS 3000 (#TY32394)

Channel	Dual, isolated
Modes	Continuous, Burst I, Burst II, Modulate I, Modulated II
Wave Form	Symmetrical Bi-Phasic Rectangular, Monophasic Rectangular
Amplitude	0-150mA peak into 500 Ω load
Pulse Rate	1 Hz-150 Hz
Pulse Width	50 µs-250 µs
Timer	Continuous, 15, 30, 45, 60, or 90 minutes
Power Source	9V alkaline battery
Dimensions	4"L x 2½"W x 1"D
Weight	0.25 lbs

#TY32394 TENS 3000\$ 49.95 ea
 #TY32 39406 6/cs252.00 cs (**\$42.00 ea!**)

Shown actual size

ProStim™ Electrodes

Long-Lasting Electrodes

- Even current distribution improves patient comfort
- Flexible cloth backing conforms to body contours
- Use for TENS, Muscle Stimulation, and Interferential
- 40/cs—4/sheet, 10 sheets/case

As low as
\$4³⁸_{sh}



#TY31705
2" x 4" Oval



#TY31704
1½" x 2½" Oval



#TY31706
2" x 3½" Rectangle



#TY31700
2" Square



#TY31703 2¾" Round



#TY31703 2" Round



#TY31703 1¼" Round

ProStim™ Electrodes		CS\$	CS\$ (4+)	CS\$ (10+)
2" Square	#TY31700	67.50 ^{cs} (6.25 ^{sh})	54.50 ^{cs} (5.45 ^{sh})	51.25 ^{cs} (5.13 ^{sh})
2" Round	#TY31701	67.50 ^{cs} (6.75 ^{sh})	54.50 ^{cs} (5.45 ^{sh})	51.25 ^{cs} (5.13 ^{sh})
2¾" Round	#TY31702	83.75 ^{cs} (8.38 ^{sh})	67.00 ^{cs} (6.70 ^{sh})	62.75 ^{cs} (6.28 ^{sh})
1¼" Round	#TY31703	58.75 ^{cs} (5.88 ^{sh})	46.75 ^{cs} (4.68 ^{sh})	43.75 ^{cs} (4.38 ^{sh})
1½" x 2½" Oval	#TY31704	74.25 ^{cs} (7.43 ^{sh})	60.00 ^{cs} (5.80 ^{sh})	56.50 ^{cs} (5.65 ^{sh})
2" x 4" Oval	#TY31705	91.25 ^{cs} (9.13 ^{sh})	72.65 ^{cs} (7.27 ^{sh})	68.00 ^{cs} (6.80 ^{sh})
2" x 3½" Rectangle	#TY31706	86.25 ^{cs} (8.63 ^{sh})	68.75 ^{cs} (6.88 ^{sh})	64.25 ^{cs} (6.43 ^{sh})

Light-Duty Electrodes

- Low-profile carbon conductive layer
- Cloth backing conforms to body contours
- Use for TENS and Muscle Stimulation
- 1.8" x 1.8" square
- 100/cs—4/sheet, 25 sheets/case

#TY32270 Light-Duty.....\$103.75 cs.....\$75.00 cs (4+) (\$3.00^{sh})



As low as
\$3⁰⁰_{sh}