



AliMed®

"Helping people help people"

Hand rest, reminder, and alarm in one

QUALCARE™ BREAKAWAY LAP CUSHION WITH ALARM

Provides patients in wheelchairs with a soft, comfortable surface to rest their hands. As a patient attempts to rise, a magnetic strap on one side of the lap cushion releases, activating the alarm. Alarm unit features on/off switch, alarm or chime tone with 85 dB or 100 dB sound levels, and a low-battery warning light. Includes one 9V battery.


#75188 BreakAway Lap Cushion

ACCESSORIES (SOLD SEPARATELY)

Replacement Lap Cushion #76808

Replacement Alarm Unit #77141

9V Batteries, 8/pk #71088



KEY BENEFITS

- ▶ Alarm sounds before patient rises
- ▶ Acts as a constant reminder for patient to ask for assistance
- ▶ Easy to clean
- ▶ Not a restraint for most patients*
- ▶ Provides a hand rest and soft work surface for patients

** If the resident is unable to remove a device without assistance, it is considered a restraint under OBRA guidelines.*



297 High Street, Dedham, MA 02026 USA
www.AliMed.com
E-mail: info@alimed.com

Order Toll Free 800.225.2610
24-Hour Fax 800.437.2966
General Info 781.329.2900