


**Be green...
pass it on**

ROUTE TO:

- Bariatric Coordinator
- Facilities Manager
- Equipment Planner
- _____

Is your facility bariatric- ready?



Retrofit your facility to meet standards for accreditation and caregiver safety:
OSHA • DPH • CARF • JCAHO

- Minimize caregiver back strain:
Use bariatric repositioning aids p10
- Reduce the risk of bed sores:
Enable patients to independently reposition themselves in bed p9
- Avoid a disaster:
Reinforce your wall-mounted toilets p4
- Stop using cushions that bottom out:
Use cushions you can rely on p12



SPECIAL!
1,200-lb. capacity
no-splash Bed Pan
p5

6617

How important is it to make your facility bariatric-ready?

Meeting the needs of obesity

As the obese population continues to grow, medical costs increase accordingly. Facilities must address both short- and long-term needs for specially designed products and equipment for obese patients. However, the majority of US healthcare facilities are under-equipped to accommodate the growing number of bariatric patients, and many cannot afford an overhaul of equipment upgrades.

Treating the whole patient

Retrofitting your facility to be bariatric-ready not only addresses the physical needs of your obese patients, but their emotional well-being as well. Obese patients often avoid treatment based on access and sensitivity to their healthcare environment. When you comply with guidelines for your bariatric patients, you treat the whole patient.

BARIATRIC-READY HOSPITAL ROOM



Product #1047 Publication #6617

Shot on location at Radius Hospital, West Roxbury, Massachusetts

What can you do?

You can take advantage of several low-cost methods throughout this catalog to make your facility bariatric-ready. Our bariatric products accommodate the individual needs common to heavier patients.

How can we help?

We can help you and your facility become bariatric-ready. Throughout this brochure we have highlighted problem areas within facilities where the proper equipment is not always readily available. We can help answer any questions you may have and supply you with products to save your facility money, make your job easier, and keep your patients' dignity and safety intact.

BARIATRIC-READY BATHROOM



Bariatric Seating

Safety Bars

SAFETY BAR

BARIATRIC-RATED TOILET SEATS

TOILET SUPPORT

BED

BARIATRIC SEATING

MOBILITY AIDS

COMMUNE

Toilet Support

Reinforce your wall-mounted toilets

The industry standard* requires that all wall-mounted toilets pass a 500-lb. static load test to avoid any breakage that could result in injuries. Make your wall-mounted toilet safer for bariatric patients. *(ASMEA112.19.2)

Solution:

Bariatric Toilet Support

Discreetly mounts underneath toilet bowl to redistribute the weight and increase weight capacity from 350 lbs. to 1,000 lbs. High-impact ABS and stain- and chemical-resistant PBC. Easy fit, easy clean. Fits under most wall-mounted toilets. Meets SRC criteria for Bariatric Centers of Excellence. Choose between Standard (height adjusts 7½"-14"H) and Short (height adjusts 5" 7½"H).



Height adjustment from 7½" to 14" #CSS71 0145



		6+	EA
#CSS71 0145R	Standard	\$175.00 ea	\$ 202.00 ea
#CSS71 1635R	Short	170.00 ea	185.00 ea

Bariatric Toilet Seat

Rated to hold up to 1,200 lbs. Wider (19"), 2" higher, and more comfort-able than standard toilet seats. Fits standard and elongated toilets. Un-breakable with extra-strong stainless steel hinges. Large bumpers grip bowl and keep seat from sliding. Available with closed or open front. White.



		6+	EA
#CSS71 0144R	Closed-seat w/Cover	\$150.00ea	\$174.00 ea
#CSS 71 1311R	Open-seat w/Cover	150.00 ea	171.00 ea
#CSS 71 1310R	Open-seat w/o Cover	145.00 ea	160.00 ea

Independent sitting and standing from toilet

Standard toilets are 14"-16" from the ground, making it difficult for heavy patients to lift themselves off the toilet.

Solution:

Easy Toilet Riser

Platform raises entire toilet by 4". Blends in with toilet. Sturdy and durable. Includes universal modification kit. Lifetime guarantee. White. 500-lb. capacity.



	6+	EA
#CSS71 1974R	\$115.00 ea	\$ 140.00 ea

Toilet Seat Riser Easy and effective way to increase height on standard toilet. Installs at seat hinges. Easy to clean. Raises seat height by 2¼". Permanent Installation. Hardware included. 500-lb. capacity.



	6+	EA
#CSS75099R	\$34.25 ea	\$44.25 ea

Easier bed pan insertion

Most bed pans are 3"-4"H, making it difficult for the caretaker to insert and uncomfortable for the patient.


Solution:

AliMed® Bed Pan is only ½" H in the back of the pan and is certified up to 1,200 lbs., holds 2 quarts of liquid, and serves a wide variety of patient needs including:

- **Bariatric patients.** Unlike traditional bed pans, this pan is ideal for bariatric patients because its low-profile design allows it to easily slide underneath patient for easy placement and removal and better comfort for patients.
- **Fracture patients.** The Universal Bed Pan's low profile reduces the need to lift the patient up high for insertion and removal, aiding patient comfort and simplifying the task for caregivers. The covered area of the pan under the flat seat better contains fluids, reducing spill and splash hazards.
- **Labor and delivery patients.** AliMed's Universal Bed Pan fills the necessity in obstetrics departments of catching the placenta: a challenge for most traditional bed pans.

- ✓ Wider seat for more comfort
- ✓ Superior splash control

**UP TO
1,200-lb.
CAPACITY**



- Bed Pan's lower profile:
- ✓ Easy placement
 - ✓ Easy removal



Holding rack

#CSS71 1255R	Bed Pan	6+	EA
		\$17.95 ea	\$24.50 ea
#CSS71 1273R	Universal Holding Rack	65.00 ea	75.00 ea

Bathroom aids for nonambulatory patients

Walking to the bathroom, especially at night, can be difficult for bariatric patients, as is standing up in the shower/bathtub. This is due to issues with their joints and stamina.

Solution:

Bariatric Steel Commode

Great for people who have difficulty walking to the bathroom. Deeper seat depth for larger individuals. Height-adjustable legs accommodate different heights. Steel. 12-qt. bucket and cover. Plastic armrests. Comes assembled. Limited lifetime warranty. **OS**

31"W x 17"-23"H • 24" between arms • Seat: 14"W x 16½"D • 21 lbs. • 650-lb. capacity

**650-lb.
CAPACITY**



#CSS70197R	6+	EA
	\$160.00 ea	\$175.00 ea

Bariatric Shower Commode

Easily transfer bariatric patients into the shower/bath. Includes built-in splash guard. Seat snaps on/off for cleaning. Sliding footrest included. **OS**

30"W x 40"H x 26"D • 4" heavy-duty wheels • 7-qt. pail • 36 lbs. • 500-lb. capacity

**500-lb.
CAPACITY**



#CSS72259R	6+	EA
	\$320.00 ea	\$357.75 ea

Specify Blue, Mauve, or Putty.

Accomodate any size visitor

Standard facility chairs and couches only hold up to 250 lbs.
As the number of heavier patients grows, so does the number of heavier visitors.

Solution:

Austin Bariatric Chairs Solid hardwood frame and extra-firm foam. Open-arm design. Antimicrobial-treated Crypton upholstery. Matching ottomans sold separately below. 1,000-lb. capacity.

		6+	EA
#CSS72146R	28"W x 33"H x 25"D	\$1,616.00 ea	\$1,667.00 ea
#CSS72147R	34"W x 34"H x 25"D	1,699.00 ea	1,756.00 ea
#CSS72148R	44"W x 34"H x 25"D	1,773.00 ea	1,844.00 ea
#CSS72009R	24"W x 17"H x 18"D	432.00 ea	473.00 ea

Frame Finish

Oak	Walnut
-----	--------

Austin Crypton Upholstery

Aqua	Blue	Champagne
Merlot	Okra	



ConvaQuip® Bariatric Arm Chair High-pressure laminated northern hardwood. Adds 5" of adjustment in 1" increments, factory installed only. Fire Rating: Cal TB 133. Warranty ensures product is free of defects in materials and workmanship for 10 years. Marina blue. 750-lb. capacity.

		6+	EA
#CSS78513R	25"W	\$668.00 ea	\$708.00 ea
#CSS78514R	25"W w/Leg Extensions	768.00 ea	808.00 ea
#CSS78515R	30"W	710.00 ea	750.00 ea
#CSS78516R	30"W w/Leg Extensions	785.00 ea	825.00 ea
#CSS78517R	40"W	768.00 ea	808.00 ea
#CSS78518R	40"W w/Leg Extensions	835.00 ea	875.00 ea



Carstone Bariatric Chairs are ideal for any medical or personal application and are some of the most durable on the market. Each chair frame has been static-load-tested to 6,000 lbs. Features a 1 1/4" tubular steel frame which is backed by a lifetime frame warranty, as well as durable Sherpa fabric and extra dense foam for extra support. 1,000-lb. capacity.

		6+	EA
#CSS71990R	20"W	\$ 981.00 ea	\$1,081.00 ea
#CSS71991R	24"W	901.00 ea	1,101.00 ea
#CSS71992R	30"W	1,053.00 ea	1,153.00 ea
#CSS71993R	40"W	1,177.00 ea	1,277.00 ea

Specify Frame and Fabric color. See chart below.

Fabric Options

Academy Blue	Hunter Green	Black	Grape	Maroon	Amber
--------------	--------------	-------	-------	--------	-------

Frame Options

Winter Grey	Pale Brown
-------------	------------

Easily manage daily care

Dignity is of great importance when caring for your obese patients. Make sure they have equipment that is large enough to accommodate them and their daily care, while still maintaining dignity.

Solution:

Blood Pressure Cuff

features chrome-plated aneroid no-pin stop gauge with lifetime calibration warranty. Bold white numerals on contrasting blue face-plate for easy reading in all light conditions. Heavy-duty nylon cuff has size, index, and range markings, artery indicator label and gauge holder. Maxiflo® deluxe air-release valve with built-in filter and extra-large knurled air release knob for precise deflation control. Double filters eliminate dust particles. Two-tube latex bladder. Measures from 20 mmHg to 300 mmHg, certified accurate to ±3 mmHg. Zippered carrying case.



	Size	Dimensions	Circum. range	6+	EA
#CSS95-186R	Adult Lg. (720X)	6.9"W x 24.9"L	24.9-50.8 cm	\$55.00 ea	\$71.00 ea
#CSS92 4978R	Thigh (700T)	8.3"W x 32.5"L	40.6-65.5 cm	65.00 ea	76.50 ea

Portable Wheelchair Scale

with sloped entrance/exit ramps to accommodate wheelchairs, straight-legged chairs, or walk-on patients. **OS**

Platform: 28"W x 32"D • 800-lb./360-kg capacity in 1/5-lb./100-g increments
 • Functions: enhanced 785C digital weighing, BMI calculator, serial connectivity to a printer or PC, push-button and keypad tare modes, preset or manual key tare weight, motion detection, weight lock-in
 • Operates on six D batteries (included)
 • 65 lbs. • White • Ships freight
 • Two-year warranty



500-lb. CAPACITY

DETECTO

#CSS93 0646R 6+ EA \$2,624.75 ea \$2,824.75 ea

Folds and latches for easy transport!

NEW! **SoftPro™ Bariatric In-Bed AFO Supports**, protects, and positions ankle to help prevent plantar-flexion contracture. Footdrop heel protection accommodates larger calf and foot. Suspended heel, anti-rotation bar, and adjustable toe post assist with skin protection and positioning. Low-profile back and plush fleece liner ensure superior comfort. Nonskid bottom for assisted transfers.
 One size: Fits calf circumference 21"-32". Footplate 10"L, 14"L w/toe post.



#CSS66400R 6+ EA \$95.00 ea \$120.00 ea

NEW! **Coloplast InterDry Ag™ Textile with Antimicrobial Silver Complex** is a skin fold management system for moisture, odor, and inflammation in skin folds and other skin-to-skin contact areas. Textile functions by wicking moisture away to keep skin dry and offers a friction-reducing surface to help minimize risk



Helps prevent skin-to-skin contact


of skin tears and pressure ulcers. Impregnated with silver to reduce the odor associated with skin folds. Textile significantly improves symptoms associated with intertrigo, maceration, denudement, inflammation, itching, and satellite lesions. 10" x 12' Roll. 10rls/box.

#CSS66360R 6+ EA \$815.00 bx \$915.00 bx

Bariatric Overbed Table Oversized table top measures 18" x 40" to reach the center of most bariatric beds. **OS**

Steel welded frame • Table top height adjusts from 30" - 45" in 1" increments • Swivel casters • Limited lifetime warranty • 70 lbs. • 500-lb. capacity

#CSS71938R 6+ EA \$575.00 \$625.00 ea



500-lb. CAPACITY

Geri-Sleeves Breathable cotton-blend LYCRA® protects against skin tears. Thumb-hole cutouts help keep arm sleeves in place. Leg sleeves help prevent skin shear at heels and calves. Protect IV and wound sites. Reusable and washable.

Bariatric Geri-Sleeves Sized to comfortably fit bariatric patients.

Sizing: Arm sleeve: S measures 3"W x 17"L; M measures 4"W x 17"L. L/Bariatric measures 5"W x 19"L; Leg sleeve measures 3 1/2"W x 17"L.

	Geri-Sleeves	Size	6+	PR
#CSS50 3357R	Arm Sleeves	L/Bariatric	\$12.00	\$14.00 pr
#CSS50 3360R	Leg Sleeves	Universal	12.00	14.00 pr


Suggested code: A9270



Optimize skin protection and comfort

When mattresses bottom out due to weight capacities being exceeded, it introduces the opportunity for skin ulcers to grow.

Solution:



Bariatric Mattresses Three-layer foam mattresses help prevent formation of pressure ulcers. High-resiliency foam top layer allows surface to mold to patient's body, middle layer adds firmer support, and bottom layer provides proper orthopedic support and prevents bottoming out. Laminated Nylex II cover offers excellent stain and fluid

resistance; heavy-gauge vinyl bottom is slip-resistant. Ticking and foam are completely antimicrobial to help prevent hospital-acquired infections and eliminate odors. Includes fire barrier. 80" L. **OS**


	Dimensions	6+	EA
800-lb. capacity			
#CSS71 1792R	39"W x 6"H	\$425.00 ea	\$495.00 ea
#CSS71 1793R	42"W x 6"H	425.00 ea	495.00 ea
#CSS71 1794R	48"W x 6"H	525.00 ea	595.00 ea
#CSS71 1795R	54"W x 6"H	525.00 ea	595.00 ea
1,000-lb. capacity			
#CSS71 1796R	39"W x 7"H	\$625.00 ea	\$695.00 ea
#CSS71 1797R	42"W x 7"H	625.00 ea	695.00 ea
#CSS71 1798R	48"W x 7"H	725.00 ea	795.00 ea
#CSS71 1799R	54"W x 7"H	725.00 ea	795.00 ea

Advantage BARI Mattress Premium pressure-relief mattress for bariatric residents. Anatomically designed to conform to body for outstanding pressure redistribution and comfort. Sculpted top layer uses cropped foam to create air channels, reducing buildup of heat and moisture.

Super-soft layered heel section for maximum pressure redistribution. Nylex top fabric combines durability, performance, and shear reduction. 80" L. **OS**

	Dimensions	6+	EA
800-lb. capacity			
#CSS71 1800R	39"W x 6"H	\$525.00 ea	\$595.00 ea
#CSS71 1801R	42"W x 6"H	525.00 ea	595.00 ea
#CSS71 1802R	48"W x 6"H	625.00 ea	695.00 ea
#CSS71 1803R	54"W x 6"H	625.00 ea	695.00 ea
1,000-lb. capacity			
#CSS71 1804R	39"W x 7"H	\$725.00 ea	\$795.00 ea
#CSS71 1805R	42"W x 7"H	725.00 ea	795.00 ea
#CSS71 1806R	48"W x 7"H	825.00 ea	895.00 ea
#CSS71 1807R	54"W x 7"H	825.00 ea	895.00 ea

Need a replacement?



#CSS71 1684 Replacement Bariatric Mattress

Replacement Bariatric Mattresses Available in either 500- or 900-lb. weight capacities, these mattresses offer bariatric patients comfort and help reduce pressure ulcer risk. Top polyurethane-coated fabric provides omni-directional stretch. Will not bottom out under patient's weight. Three-part foam core uses different foams: Base foam delivers firmness, Intermediate layer is cupped at the sides for strength and side reinforcement, and Central foam is a softer patterned foam for pressure distribution.

500-lb. capacity

Item #	W x L x Thick	6+	EA
#CSS71 1684R	39" x 80" x 7"	\$393.00 ea	\$493.00 ea
#CSS71 1685R	39" x 84" x 7"	410.00 ea	510.00 ea
#CSS71 1686R	42" x 80" x 7"	469.00 ea	569.00 ea
#CSS71 1687R	42" x 84" x 7"	490.00 ea	590.00 ea
#CSS71 1688R	48" x 80" x 7"	548.00 ea	648.00 ea
#CSS71 1689R	48" x 84" x 7"	568.00 ea	668.00 ea
#CSS71 1690R	54" x 80" x 7"	570.00ea	670.00 ea
#CSS71 1691R	54" x 84" x 7"	590.00 ea	690.00 ea

900-lb. capacity

Item #	W x L x Thick	6+	EA
#CSS71 1692R	39" x 80" x 7"	\$608.00 ea	\$708.00 ea
#CSS71 1693R	42" x 80" x 7"	647.00 ea	747.00 ea
#CSS71 1694R	42" x 84" x 7"	671.00 ea	771.00 ea
#CSS71 1695R	48" x 80" x 7"	772.00 ea	872.00 ea
#CSS71 1696R	48" x 84" x 7"	805.00 ea	905.00 ea
#CSS71 1697R	54" x 80" x 7"	858.00 ea	958.00 ea
#CSS71 1698R	54" x 84" x 7"	891.00 ea	991.00 ea

Replacement covers also available—call 800.225.2610 for details

Independence for self-repositioning

Patients need to be repositioned every two hours to avoid bed sores. Due to limited staffing in facilities, repositioning aids are available to help patients move themselves.

Solution:



1,000-lb. CAPACITY

#CSS78212

Bariatric Trapeze Lightweight, tool-free, four-piece design. Three locking overhead boom angle positions for comfortable, secure use.

Overall height: 73" • Reach: 35½" • Base width: 40" • Base length: 50½" • Product wt.: 72 lbs. • 1,000-lb. capacity • 3-year warranty

	6+	EA
#CSS78212R	\$1,390.00 ea	\$1,490.00 ea

Ships freight.

Bariatric Bed Handle mounts directly to bed frame (metal bed frame required) for extra support and safety. No tools required. Top of handle is approximately 23" above the mattress. **OS**

#CSS8639R	6+	EA
	\$135.00 ea	\$155.00 ea

Rail extends to floor for greater support



600-lb. CAPACITY

Halo Safety Ring Side Rail Circular design and vertical mounting bar reduce entrapment concerns and allow exit points for the neck, head, and chest. Rotates 90° in either direction. Durable steel construction. Includes easy-to-install mounting hardware and stay bar to secure mattress.

21"H (Bed rail clamp to top of ring) • 12" ring diameter • Fits most hospital style beds • 1,000 lb. capacity.

#CSS71 1770R	6+	EA
	\$400.00 ea	\$425.00 ea



1,000-lb. CAPACITY

Respiratory sleep aids

Overweight patients often struggle to breathe due to the obstruction and/or collapse of the upper airway (throat). Sleeping aids help them get a restful night's sleep.

Solution:

RestEZ III CPAP internal pressure transducer produces real-time measurements to ensure accurate delivery of prescribed pressure settings. Built-in Auto Compliance feature monitors usage and comes complete with USB cord, software which allows usage information to be easily downloaded to check compliance. Also has an interface alert that notifies patient every 90 days that to replace interfaces.



NEW!

Shown with optional heated humidifier

#CSS93 6054R	RestEZ III CPAP System	\$600.00ea	\$658.00 ea
#CSS93 6055R	Optional Heated Humidifier	175.00 ea	202.50 ea
#CSS98 0043R	I.Q. Adjustable Gel Mask w/Headgear, one size	147.00 ea	157.00 ea
#CSS98 0194R	CPAP Mask w/Vented Swivel, L	65.00 ea	68.25 ea
#CSS98 0197R	72" Reusable CPAP Circuit	7.00 ea	9.00 ea
#CSS98 0198R	18" CPAP Tubing, Rubber Molded Cuff	9.00 ea	11.00 ea
#CSS98 0199R	Vented Swivel Elbow Adapter (for All Contour Masks)	64.00 ea	74.75 ea

Suggested Code: E0601

CPAP Sleep Aid Pillow

improves mask comfort and freedom of movement during sleep. Multi-position for side, back, and stomach sleepers. Offers neck support and spinal alignment. Reduces mask leaks, mask pressure on face, airway resistance due to better neck support. Works with all major brands and styles of CPAP masks. Mask, hose, and straps not included.



#CSS82736R	6+	EA
	\$51.00 ea	\$61.00 ea

Safer repositioning for patients and staff

Patients often slide down their beds and need to be shifted and turned to avoid skin irritation and ulcers. Using the right shifter or booster helps reduce injury risk for both patient and staff.

Solution:

“No-Lift” Bariatric BOOSTER™

Boost patients up to 600 lbs.! Helps ease caregiver back strain by making resident easier to move toward the headboard. The BOOSTER is a flat tube of fabric with a virtually frictionless interior surface, and handles on each side. Simply position resident on top of tube just as you would with a draw sheet. Now grab the reinforced edges of top tube layer and slide it upwards. No lifting is needed—less risk of back strain! The patient’s skin integrity is never compromised since BOOSTER eliminates draw sheet drag and shear. Remove from underneath resident after each use. Comes with handles.



Flat dimensions: 46"W x 46"L • Black • Machine-washable, air-dry • Disinfectant-safe • 600-lb. capacity

	6+	EA
#CSS75582R	\$60.00 ea	\$103.00 ea

“No-Lift” TURNER™ Simply position patient on top of the tube just as you would with a drawsheet. Grab reinforced edges of the top tube layer and pull it toward you. No lifting and less back strain! Once TURNER is put in position (generally requires two people), only one staff member is needed to move most sidelying residents.



Since the TURNER eliminates drawsheet drag and shear, the resident’s skin integrity is never compromised. Comes with handles.

Flat dimensions: 29"W x 52"L • Black • Machine-washable, air-dry • Disinfectant-safe • 600-lb. capacity

	6+	EA
#CSS712177R	\$60.00 ea	\$105.00 ea

How the BOOSTER works



Place “No-Lift” BOOSTER under resident as you would a drawsheet.



Roll resident onto one side and place BOOSTER underneath.



Once BOOSTER is positioned slide resident up toward head of the bed by grabbing top layer of fabric.

How the TURNER works



“No-Lift” TURNER is positioned across bed to aid in turning schedules.



Roll resident onto side to position TURNER under patient’s hips and shoulders.




Once TURNER is positioned, patient can be safely repositioned by one staff person.



For more information, and to see a video, go to <http://video.AliMed.info> and click on “No-Lift” Booster and Turner.

Note: If patient has highly sensitive skin, BOOSTER or TURNER should be removed between uses.



Bariatric Shifter Lateral transfers for patients up to 600 lbs. Supports patient’s entire body to help reduce pinching of skin. Helps reduce caretaker back strain by minimizing friction on transfer surfaces. Elongated hand holds are ergonomically correct. Flexible plastic is great for uneven surfaces. Recommend using one caretaker per 100 lbs. 28"W x 72"L x 3/16"H. 16 lbs. 600 lb.-capacity. Purple. *

	6+	EA
#CSS93 5944R	\$220.00 ea	\$320.00 ea

***Save On Shipping** Shifter ships oversize at same rate for 1–10 units. Buy up to 10 and only pay one shipping charge.


Instantly recognize the weight capacity of your equipment

Weight Capacity Labels give quick recognition of equipment weight capacity, while maintaining patient dignity. Labels also help facilities meet criteria for SRC Centers of Excellence for Bariatric Centers.* 1½" x 1½" vinyl labels printed in large type that can be seen from 6-8 feet away. Patient dignity is improved by subtly stating weight capacity; 400-lb. capacity label shows an image of a dove and the number 40. **Wood:** Chairs, crutches, rehabilitation treatment tables, etc. **Metal:** Wheelchairs, walkers, canes, patient lifts, etc. **Plastic:** Patient shifters, transfer boards, etc. Labels can be sprayed and wiped with disinfectant

#CSS71 1511R Variety Pack Labels: 175/pk: (25 labels each of 30, 40, 50, 60, 70, 80, 100) \$60.00 pk (6+) \$69.90 pk (EA)



*Surgical Review Corporation (SRC) requires that facilities have equipment available for morbidly obese patients and that all staff members can easily recognize the weight capacity for each piece of equipment.



Higher-weight-capacity lifts and slings

Just because you have a lift in your facility does not mean the weight capacity is higher than the standard 250 lbs. Check to make sure your lift is appropriate for bariatric patients.

Solution:



Heavy-Duty Electric Patient Lift Emergency manual lowering device, universal battery, lift-capacity warning beep, six-point spreader bar and foot-pedal base opening for wheelchair/ commode.

APEX

Six-point multi-tilt spreader bar.

50"L x 48"H • Spreader bar: 22"-71"H • Base closed: 26"W (external); 21½"W (internal) • Base open: 43"W (external); 38"W (internal) • Overall base: 6"H • Front and Rear casters: 4" diam. • Lifting mechanism: 24 VDC Linear Actuator • Battery: 12V 4½ AH x 2 • Battery charger: AC 110V/220V to DC 24V 1 AMP • 137 lbs. • 600-lb. capacity

#CSS77820R 6+ EA
\$2,477.00 ea \$2,577.00 ea
Suggested code: E0630

Stella Bariatric Stand Assist Lift Quick and easy access for cleaning and changing. Less time required for lift and transfer. Secure and safe for both patient and caregiver. Doubles as a commode lift for additional flexibility. Height-adjustable leg and knee support. Detachable foot platform. Foot pedal base opening. 40"L x 68"H • Spreader bar: 33"-68"H • Base closed 261/2"W (external); 213/4"W (internal) • Base open 37"W (external); 321/2"W (internal) • Overall base: 41/2"H • Front caster w/brake: 6" diam. • Back caster w/brake: 6" diam. • Lifting mechanism: 24 VDC Actuator • Battery: AC 110V/220V to DC 24V 1 Amp • 122 lbs. • 600-lb. capacity

#CSS79852R
Suggested code: E0635

6+ EA
\$2,999.00 ea \$3,219.25 ea



AliMed Bariatric Slings

Width: 45" • Length: 55" • 600-lb. capacity



UP TO 600-lb. CAPACITY

Full-Body, Mesh Seat

Full-Body Slings	Item #	6+	EA
Mesh Seat	#CSS78258R	80.00 ea	110.00 ea
Mesh Seat, Commode	#CSS78259R	100.00 ea	140.00 ea
Universal Split-Leg Slings			
Padded Seat	#CSS72080R	146.00 ea	176.00 ea
Stand Assist Slings			
Stand Assist Sling	#CSS79853R	157.00 ea	187.00 ea
Stand Assist Buttock Strap	#CSS79854R	68.00 ea	88.25 ea

Safely ambulate overweight patients

This reinforced assistive device is used to help safely transfer a person in a weakened state, and also helps caregivers avoid extra back strain.

Solution:

Bariatric Gait Belt

Ease caregiver back strain when ambulating with bariatric patients. 96" wipeable gait belt easily wipes clean with any commonly used disinfectant. Reinforced plastic strap is easy to grasp and comfortable for the patient. Specialty buckle ensures belt stays attached. 96". Yellow. 1,000-lb. capacity.

#CSS71 1779R 96" Gait Belt \$20.00 ea (6+) \$26.00 ea (EA)

1,000-lb. CAPACITY




Wheelchair cushions that won't bottom out

Many patients can spend up to 8 hours a day in a wheelchair. Adding a cushion is a great option to optimize comfort, but it is imperative that the cushion has enough support that it does not bottom out on the wheelchair.

Solution:

AliMed® Bariatric Tri-Foam Cushion offers pressure relief for patients up to 1,000 lbs. and wheelchairs up to 30" wide. Durable cushion is made up of three distinct layers: a 2" top layer of T-Foam memory foam for exceptional pressure distribution; a 1" middle layer of open-cell foam for pressure reduction; and a 1" bottom layer of closed-cell foam to help reduce pressure and to help prevent "bottoming out". Black, washable, waterproof cover with ties and carry handle. Larger and custom sizes available.



This cushion can be
Call 800-225-2610
Customized

**UP TO
1,000-lb.
CAPACITY**

		6+	EA
#CSS1795R	20"W, 22"W, or 24"W <i>Specify 20"W x 18"D; 22"W x 18"D; 24"W x 18"D; or 24"W x 20"D.</i>	\$155.00 ea	\$185.25 ea
#CSS1796R	26"W or 28"W <i>Specify 26"W x 18"D; 26"W x 20"D; 28"W x 18"D; or 28"W x 20"D.</i>	\$165.00 ea	\$195.25 ea
#CSS1797R	30"W <i>Specify 30"W x 18"D; or 30"W x 20"D.</i>	\$175.00 ea	\$196.25 ea

Wheelchairs to fit your largest patients

Standard wheelchairs are only 18"W. Wheelchairs with much larger seats and weight capacities fit larger patients and save their dignity.

Solution:

Sentra EC Heavy Duty Wheelchair Reinforced steel gussets at all weight-bearing points. Easy transition of seat to hemi-height. Steel spoke wheels. Heavy-duty front forks. Carbon steel frame. Push-to-lock wheels. Elevating Swing-Away Legrests sold separately below. **OS**

26" Seat Width: Overall width (open-closed): 34½"-14" • Overall height: 39" • Seat-to-floor height: 17½"-19½"
• Back height: 18" • Seat depth: 20" • Overall Length (with riggings): 43" • 87 lbs. • Ship weight: 102 lbs. • 700-lb. capacity

28" Seat Width: Overall width (open-closed): 37"-14½" • Overall height: 39" • Seat-to-floor height: 17½"-19½"
• Back height: 18" • Seat depth: 20" • Overall length (with Riggings): 43" • 88 lbs. • Ship weight: 103 lbs. • 700-lb. capacity

30" Seat Width: Overall width (open-closed): 39"-15" • Overall height: 39" • Seat-to-floor height: 17½"-19½"
• Back height: 18" • Seat depth: 20" • Overall length (with riggings): 43" • 89 lbs. • Ship weight: 104 lbs. • 700-lb. capacity

AliMed #	Mfr. #		6+	EA
#CSS72338R	STD26ECDDA	26"W Seat	\$727.00 ea	\$767.00 ea
#CSS72339R	STD28ECDDA	28"W Seat	727.00 ea	767.00 ea
#CSS72340R	STD30ECDDA	30"W Seat	849.00 ea	909.00 ea
#CSS710410R	PHELR	Elevating Swing-Away Legrests	164.00 ea	184.00 ea



**700-lb.
CAPACITY**

drive

Proper wheelchair positioning

Proper wheelchair positioning helps promote good posture, enhance breathing and digestion, prevent complications like pressure sores, minimize pain, and maximize functionality.

Solution:



Skin-Guard Leg Protector guards legs against skin tears and bruises. Soft, durable neoprene pads wrap around wheelchair footrest support bars. Hook-and-loop tabs. Universal fit. Pair.

	6+	EA
#CSS70 3073R	\$24.00 pr	\$28.00 pr

Extended-Width Armrest

accommodates individuals with upper extremity edema, flaccidity, or increased tone while eliminating risk of skin tears. Individual can shift position normally while avoiding risk of slippage and damage to the arm and shoulder joint.

For petite patients, armrest can be placed in opposite configuration for upper extremity support and weight-bearing capability. Universal fit. 6"W x 17"L.



	6+	EA
#CSS77382R	\$43.00 ea	\$63.50 ea



Foot support swings away.



Swing-Away Foot Supports attach to wheelchair footrests. Cushioned side and back panels help prevent feet from slipping off the footrest. Covered in wipe-clean nylon. Fit 22"-30"W wheelchairs. Pair.

	6+	EA
#CSS70 3477R	\$115.00 pr	\$126.75 pr

Higher-weight-capacity mobility aids

Most standard mobility aids have a 200-lb. weight capacity. There are plenty of bariatric options from AliMed!

Solution:

Bariatric Large Base Quad Cane

1" aluminum tubing. Steel base. Lightweight. One-button height adjustment from 30"-39". Vinyl hand grips. Limited lifetime warranty. 500-lb. capacity

	6+	EA
#CSS82326R	\$49.00 ea	\$59.00 ea

500-lb. CAPACITY



Bariatric Sidewalker for larger individuals with use of only one hand or arm. For users who do not require a walker but need a wider base for support. Lighter than a walker and more stable than a cane. Folds easily with one hand. Height adjusts from 32"-35".

16½"W • 36"H (highest setting) • 16"D • 30" (lowest setting) • 10 lbs. • 500-lb. capacity • Limited lifetime warranty

	6+	EA
#CSS70622R	\$156.00 ea	\$176.00 ea

500-lb. CAPACITY



Drive Medical Bariatric Cane with vinyl hand grip and wrist strap. 1" steel tubing. Height-adjustable handle. 500-lb. capacity

	6+	EA
Adult: Silver, adjusts 27¾"-36¾"		
#CSS82328R	\$17.00 ea	\$22.75 ea
Tall Adult: Silver, adjusts 37"-46"		
#CSS82329R	21.00 ea	30.75 ea

500-lb. CAPACITY



Drive Medical Bariatric Forearm Crutches with vinyl-coated arm cuffs. Leg and forearm sections adjust independently. Extra large tips for added stability.

Youth: Adjusts 21½"-30½"H, fits 4'2"-5'2"
Adult: Adjusts 28"-37"H, fits 5'0"-6'2"
Tall Adult: Adjusts 33"-42"H, fits 5'10"-6'6"
 • Limited lifetime warranty • 500-lb. capacity

	6+	EA
#CSS71 1975R	\$221.00 pr	\$251.00 pr

Specify Youth, Adult, or Tall Adult.

500-lb. CAPACITY



Extra-Wide Universal Dualie

Tuff-Coat® nonslip grip. Dual 5" wheels for low strength ambulation. Rear 2" swivel casters have push-down auto brakes. Over-the-toilet design with angled front legs to ensure stability. 29"W x 24"D • 31"-41" adjustable height • 11½ lbs. • 500-lb. capacity

	6+	EA
#CSS74862R	\$226.00 ea	\$246.75 ea

500-lb. CAPACITY



Bariatric Walker Platform Attachment Cradles forearm with soft vinyl padding. Large top bracket allows platform attachment to be placed over walker hand grip. Adjustable strap secures arm in place. Limited lifetime warranty.

	6+	EA
#CSS70586R	\$99.00 ea	\$134.00 ea

500-lb. CAPACITY



To Order

- MAIL AliMed® inc.** Pub #6617
Order Department Product #1047
297 High St., Dedham, MA 02026
- PHONE 800-225-2610**
8:00 AM-7:00 PM EST, Mon.-Fri.
- 24-HOUR FAX 800-437-2966**
- International: 781-329-2900 FAX 781-329-8392**
- ONLINE CATALOG**
24-hr. web ordering www.AliMed.com
- CHARGE IT** By phone, mail, or fax. Please see our new Shipping Rate Table below for exact charges. Credit cards are only accepted at time of order entry.
- QUESTION OR RETURN?** Customer Service, x6901

MINIMUM BILLING Orders are subject to a \$25.00 minimum merchandise billing plus shipping, handling, and insurance. International orders (except Canada) require a \$300.00 (U.S.) minimum order.

TERMS Net 15 to established accounts. All prices F.O.B. shipping point. Shipment made at prevailing price at time of order. All prices in U.S. \$. Prices are professional wholesale. Patient may order directly; retail pricing applies. Prices, availability, and specs subject to change. Visit www.AliMed.com for the most current pricing. Items may be discontinued at any time. AliMed reserves the right to substitute similar products to those shown in catalogs at any time. Overdue accounts are subject to a late payment fee of 2.0% per month, resulting in an annual rate of 24%.

Quantity Discounts: Item numbers, sizes, colors, Lefts and Rights, etc. cannot be assorted for a quantity discount. One item #, one style per quantity discount order (unless otherwise specified).

SHIPPING We ship FedEx/UPS whenever possible. Charges are prepaid and added to invoice. Freight shipments are freight collect. See Shipping Rate Table for exact charges.

OS Oversized items subject to additional freight charges.

OVERSIZED CHARGES Apply to all orders including contract orders, items with free ground shipping, and free shipping promotions.

GLN #0733657000009

TAX Applies to all orders shipped in MA, SC, WA, and SD. If you are tax exempt, please forward a copy of your tax exempt or resale certificate. If your facility is not in MA, SC, WA, or SD, please submit use tax directly to your state.

INSURANCE Your shipment is 100% insured against loss or damage.

RETURNS Our return policy is consistent in maintaining the best value and highest service level for all of our customers. If the merchandise you received is defective or not as you ordered, please contact our customer service staff for quick resolution.

All other returns will be subject to a 20% restocking fee (\$10.00 minimum) except "Special Order" items that have been "drop shipped". (These items will be subject to a restocking fee that may vary at the discretion of the individual vendor.)

The guidelines below must be followed in order to receive any credit:

- Merchandise credit only
- Product must be in new, resaleable, unmarked condition and in original packaging
- Merchandise must be returned within 30 days of invoice date
- Custom products, beds, books, videos, shoes, eating utensils, and toileting products are not returnable
- Open packages cannot be returned for credit

To ensure prompt return processing, please attach the original packing list, invoice or RA number to the return. Refunds will be for merchandise value only. Please allow three to four weeks for return processing.

TRUCK SHIPMENTS Delivery Inspection Required

IMPORTANT NOTICE: All truck shipments must be fully inspected by receiver before signing for the delivery to ensure there is no concealed damage. If you cannot inspect immediately, make a note "pending inspection" with your signature. Any damaged shipment received should be refused by the receiver. These guidelines must be strictly followed in order to receive any credit on Truck shipment orders. **Overnight Truck Shipments - AliMed does NOT** provide overnight delivery on orders that require delivery by truck.

DAMAGED GOODS/SHORT SHIPMENTS Must be reported within 10 days of receipt. Be sure to note any damage or poor box conditions when signing for packages. If shipment arrives damaged or with parts missing, contact the carrier immediately. Be sure to save product packaging. All shipping charges and risk of loss are FOB shipping point and are the responsibility of the purchaser.

Bariatric aids for daily living

No one likes to ask for help!
Here are a few tools to help patients can take care of themselves.

Solution:

Inspection Mirror




	6+	EA
#CSS8255R	\$19.00 ea	\$25.50 ea

Freedom Wand™




	6+	EA
#CSS82850R Wand	\$40.00 ea	\$50.00 ea
#CSS82851R Extension	15.00 ea	20.25 ea

Easy-Slide Shoehorn



#CSS8626		#CSS8625
	6+	EA
#CSS8625R 18"	\$4.00 ea	\$6.75 ea
#CSS8626R 24"	4.50 ea	7.75 ea

Xtra Wide Sock/Stocking Aid



	6+	EA
#CSS77816R	\$ 15.00 ea	\$ 18.25 ea

Lightweight Reachers



	6+	EA
#CSS82598R 26"L, 6 oz.	\$9.00 ea	\$13.25 ea

Big Hook™ Dressing Stick



	6+	EA
#CSS8202R 28"L	\$4.00 ea	\$6.75 ea

For more ADLs visit www.AliMed.com

Shipping charges will be added to your total. **OS** Oversized items subject to additional freight charges.

SHIPPING RATE TABLE**

Includes shipping and handling charges • Temporary fuel surcharge may be added to our published rates
No small order charge • Instant freight quotations • No additional freight on backorders

Some large volume orders and/or oversized items may require additional freight charges or truck shipping.

Ground/Air (FedEx, UPS, Postal)

Merchandise value	Continental US			Outside cont. US*		
	Ground	Next Day	2nd Day	Ground	Next Day	2nd Day
\$ 0-99.99	\$11.99	\$39.99	\$25.99	\$26.99	\$59.99	\$40.99
100-199.99	16.99	49.99	29.99	35.99	69.99	48.99
200-499.99	19.99	57.99	35.99	42.99	86.99	55.99
500-799.99	28.99	64.99	42.99	49.99	93.99	63.99
800+	34.99	69.99	48.99	55.99	97.99	69.99

COD orders not available. Normally, orders are shipped within 24 hours. For same-day guaranteed rush shipment, add a \$10 special handling charge. For Saturday Delivery, add \$13 to Next Day rates. For hazardous products, add \$25.00 for each individual item to above rates. Hazardous items must ship ground. Collect shipments are subject to a \$7.50 processing fee per order and can only be processed when items are shipped direct from AliMed's warehouse. Any items shipped from an alternate site will have charges prepaid and added to the invoice as per our shipping table.

Truck Shipments

(In addition to any applicable Ground Shipments.)

Merchandise value	Shipping & handling Continental US only
\$ 00-299.99	\$145.00
300-499.99	195.00
500-999.99	225.00
1000-1999.99	280.00
2000-3599.99	380.00
3600+	10% of order value

COD orders not available. Products are delivered to your receiving dock. Add \$75 for inside delivery and \$125 for Liftgates to above rates.

*Shipments to HI, AK, and Puerto Rico are sent by 2nd Day Air unless otherwise requested. Shipments to Canada may incur additional charges, including, without limitation, duties, taxes, or country-specific fees.

**Shipping Rate Table does not apply to international shipments. International shipments (except Canada) require a freight quote. Contact our International Department for assistance.

Everything you need to make your facility bariatric-ready!


For our full list of over 1,000 bariatric products, visit www.AliMed.com



297 High Street, Dedham, MA 02026

OR CURRENT OCCUPANT

PRESORTED
STANDARD
U.S. POSTAGE
PAID
AliMed® inc.



©2011 AliMed® inc. All rights reserved. Printed in U.S.A. Specifications, prices, and availability subject to change without notice. \$1.50

Supports 1,200 lbs.

Bariatric Bed Pan

- ✓ Low profile ($\frac{1}{2}$ "H in back)
- ✓ Easier to insert for caregivers
- ✓ Better patient comfort
- ✓ Low-splash guard

p5



Instantly recognize the weight capacity of your equipment

- ✓ Satisfies SRC requirement for certification
- ✓ Adheres to wood, metal, plastic
- ✓ Preserves patient dignity

p11



$\frac{1}{2}$ " x $\frac{1}{2}$ "
Shown actual size



6617

ORDER NOW! [CLICK www.AliMed.com](http://www.AliMed.com) CALL 800-225-2610 FAX 800-437-2966