



How to go lift-free... and everything you need to do it!

INSIDE!
*Hundreds of products
to increase safety
and reduce costs*



See how to save money with safe patient handling p2



**Sliding
Transfer Boards**
p18



**Gait Belts and
Handles**
p22



**Lateral
Transfer**
p12



Lifters
p4



Repositioning
p11



Slings
p8

No-Lift makes \$ense

What safer patient handling means for your budget

Better equipping your facility for safer patient handling has quantifiable return on investment (ROI).

It reduces staff injury potential and resulting workers' compensation, enhances patient safety and comfort, and saves money for the facility. Simply, it makes good business sense.

Below are results* of longitudinal research based on before and after implementation of safe patient handling protocols within nursing home facilities:

- Up to 60% reduction in work injury claims
- Up to 66% reduction in lost days
- Up to 51% reduction in workers' compensation costs
- Dramatic reductions in staff turnover
- 18 to 1 paybacks over five years while increasing patient safety and comfort.

*Source: (NIOSH longitudinal study of nursing home back strain injuries, before and after a safe mechanical lifting program, involving six facilities, six years, and 2,646 nursing staff.)

**For every \$1 spent,
...you save up to \$18*
That's about \$3,900
...per bed year**



Over 5 years

- Invest: \$10,880
- Save: \$195,840
- Financial gain: \$184,960

Legislation

More states now require "No-Lift" patient handling. Avoid citations due to improper equipment and patient handling devices.

Safer and more comfortable

With the implementation of better, safer patient handling protocols, your facility will see safer transfers, fewer slips and falls, and better overall patient comfort.

Think about the savings!

For every \$1 spent on safe patient handling products, you can save up to \$18. That means a savings of up to \$3,900 per bed year. In five years if you've invested \$10,880. That's like saving \$195,840, a financial gain of \$184,960!

Want to save like this? Let AliMed show you how!

Download our free safe patient handling white paper and ROI calculator at <http://safetycalculator.AliMed.info> and understand the financial benefits of practicing good healthcare ergonomics.



The inexpensive sit-to-stand lift

Stella Stand Assist Lift Quick and easy access to cleaning and changing. Less time required for lift and transfer. Secure and safe for both patient and caregiver. Doubles as a commode lift for additional flexibility. Assists with supported-walking and weight-bearing exercises. Ergonomic design and easy maneuverability. Height-adjustable leg and knee support. Detachable foot platform. Foot pedal base opening.

40" L x 68" H • Spreader bar: 33"-68" H • Base closed: 26½" W (external); 21¼" W (internal) • Base open: 37" W (external); 32½" W (internal) • Overall base: 4½" H • Front caster w/brake: 6" diam. • Back caster w/brake: 6" diam. • Lifting mechanism: 24 VDC Actuator • Battery: AC 110V/220V to DC 24V 1 Amp • **Stella:** 117 lbs. • 400-lb. capacity • **Stella Bariatric:** 122 lbs. • 600-lb. capacity

#CXX77747 Stella Lift\$1,868.25 ea

#CXX79852 Stella Bariatric Lift\$3,319.25 ea

Suggested code: E0635

CAPACITY CAPACITY
UP TO
**600
lbs.**
CAPACITY CAPACITY



Compare and Save!

Stand Assist Slings				
Description	Size	Capacity	Item #	Price \$
Deluxe, Padded Sling	S/M/L*	450 lbs.	#CXX72081	197.00 ^{ea}
Stand Assist Buttock Strap	Adjustable	450 lbs.	#CXX77769	66.75 ^{ea}
Deluxe Padded Sling	Bariatric	600 lbs.	#CXX79853	203.25 ^{ea}
Stand Assist Buttock Strap	Bariatric	600 lbs.	#CXX79854	98.25 ^{ea}

*Specify size.



Replacement Standing Sling

Fleece cover guards against skin shear. Sturdy metal clips easily attach for sitting or standing and offer added lifting leverage and protection. Safety belt secures patient with locking belt buckle. Washable. Latex-free. Compatible with Arjo's Sarita™ lift.

#CXX77666\$197.00 ea



Spryte helps keep the resident active and engaged in the lifting process, with minimal caregiver assistance. Padded split seat swings out to allow loading or unloading and swings back to form a comfortable and secure seat for transport. Dual knee pads.

Max height: 43" • Max width: 25½" • Overall: 32" L
• Base: 4½" H • 6" casters • 60 lbs. • 400-lb. capacity

#CXX71 0090\$1,343.50 ea.....\$1,025.64 ea (2+)
Suggested code: E0635

Note: Spryte is for residents who have the strength and stability to lift and support themselves. For residents with limited strength and stability, an electric stand assist is recommended.



NEW!

CAPACITY CAPACITY
**400
lbs.**
CAPACITY CAPACITY



Convenient stand and transport aid

Genesis 400H Patient Lift Hydraulic pump. Easy assembly; requires no tools. Six-point sling bar for safe and comfortable lifting. Foot pedal opens base for chair/commode access. Fits into car trunk when disabled.

41"L x 46"H • Spreader bar: 21"-66"H • Base closed 26"W (external); 21"W (internal) • Base open: 35½"W (external); 30½"W (internal) • Overall base: 4½"H • Front caster: 3" diam. • Back caster w/brake: 4" diam. • 75 lbs. • 400-lb. capacity

#CXX79466\$784.00 ea
Suggested code: E0630

Slings sold separately on pp. 8, 9.



Hydraulic Standard Patient Lift Safely raises or lowers individual from any stationary position. High-performance hydraulics with easy-to-operate caster brakes. Adjustable-width base. Six-point cradle design.

41"L x 46"H • Spreader bar: 21"-72"H • Base closed: 23"W • Base open: 43"W • Overall base: 5½"H • Front and rear casters: 4" diam. • 87 lbs. • 400-lb. capacity

#CXX13010 Manual Patient Lift.....\$ 701.00 ea
#CXX70137 Digital Scale.....1,030.00 ea
Suggested code: E0630

Slings sold separately on pp. 8, 9.




Invacare 9805 Hydraulic Lift for safe and affordable transfers. Durable steel construction. Includes a padded swivel bar and push handle. Six-point swivel bar easily adapts to all sling styles. Elongated handle design serves a broad range of operator heights and offers improved handling ability.

46¾"L x 73"H • Spreader bar: 20"-64"H • Base closed: 22"W • Base open: 42½"W • Overall base: 6½"H • Front and rear casters (no brakes): 5" diam. • 68 lbs. • 450-lb. capacity

#CXX77531.....\$896.50 ea
Suggested code: E0630


Slings sold separately on pp. 8, 9.

Easily know weight capacity!

 **Weight Capacity Labels** for quick recognition of the weight capacity of equipment.

These labels help facilities meet the criteria for SRC Centers of Excellence for Bariatric Surgery.* These 1½" x 1½" vinyl labels are printed in large type that can be easily seen from 6-8 feet away. Patient dignity is improved by subtly stating the weight capacity; 400-lb. capacity label shows an image of a bird and the number 40. Labels safely adhere to flat or curved surfaces on: **Wood:** Chairs, crutches, rehabilitation treatment tables, etc. **Metal:** Wheelchairs, walkers, canes, patient lifts, etc. **Plastic:** Patient shifters, transfer boards, etc. **Vinyl:** Wheelchairs, geri-chairs, etc. **Can be sprayed and wiped with disinfectant.**

*Surgical Review Corporation (SRC) requires that facilities have equipment available for morbidly obese patients and that all staff members can easily recognize the weight capacity for each piece of equipment. See: <http://www.surgicalreview.org/SAG/v4eBSCOE.aspx#Equipment>



Adheres well to curves!

#CXX71 1511 Variety Pack Labels: 175/pk: (25 labels each of 30, 40, 50, 60, 70, 80, 100)\$70.00 pk....\$61.00 pk (3+)

450E Electric Patient Lift Suited for nursing homes and in-home care. Six-point spreader bar offers additional safety and comfort. Unique foot pedal base opening for wheelchair or commode. Low base height for maximum stability and ability to fit under lower beds. Compact lift has backup manual crank to allow continuous use during emergencies.

41" L x 48" H • Spreader bar: 15½"-64" H • Base closed: 26½"W (external); 21¾"W (internal) • Base open: 38"W (external); 33½"W (internal) • Overall base: 4½" H • Front caster w/brake: 3" diam. • Back caster w/brake: 4" diam. • Lifting mechanism: 24 VDC Linear Actuator • Battery: 12V 4½ AH x 2 • Battery Charger: AC 110/220V to DC 24V 1 AMP • 114½ lbs. • 400-lb. capacity

#CXX77745\$1,468.75 ea

Suggested code: E0635

Slings sold separately on pp. 8, 9.



Easy to maneuver

400 lbs.

Electric Patient Lift Jumbo actuator pump with emergency button can switch to manual mode. Easy charge in AC outlet. Audible warning when battery power falls below 50%. Foldable mast and base for easy transportation. Holds 4- and 2-point slings.

41" L x 76" H • Spreader bar: 27½"-76" H • Base closed: 15"W • Base open: 41"W • Overall base: 2"H • Front/Rear caster: 4" diam. • Lifting mechanism: 24 VDC Linear Actuator • Battery: 12V 4½" AH x 2 • Battery charger: AC 110V/220V to DC 24V 1 AMP • 106 lbs. • 450-lb. capacity

Electric Bariatric Lift: 41" L x 76" H • Spreader bar: 32"-74½" H • Base closed: 25¾"W • Base open: 45¾"W • Overall base: 2"H • Front casters: 3" diam. • Rear casters: 4" diam. • Lifting mechanism: 24 VDC Linear Actuator • Battery: 12V 4½" AH x 2 • Battery charger: AC 110V/220V to DC 24V 1 AMP • 127 lbs. • 600-lb. capacity

#CXX13240 Electric Lift.....\$1,855.00 ea

#CXX70030 Bariatric Lift.....2,968.00 ea

#CXX70137 Digital Scale.....1,030.00 ea

Suggested code: E0635

Slings sold separately on pp. 8, 9.



drive

600 lbs.



#CXX70137

Invacare Reliant Battery-Powered Lift

Removable battery easily charges overnight on wall-mounted charging station. Equipped with a quiet long-life motor. Wide lift range for easier transfer. Adjustable wide base for maximum stability. Six-point sling.

50" L x 48" H • Spreader bar: 26"-76" H • Base closed: 26½"W • Base open: 41"W • Overall base: 6¾"H • Front and Rear casters: 5" diam. • Lifting mechanism: 24 VDC Linear Actuator • Battery: 12V 4½ AH x 2 • Battery charger: Output 24V DC • 106 lbs. • 450-lb. capacity

#CXX77532\$2,231.75 ea

Suggested code: E0630

Slings sold separately on pp. 8, 9.



450 lbs.

Heavy-Duty Electric Patient Lift

Emergency manual lowering device, universal battery, lift-capacity warning beep, six-point spreader bar, and foot-pedal base opening for wheelchair/commode. Six-point multi-tilt spreader bar.

50" L x 48" H • Spreader bar: 22"-71" H • Base closed: 26"W (external); 21½"W (internal) • Base open: 43"W (external); 38"W (internal) • Overall base: 6"H • Front and rear casters: 4" diam. • Lifting mechanism: 24 VDC Linear Actuator • Battery: 12V 4½ AH x 2 • Battery charger: AC 110V/220V to DC 24V 1 AMP • 137 lbs. • 600-lb. capacity

#CXX77820\$2,739.75 ea

Suggested code: E0630

Slings sold separately on pp. 8, 9.



600 lbs.

Molift patient lifts and slings

Molift Smart Portable lift that folds to fit into a closet or car trunk, but has big lift capabilities. Versatile locations include: facilities, home care, home, schools, and trips. Truly a versatile and compact lift—lightweight, portable, and easy to carry and move. Folds without tools and is easy to store. Uses standard 4-point lift slings. Electric with two lifting speeds, soft start/stop and standard battery power pack. Hard-Sided Transportation Case sold separately below.



Fits into a closet or car trunk!



26"W x 45¾"L x 70"H • Spreader bar: 19¼" - 65"H • Base external: 26"W • Base internal: 21½"W • Overall base: 4"H • Front caster: 3" diam. • Folded: 26"W x 45¾"L x 14"H • 55 lbs. (biggest component only 33 lbs.) • 330-lb. capacity

- #CXX71 1395 Molift Smart\$3,510.00 ea
- #CXX71 1398 Hard-Sided Transportation Case, 20"W x 50"L x 15"H.....365.75 ea

Molift Quick Raiser is ergonomically correct for sit-to-stand lifting and moving patients who can provide some assistance. It uses the natural movement pattern of bending forward when standing.

Its unique angled column raises the user in an up-and-forward movement that imitates the natural movement pattern to make the lift easier on the patient. Its extra low base and unique wheel pattern make it highly maneuverable and easy to position. The Quick Raiser is available in two versions.

Quick Raiser 1 is the standard version. **Quick Raiser 2+** includes electrical adjustable legs so the patient can be closer to the lift and self-adjusting knee pads. Has redesigned lifting arms to increase lifting capacity and provide more space for the user.



CAPACITY CAPACITY
UP TO
350
lbs.
CAPACITY CAPACITY

Quick Raiser 1: 13¼"W x 42"L x 68¾"H • Spreader bar: 39¾"-68¾"H • Base external: 13¼"W • Base internal: 10¼"W • Overall base: 2½"H • Front caster: 2" diam. • Side caster: 4" diam. • Back caster: 4½" diam. • 60 lbs. • 350-lb. capacity

Quick Raiser 2+: 14¼"W x 50"L x 50½"H • Spreader bar: 39"-68"H • Base closed external: 14¼"W • Base closed internal: 9½"W • Base open external: 47¼"W • Base open internal: 42½"W • Overall base: 4"H • Front caster: 3" diam • Back caster: 4½" diam • 81½ lbs. • 440-lb. capacity

- #CXX71 1399 Quick Raiser 1\$3,661.75 ea
- #CXX71 1401 Quick Raiser 2+5,083.50 ea



CAPACITY CAPACITY
UP TO
440
lbs.
CAPACITY CAPACITY

Battery accessories

Molift is the first manufacturer in the industry that equips all patient lifters with environmentally friendly batteries. The "green" batteries are based on the Nickel Metal Hydride System (NiMH) and contain no toxic cadmium mercury and lead, as do traditional batteries in the industry. They are particularly advantageous because of their low impact on the environment, and because of their strong and long-lasting performance. Their ability to manage many lifts before recharging and the ability to recharge hundreds of times can save a lot of money and resources.


- #CXX71 1396 Battery, all Molifts\$266.75 ea
- #CXX71 1397 Battery Charger.....233.50 ea



Molift Sling Guide

Choosing a sling depends on both the weight, height, and body type of the patient:

- A. Shoulder Width**
- B. Hip Width**
- C. Back Height**



Molift Sling Specifications

Size	Color Code	Max Width (Shoulders) A	Max Width (Hips) B	Height (Back) C	Recommended Weight	Safe Working Load
S	Light Blue	20"	30"	23½"	55-110 lbs.	350 lbs.
M	Red	22"	33¾"	24½"	98-200 lbs.	350 lbs.
L	Yellow	26"	40"	25"	200-350 lbs.	350 lbs.
XL	Green	30"	45½"	26"	350-525 lbs.	650 lbs.
XXL	Blue	33¾"	51½"	26"	500-650 lbs.	650 lbs.



1. Easy Sling

#CXX71 1404	S	\$206.75 ea
#CXX71 1405	M	206.75 ea
#CXX71 1406	L	190.75 ea
#CXX71 1407	XL	405.25 ea
#CXX71 1408	XXL	405.25 ea

2. Easy Sling w/Head Support

#CXX71 1409*	S, M, L	\$259.75 ea
#CXX71 2025	XL	325.00 ea
#CXX71 2026	XXL	489.00 ea

*Specify Size.



3. Easy Raiser Sling

#CXX71 1422	S	\$163.75 ea
#CXX71 1423	M	163.75 ea
#CXX71 1424	L	179.00 ea
#CXX71 1425	XL	255.75 ea

4. Easy Raiser Sling w/Support

#CXX71 1427	S	\$248.75 ea
#CXX71 1428	M	248.75 ea
#CXX71 1429	L	276.75 ea
#CXX71 1430	XL	490.00 ea

Shop **AliMed.com**

- :: See our newest products
- :: Up-to-date pricing
- :: Real-time inventory
- :: Order tracking

40% savings
everyday of the year!
Fit Hoyer®, Invacare®, Apex,
and most other lifters



Full-Body, Solid Seat



Commode Opening



Split Leg Mesh Seat Sling



Full-Body, Mesh Seat



Universal Padded Sling, Bariatric

Compare and \$ave!

FULL-BODY SLINGS				
Description	Size	Capacity	Item #	Price \$
Mesh Seat	M	400 lbs.	#CXX78250	103.25 ^{ea}
Mesh Seat	L	400 lbs.	#CXX78251	103.25 ^{ea}
Mesh Seat, Commode	M	400 lbs.	#CXX78254	110.25 ^{ea}
Mesh Seat, Commode	L	400 lbs.	#CXX78255	110.25 ^{ea}
Mesh Seat, Commode	XL	400 lbs.	#CXX78256	116.75 ^{ea}
Solid Seat	M	400 lbs.	#CXX78252	102.25 ^{ea}
Solid Seat	L	400 lbs.	#CXX78253	106.25 ^{ea}
Solid Seat	XL	400 lbs.	#CXX78257	113.50 ^{ea}
Mesh Seat	Bariatric	600 lbs.	#CXX78258	119.75 ^{ea}
Mesh Seat, Commode	Bariatric	600 lbs.	#CXX78259	127.00 ^{ea}
UNIVERSAL SPLIT-LEG SLINGS				
Description	Size	Capacity	Item #	Price \$
Mesh Seat	S/M*	400 lbs.	#CXX72077	130.00 ^{ea}
Mesh Seat	L	400 lbs.	#CXX72078	143.50 ^{ea}
Padded Seat	S/M/L*	400 lbs.	#CXX72079	179.25 ^{ea}
Padded Seat	Bariatric	600 lbs.	#CXX72080	191.75 ^{ea}
STAND-ASSIST SLINGS				
Description	Size	Capacity	Item #	Price \$
Deluxe, Padded Seat	S/M/L*	450 lbs.	#CXX72081	197.00 ^{ea}
Stand-Assist Buttocks Strap	—	450 lbs.	#CXX77769	66.75 ^{ea}
Stand-Assist Sling	Bariatric	600 lbs.	#CXX79853	203.25 ^{ea}
Stand-Assist Buttocks Strap	Bariatric	—	#CXX79854	98.25 ^{ea}

*Specify size.

For most 4- and 2-hook lifts including Apex, Invacare, Hoyer, EZWay, and Vancare.

Sling Specifications

	M	L	XL	2XL
Width	25"	30"	30"	45"
Length	50"	55"	57"	57"
Cap. (lbs.)	450	450	450	600

Arjo replacement slings



Replacement Slings for Arjo lifts are made to the same specifications and quality in the USA as imported Arjo-brand slings. Include rear pocket to help caregiver push sling under patient. Heavy-duty nylon, latex-free, washable.

#CXX71 2216\$195.00 ea.....\$182.75 ea (10+)
Specify S, M, L, XL.

Call! **VOLUME PRICING**

Best for infection control—single-patient-use slings

Can be allocated to a patient and cost only a fraction of a nylon re-use sling. They can be billed to a patient and become revenue instead of a cost item. They reduce cross-contamination and are disposable—no laundry costs, and no waiting for the laundry that never returns. May be used with most 4- or 6-hook patient lifts. 600-lb. capacity. Individually polywrapped. 10/cs.



SAFE PATIENT HANDLING

Patient Slings

Single-Patient-Use Slings*


Size	Wt. Range	Case #	Price \$	(3+) \$
W/O HEAD SUPPORT				
M	140-200 lbs.	#CXX70130	297.25 ^{cs}	281.35 ^{cs} (28.14 ^{ea})
L	200-400 lbs.	#CXX70131	317.75 ^{cs}	299.69 ^{cs} (29.97 ^{ea})
W/HEAD SUPPORT				
S	85-140 lbs.	#CXX70132	291.50 ^{cs}	276.05 ^{cs} (27.60 ^{ea})
M	140-200 lbs.	#CXX70133	333.00 ^{cs}	314.72 ^{cs} (31.47 ^{ea})
L	200-400 lbs.	#CXX70134	343.00 ^{cs}	324.92 ^{cs} (32.49 ^{ea})
Bariatric	400-600 lbs.	#CXX71 0088	379.00 ^{cs}	335.44 ^{cs} (33.54 ^{ea})
FOR STAND-ASSIST LIFTS #CXX77747 AND #CXX79852, 10/PK				
M	140-200 lbs.	#CXX71 1172	525.25 ^{pk}	492.56 ^{pk} (49.26 ^{ea})
L	200-400 lbs.	#CXX71 1173	561.25 ^{pk}	528.62 ^{pk} (52.86 ^{ea})

* Fit most Invacare, Apex, Hoyer, EZWay, Vancare and other 4-hook lifts.

Bestcare Apex® Dynamics Slings

Lifts, pp. 4, 5

Description	Item #	Size	Capacity	Regular \$	Special AliMed® \$
Universal Sling, Mesh	#CXX77753	M	400 lbs.	145.00 ^{ea}	116.75^{ea}
	#CXX77754	L	400 lbs.	150.00 ^{ea}	129.75^{ea}
Universal Sling, Padded	#CXX77756	M	400 lbs.	250.00 ^{ea}	164.50^{ea}
	#CXX77757	L	400 lbs.	250.00 ^{ea}	184.50^{ea}
	#CXX77758	XL	600 lbs.	300.00 ^{ea}	190.75^{ea}
Hammock Sling, Padded	#CXX77762	M	400 lbs.	280.00 ^{ea}	178.25^{ea}
	#CXX77763	L	400 lbs.	280.00 ^{ea}	190.00^{ea}
Bariatric Stand-Assist Sling	#CXX79853		600 lbs.	380.00 ^{ea}	203.25^{ea}
Stand-Assist Buttocks Strap	#CXX77769		400 lbs.	100.00 ^{ea}	66.75^{ea}
Bariatric Buttocks Strap	#CXX79854		600 lbs.	130.00 ^{ea}	98.25^{ea}
Toileting Sling Padded	#CXX77771	M	400 lbs.	360.00 ^{ea}	187.00^{ea}
	#CXX77772	L	400 lbs.	360.00 ^{ea}	197.25^{ea}




Invacare® Slings

Lifts, pp. 4, 5

Fully padded polyester/nylon. Can be used in wet or dry environments.

Sling Specifications	Full-Body			Full-Body w/Commode			Divided-Leg w/Headrest		
	M	L	XL	M	L	XL	M	L	XL
Size									
W	41½"	45½"	45½"	41½"	45½"	45½"	37½"	41"	44½"
L	54⅞"	60½"	65⅜"	54⅞"	60½"	65⅜"	62¼"	67⅞"	72⅜"
Commode Open.	—	—	—	8"x11"	8"x11"	9"x13"	—	—	—
Capacity (lbs.)	190	285	400	190	285	400	190	285	400

Description	M	L	Price \$	XL	Price \$
Mesh Full-Body	#CXX77521	#CXX77524	161.25 ^{ea}	—	—
Solid Full-Body	—	—	—	#CXX77527	180.00 ^{ea}
Mesh Full-Body w/Commode Opening	#CXX77522	#CXX77525	173.75 ^{ea}	#CXX77528	180.00 ^{ea}
Solid Divided-Leg w/Headrest	#CXX77523	#CXX77526	211.25 ^{ea}	#CXX77529	217.50 ^{ea}



Disposable transfer aids

Single-Patient-Use Slings cost only a fraction of nylon reusable slings, save money, and help reduce cross-contamination. Simply allocate to one patient for the entire visit to reduce spreading germs, and bill the patient rather than eat the cost. Plus, save with one less thing for the laundry. May be used with most 4- or 6-hook patient lifts. 600-lb. capacity. Individually polywrapped. 10/cs.



Single-Patient-Use Slings*

Size	Wt. Range	Case #	Price \$	(3+) \$
W/O HEAD SUPPORT				
M	140-200 lbs.	#CXX70130	297.25 ^{cs}	281.35 ^{cs} (28.13 ^{ea})
L	200-400 lbs.	#CXX70131	317.75 ^{cs}	299.69 ^{cs} (29.97 ^{ea})
W/HEAD SUPPORT				
S	85-140 lbs.	#CXX70132	291.50 ^{cs}	276.05 ^{cs} (27.61 ^{ea})
M	140-200 lbs.	#CXX70133	333.00 ^{cs}	314.72 ^{cs} (31.47 ^{ea})
L	200-400 lbs.	#CXX70134	343.00 ^{cs}	324.92 ^{cs} (32.50 ^{ea})
Bariatric	400-600 lbs.	#CXX71 0088	379.00 ^{cs}	335.44 ^{cs} (33.54 ^{ea})

* For Apex, Invacare, Hoyer, EZway, Vancare and most 4- or 6-hook lifts



Single-Patient-Use Gait Belts Disposable cotton belt you can allocate to the patient throughout their stay and for them to take home after their discharge. 36/cs.

Metal Buckle

#CXX71 0126	54"	\$242.00 cs194.00 cs (5+) (5.39 ^{ea})
#CXX71 0132	70"	242.00 cs194.00 cs (5+) (5.39 ^{ea})

Plastic Buckle

#CXX71 0129	54"	\$242.00 cs194.00 cs (5+) (5.39 ^{ea})
#CXX71 0135	70"	242.00 cs194.00 cs (5+) (5.39 ^{ea})

Cloth transfer devices

SafetySure™ Mobility Sheet For easier boosts, lifts, and transfers. Plastic insert gives some rigidity for sliding under patient on flat gurney. Model without plastic insert recommended when head is raised.

S: 16"W x 39"L; 4 handles; 1 lb.; 250-lb. capacity • **M:** 36"W x 39"L; 4 handles; 2 lbs.; 350-lb. capacity • **L:** 34"W x 60"L; 8 handles; 4 lbs.; 500-lb. capacity • Two-year warranty



With Plastic Insert

#CXX78592	S	\$121.00 ea
#CXX78593	L	261.00 ea

Without Plastic Insert

#CXX78594	S	\$ 86.50 ea
#CXX78595	M	106.00 ea
#CXX78596	L	127.00 ea



SafetySure™ Super Slide Sliding board with handles and an upholstered mobility sheet attached underneath. Roll patient onto side to place device. Grasp handles and easily pull patient to targeted surface. 250-lb. capacity.

#CXX70515\$325.00 ea.....\$304.00 ea (6+)





Sliderz Disposable Transfer Sheets

Easy-to-use, friction-reducing sheets minimize the force required by caregivers to reposition patients. Dispense from a continuous roll with perforated edges for optimum use. Eliminate heavy manual lifting while minimizing costs and infection risk found with shared transfer devices. For single-patient use; may be used several times before discarding. One roll of 50 perforated sheets. Latex-free. Made with recycled material. Multi-directional sliding capability. Do not wash. Discard if soiled or damaged.

#CXX71 1175\$86.25 pk.....\$78.25 pk (3+)

Single patient use—help stop the spread of infection

Single Patient Use Shifter/Booster

Hailed by clinicians as an economical transfer aid that helps reduce cross-contamination, this boosting aid is issued to one patient and should remain with that patient regardless of what department he/she has been transported to—making it always available for repositioning or transfers when needed. Works great on Geri-chairs too. Folds for storage. Dispose of when patient is discharged, or if wear is visible. Nonwoven fabric with "Slide Glide" coated bottom surface. Five reinforced handholds on each side. 33 $\frac{3}{4}$ "W x 77"L. Individually polybagged. 10/cs.



"Slide Glide" coating on bottom minimizes friction.



Works great on Geri-chairs too.

Store Booster in mesh hanging bag.




#CXX10802 Shifter/Booster.....\$47.75 cs
.....\$39.50 cs (10+) (3.95 ea)
#CXX71444 Mesh Hanging Bag.....7.75 ea
.....6.25 ea (6+)

In-Bed Patient Positioning System

Includes pad and two 34" wedges. Bolsters support sidelying resident. Hook-and-loop straps hold bolsters in place. Two strong handles for easier turning. Launderable.



#CXX55 5036\$118.00 set



Bariatric Turner allows up to six attendants to safely transfer, turn, reposition, or lift bariatric patients. 12 handles. Wipe clean or launder (temperature of water not to exceed 180°). 30"W x 72"L. 400-lb. capacity.

#CXX77806\$129.00 ea



Half-length and full-length shifters

Comfortably transfer patients in any position—even those with tubes and lines

Half Shifter helps reduce injury risk during transfers. Transfer patients in any position—upright, lying, or those with epidural catheters or spinal anesthesia lines.

Six built-in handholds also function as hanging holes. Bridges small gaps for uneven transfers. Radiolucent. Warp-resistant for a smooth sliding surface. Store a Half Shifter in every room for staff safety and convenience. **[OS]**


23"W x 31"H x 0.08" thick • 1.7 lbs. • 250-lb. capacity

#CXX93 1827\$114.00 ea\$93.50 ea (10+)



Bariatric Shifter Lateral transfers for patients up to 600 lbs. Reduce caretaker back strain by minimizing friction on transfer surfaces. Elongated shape handholds are ergonomically correct. Flexible plastic is great for uneven surfaces. 28"W x 72"L x $\frac{3}{16}$ "H. 600 lb.-capacity. Purple. **[OS]**

#CXX93 5944\$325.00 ea\$274.00 ea (3+)



Flexible half shifter!

Ergonomic Half Shifter helps reduce back strain by minimizing friction on transfer surfaces. Elongated shape brings handles to an easy-to-reach distance. Flexible plastic is great for uneven surfaces.

23"W x 33"L x $\frac{1}{2}$ " thick • 1.6 lbs. • 250-lb. capacity

#CXX92 7349\$47.50 ea
Specify Orange, Royal Blue, White, Gold, or Teal.



Flexible for easier transfers; easy-to-reach handles



Great for uneven surfaces

Save On Shipping

Shifter ships oversized with the same rate for 1–10 units. Buy up to 10 and only pay one shipping charge.

Store off the floor!


Wall Hanger Adjustable pegs accommodate most shifters. Solid oak, clear finish. Mounting hardware included. Vertical is 5"W x 18"L x $\frac{3}{4}$ " thick. Horizontal is $5\frac{1}{2}$ "W x 32"L x 1" thick.

#CXX9-704 Vertical\$71.50 ea

#CXX75134 Horizontal.....\$96.00 ea

Patient slides smoothly

Conductive can't shock!




Antistat Shifter slides easily on bed sheet or stretcher pad. Slippery coating protects patient from sticking to mover and shields against static buildup. Thin, warp-resistant, smooth polyethylene board. 10 handholds. Semirigid board—bridges uneven gaps.

Conductive Shifter Surprise “zaps” are painful and dangerous for patients and can jeopardize electronic equipment. With this shifter, static charges can’t build up. No more painful shocks, ever—guaranteed. Carbon-filled plastic actually conducts electricity (not a surface treatment). 10 built-in handles.

3/16" thick • 250-lb. capacity Standard Width: 22"W x 72"L for shifting patients from bed to stretcher to table • 9 lbs. Narrow Width: 18" x 72"L for close quarters. • 7 3/4 lbs. **OS**

Shifters				
Item #	Item #	Color	ea \$	(3+) \$
BASIC SHIFTER				
Standard Antistat	#CXX9-076	Aqua	235.00 ^{ea}	197.00 ^{ea}
Standard Antistat	#CXX9-719	White	238.75 ^{ea}	200.75 ^{ea}
Narrow Antistat	#CXX9-077	Aqua	189.00 ^{ea}	158.25 ^{ea}
CONDUCTIVE SHIFTER				
Standard	#CXX9-182	Black	325.00 ^{ea}	274.00 ^{ea}



Horizontal



Vertical

MoveMaster® Patient Shifter has a contoured shape for easier patient loading and unloading, especially in tight quarters. Built-in Agion® ionic silver antimicrobial agent naturally and continuously controls odor for the life of the board. Antistatic inhibitors, X-ray translucent. Anti-pinch handholds for improved, painless flat-surface transfers. Ergonomic design to help reduce healthcare workers’ stress and injuries. The MoveMaster is made from high-density, stress-relieved polyethylene. 22"W x 72" L x 3/16" thick. Shipping weight: 13 lbs. **OS**

- #CXX9-183 Transfer Board\$254.00 ea
- #CXX9-18305 Transfer Board, 5/pk608.00 pk

Transfer without lifting

- Reduces risk of back injuries
- Perfect for O.R. transfers
- Ball-bearing rollers
- Durable cover




Patient Roller Transfer Device lets one or two staff members move a patient from one surface to another with very little effort—no lifting or pulling—and without disturbing the sheets. Transfers the patient between slightly uneven surfaces (up to 3" height difference) without lifting. Long-lasting, trouble-free aluminum rollers with steel ball bearings and heavy-duty black conductive vinyl cover. Replacement covers available, see below.

440-lb. capacity • Cleans with soap and water or disinfectant • Aluminum rollers
 • **Short:** 14"W x 25"L • 5.8 lbs. **Medium:** 14"W x 30"L • 6.3 lbs. **Long:** 14"W x 67"L • 13 lbs.

Replacement Covers fit every roller on the market except those made by Spenco®. Make your transfer device look brand new with fresh covers! Made of black conductive vinyl.

Short: 14"W x 25" L • **Medium:** 14"W x 30" L

• **Long:** 14"W x 67" L



Patient Transfer Board allows transfer of heaviest patients without lifting. This roller-free transfer cover rotates around board during transfer. Long Board folds for easy storage. Allows transfer of patients into a Fowler position. Lightweight and radiolucent. Disposable covers for the Long Board help prevent cross-contamination.

500-lb. capacity • **Short:** 16"W x 39"L • 3 lbs. **Long:** 20"W x 71"L • 6.3 lbs.

#CXX93 3279 Long\$922.25 ea
 #CXX93 3280 Short614.00 ea
 #CXX93 3285 Disposable Covers for Long, 50/cs
128.00 cs

Transfer Device

Rack Easy to install and adjustable to fit any length roller. Your transfer device will last longer, stay cleaner, and take up less space when it's conveniently stored on the wall. Top bracket: 7" L x 1 5/8" thick. Bottom bracket 9 1/4" L x 1 3/4" thick.



Fits any manufacturer's model

#CXX9-721.....\$57.00 set.....\$52.00 set (3+)

Patient Roller Transfer Device **BUY A KIT AND SAVE \$35.25!**

Description	Roller		Roller/Rack Kit		Replacement Covers	
	Item #	Price \$	Item #	Price \$	Item #	Price \$
Short	#CXX9-254	225.00 ^{ea}	#CXX92 0889	246.00 ^{kit}	#CXX9-209	84.25 ^{ea}
Medium	#CXX9-729	225.00 ^{ea}	#CXX92 0890	246.00 ^{kit}	#CXX9-478	84.25 ^{ea}
Long	#CXX9-728	315.00 ^{ea}	#CXX92 0891	331.00 ^{kit}	#CXX9-477	125.00 ^{ea}

Air-assisted transfer mats

Slides on Air Inflated mat supports patient in a horizontal position allowing patient to be transferred or repositioned with little to no effort by the caregiver and without discomfort or trauma to the patient. Cushions patient during positioning and transfers. Extra-long handles offer proper ergonomics for safe patient handling. No weight limit—ideal for bariatric patients. Double layered top surface creates a smooth top surface for extra shear protection. Sani-liner included with transfer mat. Air supply pump sold separately. Specify Launderable or Wipe-clean Mat. Shipped freight.

70"L (Short model is 40"L) • Five-year warranty on Transfer Mat and Sani-liner • Air supply: 1,100 Watts • Grounded 15' cord

#CXX71 1780	26" W.....	\$4,465.00 ea
#CXX71 1781	28" W.....	4,465.00 ea
#CXX71 1782	32" W.....	4,465.00 ea
#CXX71 1783	39" W.....	4,675.00 ea
#CXX71 1784	50" W.....	4,680.00 ea
#CXX71 1785	32" W, Short	3,195.00 ea
Specify Launderable or Wipe-clean Mat.		
#CXX71 1786	Air Supply Pump	2,019.00 ea

Call! **VOLUME PRICING**



Transfer Mat

Air Supply Pump

Air lifting cushions

CAMEL and ELK are inflatable cushions that lift a fallen person from the floor. Easy to store. When inflated and placed underneath a supine person and inflated, the **CAMEL** gently lifts the person into a raised, seated position ready to stand. Comes with two AirFlo Compressors, slide sheet, transfer belt, carry bag, and battery charger. 27½"W x 44"L. Seat is 18½"D. 700-lb. capacity.

The **ELK** lifts the person up in a supine flat position. Inflated height is 22". Comes with an AirFlo Compressor, lifting cushion, and battery charger. 20"W x 20"L. 8 lbs. 1,000-lb. capacity.

#CXX71998	Camel.....	\$3,708.00 ea
#CXX71999	Elk	2,467.00 ea
#CXX72000	AirFlo Compressor.....	837.00 ea

Call! **VOLUME PRICING**



The CAMEL

CAPACITY CAPACITY
700 lbs.
CAPACITY CAPACITY



The ELK

CAPACITY CAPACITY
1000 lbs.
CAPACITY CAPACITY

Take control of bed repositioning

"No-Lift" TURNER™



**One person
can turn
patient!**

**Turn patient from
side to side
without lifting!**

Bed Repositioning

HOW THE TURNER WORKS...



"No-Lift" TURNER is positioned across bed to aid in turning schedules.



Roll resident onto side to position TURNER under patient's hips and shoulders.



Once TURNER is positioned, patient can be safely repositioned by one staff person.



For more information, and to see a video, go to

<http://video.AliMed.info>
and click on
"No Lift" BOOSTER and TURNER

"No-Lift" TURNER™ Simply position resident on top of the tube just as you would with a drawsheet. Grab reinforced edges of the top tube layer and pull it toward you. No lifting and less back strain! Once TURNER is put in position, which generally requires two people, only one staff member is needed to move most sidelying residents up to 250 lbs. Since the TURNER eliminates drawsheet drag and shear, the resident's skin integrity is never compromised.

Flat dimensions: 29"W x 52"L • Black • Machine-wash, air-dry • Disinfectant-safe

#CXX72124	TURNER w/Handygrips.....	\$ 83.00 ea\$ 72.75 ea (3+)
#CXX71670	10 TURNERs w/o Handles w/10 Free Hanging Bags592.00 cs (\$59.20 ^{set!})	
#CXX71671	25 TURNERs w/o Handles w/25 Free Hanging Bags1,292.00 cs (\$51.68 ^{set!})	
#CXX71 2177	Bariatric TURNER	105.00 cs	
#CXX71444	TURNER Hanging Bags.....	7.75 ea6.25 ea (6+)

Note: If patient has highly sensitive skin, TURNER should be removed between uses.

"No-Lift" BOOSTER™



Boosts patients up the bed without lifting

HOW THE BOOSTER WORKS...



Place "No-Lift" BOOSTER under resident as you would a drawsheet.




Roll resident onto side and place BOOSTER underneath.



Slide resident up toward head of the bed by grabbing top layer of fabric.



Roll resident onto side and slide BOOSTER out from underneath.



Store TURNER or BOOSTER in mesh hanging bag.

"No-Lift" BOOSTER™ Eases caregiver back strain by making the resident easier to move toward the headboard. The BOOSTER is a flat tube of fabric with a virtually frictionless interior surface. Simply position the resident on top of the tube just as you would with a draw sheet. Now grab the reinforced edges of the top tube layer and slide it upwards. No lifting is needed—less risk of back strain! The patient's skin integrity is never compromised since the BOOSTER eliminates draw sheet drag and shear. Remove it from underneath the resident after each use.

The Standard BOOSTER has handles and will fit most patients. The Bariatric BOOSTER, which comes with handles, is ideal for larger residents. Our research team of four has moved residents up to 600 lbs.

Standard: (flat dimensions), 29"W x 46"L • 250-lb. capacity; **Bariatric,** 46"W x 46"L • 600-lb. capacity • Black • Machine-washable, air-dry • Disinfectant-safe

#CXX72123	Standard	\$ 72.50 ea	\$62.50 ea (3+)
#CXX72 12312	Standard, 12/cs	546.00 cs	(\$45.50 ea!)
#CXX75582	Bariatric	105.00 ea	94.50 ea (3+)
#CXX75 58212	Bariatric, 12/cs	720.00 cs	(\$60.00 ea!)
#CXX71444	BOOSTER/TURNER Hanging Bag	7.75 ea	6.25 ea (6+)



For more information, and to see a video, go to

<http://video.AliMed.info>

and click on

"No Lift" BOOSTER and TURNER



A. Double Notched Transfer Board fits around the upright of the wheelchair arm support. 24" has a single handhold; 29" has two. Plywood. Stress-tested to 400 lbs.

#CXX8501 24"\$72.75 ea
#CXX8503 29"93.75 ea

B. Multi-Ply™ Flat Transfer Board
Polyclad, splinter-free birch plywood. Stress-tested to 400 lbs. 8"W x 24"L. 1/2" thick at center; 1/4" at both ends.

#CXX8500\$57.50 ea

C. Maple Transfer Board 8"W and 3/4" thick in the center, tapering to 1/2" at edges. Stress-tested to 300 lbs.

#CXX8502 Solid, 24"\$48.50 ea
#CXX8146 Slots, 24"58.50 ea
#CXX8047 Slots, 30"64.50 ea

D. Flat Plastic Transfer Board High-density, linear plastic. Satin-surface is sanded smooth. 1/2" thick. 9"W. Stress-tested to 400 lbs.

#CXX8510 22"\$38.00 ea
#CXX8512 29"49.25 ea



SafetySure™ Transfer Boards bridge the gap between two transfer points. Butterfly model features backside protection panels. Two-year warranty.

Transfer Board, 12 1/2"W, 250-lb. cap

#CXX78599 23"\$62.00 ea
#CXX78600 30"70.00 ea

Butterfly Board, 12 1/2"W, 250-lb. cap

#CXX78601 23"\$64.25 ea
#CXX78602 30"82.75 ea

Shop **AliMed.com**

- :: See our newest products
- :: Up-to-date pricing
- :: Real-time inventory
- :: Order tracking

Positioning devices



Swivel Cushion Swivel in and out of cars with this well-padded cushion. Acts like a Lazy Susan. For those with leg and/or trunk weakness.

#CXX79856\$29.75 ea\$28.65 ea (3+)

Sit-Rite™ One-way slide allows patient to boost, or be boosted, back in chair with minimal friction, helping to prevent patient and staff injuries. Improves seating posture by helping to prevent the patient from sliding forward

in chair, reducing the potential for ulcer-causing sacral sitting.

#CXX78529\$121.00 ea



Easier to slide back



Helps prevent sliding forward



#CXX77415



#CXX77414

BeasyTrans™ Transfer a patient without lifting. Transfer from a wheelchair to an automobile, bed, commode, shower, or bench. Created from a combination of super-strong, space-age polymers. Does not cause friction on the user's skin. Circular seat rotates 360° to easily turn the user to the exact angle for comfortable placement. Absolutely no lifting involved!
 40" L x 1" H • Seat: ½" H x 12" D • 6½ lbs. • 400-lb. capacity

Beasy II™ More compact than the BeasyTrans. Great for car transfers.

27½" L x 1" H • Seat: ½" H x 12" D • 4½ lbs. • 400-lb. capacity

BeasyGlyder™ Crescent-shaped to help maneuver around fixed wheelchair arms and wheel obstructions.

32" L x ¾" H • Seat: ¾" H x 10" D • 4¾ lbs. • 400-lb. capacity

#CXX77413 BeasyTrans.....\$274.00 ea

#CXX77414 Beasy II263.00 ea

#CXX77415 BeasyGlyder269.00 ea



Like a Lazy Susan for feet!



SafetySure™ Pivot Disc Makes pivot transfers safer and easier. Reduces backstrain on caregiver. Nonslip surface. Two-year warranty.

#CXX78597 13"\$ 96.50 ea

#CXX78598 15"111.00 ea

#CXX70479 18"141.00 ea


Lift-Easy Transfer Belt offers additional lift leverage when thigh straps are wrapped in place to make lifting easier. Less force is exerted by caregiver, resulting in a heightened sense of security for the patient. Two side loops facilitate use as a gait belt. Hook-and-loop enclosures.

Sizing: Measure waist in inches.

- #CXX74093 S, 24"-30" ...\$53.75 ea ...\$41.50 ea (3+)
 #CXX74094 M, 30"-44" ...\$53.75 ea\$41.50 ea (3+)
 #CXX74095 L, 44"-58"\$53.75 ea\$41.50 ea (3+)



Ambulation belt w/handles on clinician gives patient secure grip and less strain on clinician.



Easy-grip handles!

Ambulation Belt is wider than typical gait belts. Nylon web handles and metal buckle closure. 4"W. Fits 30"-50" waist.

#CXX77905\$45.25 ea\$40.00 ea (3+)



Ergonomic Ambulation Belt is the widest, most supportive gait belt offered. Cone-shaped, reinforced handles, and front hook-and-loop closure with two plastic buckles. 4¾"W.

#CXX77906 Fits waists 36"-46"
\$44.25 ea\$39.00 ea (3+)
 #CXX77907 Fits waists 44"-52"
\$45.25 ea\$40.00 ea (3+)

NEW!



Lift Vest Transfer device is beneficial to both the patient and caregiver. Patented design spreads the lifting pressure evenly over the whole torso, while the high-strength webbing harness secures the upper body and prevents injuries. Promotes safer lifting techniques for the caregiver and allows them to lift from any direction.

Sizing: XXS (Chest: 33"-39"/Waist: 29"-35"), XS (Chest: 35"-41"/Waist: 31"-38"), S (Chest: 37"-43"/Waist: 33"-41"), M (Chest: 39"-45"/Waist: 35"-44"), L (Chest: 41"-47"/Waist: 37"-47"), XL (Chest: 43"-49"/Waist: 39"-50"), XXL (Chest: 45"-51"/Waist: 41"-52").

#CXX71 1938 Jean Blue/Khaki Poncho Vest*\$264.00 ea
 #CXX71 1939 Bath Poncho\$264.00 ea
 #CXX71 1940 Jean Blue/Khaki Autism Poncho*\$279.00 ea
 Specify XXS, XS, S, M, L, XL, or XXL.
 *Specify color

Easier rising aids

Portable Easy Up™ for post-hip surgery or for those with arthritis of the hips and knees. Lightweight chair lift adds 3" or 4" of height to any chair. Positions hips and knees at 90°. Handy carry strap.

Extra-firm foam • 15"W x 15"L sitting surface
 • Zipper closure • Machine-washable, navy blue fabric • Weighs less than 1 lb.

#CXX1220 3"H\$55.00 ea
 #CXX1221 4"H\$61.50 ea



Portable easy-up raises seat height.

Makes sit to stand easier.

Comes with convenient carry strap.



Post-Hip Surgery Cushion High-density foam for comfort and pressure reduction. Wipe-clean vinyl cover. 6" thick. 18"W x 16"L.

#CXX70 7052\$43.00 ea

Transfer slings



SafetySure™ Transfer Sling assists in repositioning patients while easing caregiver back strain. Anti-slip. Two-year warranty.

#CXX78591\$44.25 ea



Patient Handler Strong, flexible polyurethane ensures safe and secure patient transfer that is less stressful on the attendant's back.

#CXX80588\$56.50 ea

 **Order online**
www.AliMed.com



Bariatric Gait Belt eases caretaker back strain when ambulating with bariatric patients. 96" wipeable gait belt easily wipes clean with any commonly used disinfectant. Reinforced plastic strap is easy to grasp and comfortable for the patient. Specialty buckle ensures belt stays attached. 96". Yellow. 1,000-lb. capacity.



#CXX71 1779 Gait Belt\$ 29.25 ea
#CXX71 177920 20/cs.....513.00 cs (**\$25.65** ea!)



#CXX7696
#CXX7695

VELCRO® brand Walker Belts with strong VELCRO® brand hook-and-loop closures.

Quick Grab Strap model is easier for the attendant to grab when the patient needs assistance. Blue. **Two-Attendant Strap** model has two sturdy handles sewn securely into walker belt so two attendants can walk large patients with confidence. Black.

Quick Grab Strap
#CXX7695 3" x 54"\$25.00 ea ...\$20.25 ea (6+)
Two-Attendant Strap
#CXX7696 4" x 50"\$30.75 ea ...\$24.50 ea (6+)
#CXX77804 4" x 64"\$33.25 ea



SafetySure™ Transfer Belt Multiple hand grips. Anti-slip interior keeps belt in place. Two-year warranty.

Sizing: S fits 23"-36" waist; M fits 32"-48" waist; and L fits 42"-60" waist.

#CXX78578 XS.....\$42.00 ea
#CXX78579 S.....44.25 ea
#CXX78580 M.....46.25 ea
#CXX78581 L.....48.25 ea
With Synthetic Sheepskin Lining
#CXX78582 S.....\$61.25 ea
#CXX78583 M.....63.25 ea
#CXX78584 L.....65.25 ea

Comfort Take-Along™ Cushion 1" of T-Foam™ for comfort and shock absorbency. 1" of high-density foam for support and resilience.

Knit cover • Carry strap • 15"W x 15"D x 2"H

#CXX1270\$70.25 ea



Uplift Seat Assist™ This comfortable, easy-to-use cushion has a spring to help the user stand up. Memory foam with water-resistant coating helps protect against pressure sores. Waterproof cover is machine-washable in cold water. Built-in handle. Available in Standard or Plus size. 17"W x 19"D; 2"H closed.

#CXX1069 Standard: 17 lbs. • 80-230 lb. capacity
#CXX1961 Plus: 21 lbs. • 195-350 lb. capacity
#CXX1069 Standard.....\$180.50 ea
#CXX1961 Plus.....180.50 ea



Economical gait belts

Gait Belts 2"W reinforced cotton belt. Used to assist with ambulation, as well as chair and bed transfers. Machine-washable. Choose between a nickel-plated, easy-release metal buckle or a sturdy plastic buckle. 3-year warranty.



Call! **VOLUME PRICING**



Pinstripe Gait Belts

METAL BUCKLE			Cases of 20		
54"	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX7414	12.50 ^{ea}	#CXX74 1420	187.00 ^{cs}	167.00 ^{cs}
70"	#CXX7083	12.50 ^{ea}	#CXX70 8320	187.00 ^{cs}	167.00 ^{cs}
PLASTIC BUCKLE			Cases of 20		
54"	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX70868	12.50 ^{ea}	#CXX70 86820	187.00 ^{cs}	167.00 ^{cs}
70"	#CXX70869	12.50 ^{ea}	#CXX70 86920	187.00 ^{cs}	167.00 ^{cs}
				Only \$9.35^{ea}!	Only \$8.35^{ea}!



Red/White/Blue Gait Belts

METAL BUCKLE			Cases of 20		
54"	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX77175	12.50 ^{ea}	#CXX77 17520	187.00 ^{cs}	167.00 ^{cs}
70"	#CXX77178	12.50 ^{ea}	#CXX77 17820	187.00 ^{cs}	167.00 ^{cs}
PLASTIC BUCKLE			Cases of 20		
54"	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX77176	12.50 ^{ea}	#CXX77 17620	187.00 ^{cs}	167.00 ^{cs}
70"	#CXX77177	12.50 ^{ea}	#CXX77 17720	187.00 ^{cs}	167.00 ^{cs}
				Only \$9.35^{ea}!	Only \$8.35^{ea}!



Pastel Gait Belts

METAL BUCKLE			Cases of 20		
54"	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX77171	12.50 ^{ea}	#CXX77 17120	187.00 ^{cs}	167.00 ^{cs}
70"	#CXX77172	12.50 ^{ea}	#CXX77 17220	187.00 ^{cs}	167.00 ^{cs}
PLASTIC BUCKLE			Cases of 20		
54"	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX77173	12.50 ^{ea}	#CXX77 17320	187.00 ^{cs}	167.00 ^{cs}
70"	#CXX77174	12.50 ^{ea}	#CXX77 17420	187.00 ^{cs}	167.00 ^{cs}
				Only \$9.35^{ea}!	Only \$8.35^{ea}!



- MAIL AliMed® inc.**
Order Department, 297 High St., Dedham, MA 02026
- PHONE 800-225-2610** 8:00 AM–7:00 PM EST, Mon.–Fri.
- 24-HOUR FAX 800-437-2966**
- International:** 781-329-2900 FAX 781-329-8392
- ONLINE CATALOG** 24-hr. web ordering www.AliMed.com
- CHARGE IT**
By phone, mail, or fax. Please see our new Shipping Rate Table at right for exact charges. Credit cards are only accepted at time of order entry.
- QUESTION OR RETURN?** Customer Service, x6901

MINIMUM BILLING Orders are subject to a \$25.00 minimum merchandise billing plus shipping, handling, and insurance. International orders (except Canada) require a \$300.00 (U.S.) minimum order.

TERMS Net 15 to established accounts. All prices F.O.B. shipping point. Shipment made at prevailing price at time of order. All prices in U.S. \$. Prices are professional wholesale. Patient may order directly; retail pricing applies. Prices, availability, and specs subject to change. Visit www.AliMed.com for the most current pricing. Items may be discontinued at any time. AliMed reserves the right to substitute similar products to those shown in catalogs at any time. Overdue accounts are subject to a late payment fee of 2.0% per month, resulting in an annual rate of 24%.

Quantity Discounts: Item numbers, sizes, colors, Lefts and Rights, etc. cannot be assorted for a quantity discount. One item #, one style per quantity discount order (unless otherwise specified).

SHIPPING We ship FedEx/UPS whenever possible. Charges are prepaid and added to invoice. Freight shipments are freight collect. See Shipping Rate Table for exact charges.

OS Oversized items subject to additional freight charges.

OVERSIZED CHARGES Apply to all orders including contract orders, items with free ground shipping, and free shipping promotions.

GLN #0733657000009

TAX Applies to all orders shipped in MA, SC, WA, and SD. If you are tax exempt, please forward a copy of your tax exempt or resale certificate. If your facility is not in MA, SC, WA, or SD, please submit use tax directly to your state.

INSURANCE Your shipment is 100% insured against loss or damage.

RETURNS Our return policy is consistent in maintaining the best value and highest service level for all of our customers. If the merchandise you received is defective or not as you ordered, please contact our customer service staff for quick resolution.

All other returns will be subject to a 20% restocking fee (\$10.00 minimum) except "Special Order" items that have been "drop shipped". (These items will be subject to a restocking fee that may vary at the discretion of the individual vendor.) The guidelines below must be followed in order to receive any credit:

- Merchandise credit only
- Product must be in new, resaleable, unmarked condition and in original packaging
- Merchandise must be returned within 30 days of invoice date
- Custom products, beds, books, videos, shoes, eating utensils, and toileting products are not returnable
- Open packages cannot be returned for credit

Infection control gait belts

Wipeable Gait Belts Smooth plastic strap easily wipes clean with any commonly used disinfectant. We recommend using PureGreen24™, a surface sanitizer that remains effective for up to 24 hours (see www.AliMed.com).



Wipeable Gait Belts

Cases of 20

YELLOW	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX71 1274	18.25 ^{ea}	#CXX71 127420	320.00 ^{cs} (16.00 ^{ea})	279.00 ^{cs} (13.95 ^{ea})
70"	#CXX71 1275	19.75 ^{ea}	#CXX71 127520	350.00 ^{cs} (17.50 ^{ea})	315.00 ^{cs} (15.75 ^{ea})
NEW! 96" (Bariatric)	#CXX71 1779	29.25 ^{ea}	#CXX71 177920	513.00 ^{cs} (25.65 ^{ea})	483.81 ^{cs} (24.19 ^{ea})
ORANGE					
54"	#CXX71 1020	18.25 ^{ea}	#CXX71 1022	340.75 ^{cs} (17.03 ^{ea})	329.93 ^{cs} (16.50 ^{ea})
70"	#CXX71 1023	19.75 ^{ea}	#CXX71 1025	730.75 ^{cs} (36.54 ^{ea})	709.87 ^{cs} (35.49 ^{ea})
BLUE					
54"	#CXX71 1026	18.25 ^{ea}	#CXX71 1028	324.75 ^{cs} (16.23 ^{ea})	314.44 ^{cs} (15.72 ^{ea})
70"	#CXX71 1029	19.75 ^{ea}	#CXX71 1031	355.50 ^{cs} (17.78 ^{ea})	345.20 ^{cs} (17.26 ^{ea})



Gait Belt w/Antimicrobial Agent

Cases of 20

METAL BUCKLE	ea #	ea \$	case #	1 case \$	5 cases \$
54"	#CXX70257	15.00 ^{ea}	#CXX70 25720	189.00 ^{cs} (9.45 ^{ea})	141.00 ^{cs} (7.05 ^{ea})
70"	#CXX70256	16.00 ^{ea}	#CXX70260	210.00 ^{cs} (10.50 ^{ea})	147.00 ^{cs} (7.35 ^{ea})
PLASTIC BUCKLE					
54"	#CXX72108	15.00 ^{ea}	#CXX72 10820	189.00 ^{cs} (9.45 ^{ea})	141.00 ^{cs} (7.05 ^{ea})
70"	#CXX72109	16.00 ^{ea}	#CXX72 10920	210.00 ^{cs} (10.50 ^{ea})	147.00 ^{cs} (7.35 ^{ea})



Single-Patient Use Gait Belts Disposable cotton belt you can allocate to the patient throughout their stay and for them to take home after their discharge.

Single Patient Use Gait Belts

Cases of 36

METAL BUCKLE	case #	1 case \$	5 cases \$
54"	#CXX71 0126	242.00 ^{cs} (6.72 ^{ea})	194.00 ^{cs} (5.39 ^{ea})
70"	#CXX71 0132	242.00 ^{cs} (6.72 ^{ea})	194.00 ^{cs} (5.39 ^{ea})
PLASTIC BUCKLE			
54"	#CXX71 0129	242.00 ^{cs} (6.72 ^{ea})	194.00 ^{cs} (5.39 ^{ea})
70"	#CXX71 0135	242.00 ^{cs} (6.72 ^{ea})	194.00 ^{cs} (5.39 ^{ea})



To ensure prompt return processing, please attach the original packing list, invoice or RA number to the return. Refunds will be for merchandise value only. Please allow three to four weeks for return processing.

TRUCK SHIPMENTS Delivery Inspection Required

IMPORTANT NOTICE: All truck shipments must be fully inspected by receiver before signing for the delivery to ensure there is no concealed damage. If you cannot inspect immediately, make a note "pending inspection" with your signature. Any damaged shipment received should be refused by the receiver. These guidelines must be strictly followed in order to receive any credit on Truck shipment orders. **Overnight Truck Shipments - AliMed does NOT provide overnight delivery on orders that require delivery by truck.**

DAMAGED GOODS/SHORT SHIPMENTS Must be reported within 10 days of receipt. Be sure to note any damage or poor box conditions when signing for packages. If shipment arrives damaged or with parts missing, contact the carrier immediately. Be sure to save product packaging. All shipping charges and risk of loss are FOB shipping point and are the responsibility of the purchaser.

Shipping Rate Table**

Includes shipping and handling charges • Temporary fuel surcharge may be added to our published rates • No small order charge • Instant freight quotations • No additional freight or backorders
Some large volume orders and/or oversized items may require additional freight charges or truck shipping.

Ground/Air (FedEx, UPS, Postal) Shipping, handling and insurance

Merchandise value	Inside continental US	Outside continental US* and Canada
\$ 00-99.99	\$11.99	\$26.99
100-199.99	16.99	35.99
200-499.99	19.99	42.99
500-799.99	28.99	49.99
800+	34.99	55.99

COD orders not available. For 2nd Day, Next Day Air, or Saturday Delivery, additional charges apply to above rates. Normally, orders are shipped within 24 hours. For same-day guarantee rush shipment, add a \$10 special handling charge. For hazardous products, add \$25 for each individual item to above rates. Hazardous items must ship ground. Collect shipments are subject to a \$7.50 processing fee per order and can only be processed when items are shipped direct from AliMed's warehouse. Any items shipped from an alternate site will have charges prepaid and added to the invoice as per our shipping table.

Truck Shipments (In addition to any applicable Ground Shipments.)

Merchandise value	Shipping and handling Continental US only
\$ 00-299.99	\$145.00
300-499.99	195.00
500-999.99	225.00
1000-1999.99	280.00
2000-3599.99	380.00
3600+	10% of order value

No COD orders available. Products are delivered to your receiving dock. Inside delivery and liftgates are additional. International shipping prices, including Canada, reflect the base shipping charge, not including possible customs duties, taxes, and other country specific fees.

* Shipments to HI, AK, and Puerto Rico are sent by 2nd Day Air, unless otherwise requested. Shipments to Canada may incur additional charges, including without limitation, duties, taxes or country specific items.

** Shipping rate table does not apply to International shipments. International shipments (except Canada) require a freight quote. Contact our International department for assistance.

AliMed's new website makes purchasing faster and easier with even more useful information and resources!



Quick access with QR codes!
Just point and scan for instant connection to:

<http://www.alimed.com/alimed/patient-safety-fall-management/>



297 High Street, Dedham, MA 02026 OR CURRENT OCCUPANT

PRESORTED
STANDARD
U.S. POSTAGE
PAID
AliMed® inc.



©2011 AliMed® inc. All rights reserved. Printed in U.S.A. Specifications, prices, and availability subject to change without notice.

Never overload DME again

Weight Capacity Labels

Identify equipment weight capacity at a glance. Large type can be easily seen from 6-8 feet away. Preserves patient dignity by subtly stating the weight capacity.

#DNN71 1511 Variety Pack Labels:
175/pk: (25 labels each of 30, 40, 50,
60, 70, 80, 100)\$70.00 pk
.....\$61.00 pk (3+)

Details, p4



Save 40% on sling replacement

Replacement Slings

fit Hoyer®, Invacare®, Apex®, and most other lifters.

#DNN78252 Medium\$102.25 ea
#DNN78253 Large106.25 ea
#DNN78257 X-Large.....113.50 ea

Details, p8



Full-Body, Solid Seat

Portable lift for home health care

Molift Smart Folding Lift

fits into a closet or car trunk, but has big lift capabilities. Versatile, compact, lightweight, portable, and easy to carry and move. Folds without tools and is easy to store.

#DNN71 1395\$3,510.00 ea

Details, p6



Bed positioning made easy

"No-Lift" BOOSTER™ and TURNER™ enable caregivers to safely handle residents without risking injury to themselves or to the patient.

BOOSTER™ lets the caregiver boost patients in bed (or geri-chair) without risk of back strain, and with fewer staff. Protect from common injuries associated with boosting.

TURNER™ allows one staff member to safely turn a patient in bed without lifting.

Details, p16

