

## One alarm, two functions



*Pull-Cord or Sensor option*

- **Dual functionality—**  
choose sensor or pull-cord
- **Compact and easy to use**
- **Extended battery life**

### **IQ Duo Plus Alarm**

The dual-functioning IQ Duo Plus Alarm operates as either a sensor or pull-cord alarm, enabling staff to choose the best option for their patients' needs. Durable construction and an easy-to-use, compact design make the IQ Duo Plus Alarm an integral part of any Fall Management program.


**Once activated, the alarm's protected settings** and "Always On" technology eliminate the need for reconfiguration, and prevent staff from accidentally switching off the power.

**The alarm's extended battery life** means less frequent battery replacement and fewer occurrences of low-battery warnings.

**Compatible with AliMed sensors,** the IQ Duo Plus works in tandem with pressure and ultra-thin sensor pads, roll-control bed belts, seatbelts, or cushions to create a complete fall monitoring system.

**To order, see reverse side.**

# IQ Duo Plus Alarm



## Alarm Specifications:

<b>Functionality</b>	Pull-Cord or Sensor	Pull-cord is included (cord adjusts to 28" -58") For use with a sensor, see below
<b>Low Battery Indicator</b>	Yes	Flashes red when batteries need to be changed
<b>Reset Button</b>	Yes	Reset button is located on the front of the unit
<b>Nurse-Call Integration</b>	Yes	Compatible with most nurse-call systems
<b>Number of Tones</b>	2	Chime and alarm
<b>Silent Setting</b>	Yes	For use with nurse-call
<b>Number of Delays</b>	2	0 sec., 2 sec.
<b>Volume</b>	High/Low	100dB/85dB
<b>Size</b>	Compact	2 1/2" W x 3 3/4" H x 1" D
<b>Power Source</b>	3 AAA Batteries	Included, average 6 months between battery changes*
<b>Mounting Hardware</b>	Yes	Holster with clip, strap
<b>Warranty</b>	Yes	1 year

\*Results may vary, based on alarm use.

## Compatible with AliMed sensors:

### Pressure Sensor Pads

- 30-day bed and chair pads
- 6-month bed and chair pads
- 12-month bed and chair pads

### Ultra-Thin™, Touch-Sensitive Sensor Pads

- 30-day pads
- 6-month pads

### Belts

- Roll-control belts for beds
- Seatbelts for chairs

### Cushions

- BreakAway Lap cushion
- Sit-Straight™ wedge cushion
- Sensor cushion

**Ordering Information:** #711290 IQ Duo Plus Alarm

Place your order before 3 p.m. and, if in stock, it ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #6351REV0918 Product #11141

©2018 AliMed® Inc. All rights reserved. • Printed in USA



AliMed.com/fm

Call: 800.225.2610

Fax: 800.437.2966

