

Simple. Disposable. Economical.

- **Single-patient-use helps prevent cross-contamination**
- **Easy-to-use with no setup**
- **30-day permanent battery**

Worry-Free[®] Pull-Cord Alarm

With the simplest and easiest functionality, the Worry-Free Pull-Cord Alarm requires no setup—simply switch it on and it's ready for immediate use providing fast, dependable fall monitoring.


The alarm's permanent battery lasts up to 30 days under typical usage, making it an economical and disposable option for single-patient-use needs.

Equipped with an adjustable clip for easy mounting, the Worry-Free Pull-Cord Alarm is a perfect local alarm for use with wheelchairs, chairs, or beds not connected to nurse-call systems.

To order, see reverse side.



Worry-Free® Pull-Cord Alarm



Pull-Cord Alarm Specifications:

Functionality	Pull-Cord	Adjustable length pull-cord
Volume	Set at High	100dB
Size	Compact	1½" W x 3" L x 1½" D
Power Source	Permanent Battery	30 day battery life under typical usage*
Mounting Hardware	Yes	Clip
Warranty	Yes	30 days

*Results may vary, based on alarm use.

Ordering Information: #72140 Worry-Free® Pull-Cord Alarm

Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

The AliMed Advantage

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items**
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

** Place your order before 3 p.m. and, if in stock, it ships the same day!
Prices, product availability, and specs are subject to change.

Pub. #6354REV0918 Product #11146

©2018 AliMed® inc. All rights reserved. • Printed in USA



AliMed.com/fm

Call: 800.225.2610

Fax: 800.437.2966

