


**INTRODUCING!**  
New AliGuard model

## Premium features at an economical price

**Free of powder, latex\* and lead**

Our **NEW** premium AliGuard Glove provides high attenuation and excellent dexterity at a lower price than leading models



\*Not made with natural rubber latex.

- Ultra-thin formulation preserves dexterity, sensitivity, and flexibility

- Curved fingers offer greater precision and less fatigue

- Choose premium AliGuard or economical Original

### AliMed Radiation Attenuation Glove Family

Both gloves provide high attenuation and excellent tactile sensitivity without harmful skin irritants such as powder and natural rubber latex, reducing allergens for complete hand protection. Curved fingers help lessen fatigue during procedures and a coated interior makes them easy to don. Lead-free material requires no special disposal for a cost-effective and environmentally friendly alternative to traditional leaded models.

**AliGuard** Our new premium glove is our highest attenuation model, offering superior protection with finely textured fingers and palm for a better grip on delicate instruments—all at a lower price than other comparable models.

**Original** Our tried and true, the Original is our most economical glove without sacrificing attenuation or tactile sensitivity.

See reverse for details and ordering.

Try AliMed Gloves  
**FREE!**



Request your sample and watch the video at:  
[AliMed.com/glove-sample](http://AliMed.com/glove-sample)

# AliMed® Radiation Attenuation Glove Family

- ✓ Lead-free, powder-free, latex-free\*
- ✓ Ultra-thin for greater tactile sensitivity
- ✓ Curved fingers for reduced fatigue

**NEW!**

## Premium features at a lower cost

### A. AliMed AliGuard Gloves

Premium option with 52% attenuation and finely textured fingers and palm for a better grip on delicate instruments.

#960842	SIZE 6	#960846	SIZE 8
#960843	SIZE 6½	#960847	SIZE 8½
#960844	SIZE 7	#960848	SIZE 9
#960845	SIZE 7½	#960849	SIZE 9½

See glove comparison chart below.

MD



AliGuard's finely textured fingers and palm offer advanced tactile sensitivity and dexterity

## Most economical

### B. AliMed Original Gloves

Most economical option without sacrificing attenuation, tactile sensitivity, or dexterity.

#95-493	SIZE 6½	#95-496	SIZE 8
#95-494	SIZE 7	#95-497	SIZE 8½
#95-495	SIZE 7½	#95-498	SIZE 9

See glove comparison chart below.

MD



Economical Original for the cost-conscious customer

## Attenuation Gloves Comparison Chart

Compare	Attenuation at 60 kVp <sup>†</sup>	Attenuation at 120 kVp <sup>†</sup>	Latex-Free*	Lead-Free	Powder-Free	Pb-Equivalent at 60 kVp <sup>†</sup>	Beaded Cuff	Length	Shelf Life	Thickness	Color	Qty/Box
<b>A. AliMed AliGuard</b> <span style="color:red">NEW!</span>	52%	36%	✓	✓	✓	0.050 mm		11"	3 years	0.25 mm	Grey	5 pairs
<b>B. AliMed Original</b>	47%	27%	✓	✓	✓	0.018 mm		11¾"	3 years	0.23 mm	Yellow	5 pairs
<b>C. XGuard RR1</b>	45%	23%			✓	0.015 mm		12"	2 years	0.22 mm	Brown	5 pairs
XGuard RR2	55%	31%			✓	0.021 mm		12"	2 years	0.30 mm	Brown	5 pairs
<b>D. Secure Touch XR1</b>	46%	24%		✓	✓	0.020 mm	✓	12"	3 years	0.20 mm	Khaki	5 pairs
Secure Touch XR2	64%	42%		✓	✓	0.040 mm	✓	12"	3 years	0.35 mm	Khaki	5 pairs
<b>E. FreeGuard 1</b>	42%	25%	✓	✓	✓	0.024 mm		12"	2 years	0.22 mm	Black	5 pairs
FreeGuard 2	50%	31%	✓	✓	✓	0.032 mm		12"	2 years	0.30 mm	Black	5 pairs

<sup>†</sup>At fingertip. \*Not made with natural rubber latex.

## Try AliMed Gloves

**FREE!**



Request your sample and watch the video at:  
[AliMed.com/glove-sample](http://AliMed.com/glove-sample)

## More Radiation Attenuation Gloves

The following models are offered in two levels of attenuation and thickness to meet your procedural needs.

### C. XGuard® Attenuation Gloves

- Lead-lined, powder-free
- Thinner, texture-grained fingertips for enhanced tactile sensitivity

See glove comparison chart above.

#925920	RR1
#925921	RR2

Specify size: 6, 6½, 7, 7½, 8, 8½, or 9.

MD



### D. Secure Touch® Attenuation Gloves

- Lead- and powder-free with hypoallergenic liner
- Textured fingers and palms

See glove comparison chart above.

#937947	XR1 ULTRA-THIN
#937948	XR2

Specify size: 6½, 7, 7½, 8, 8½, or 9.

MD



### E. FreeGuard® 1 & 2 Attenuation Gloves

- Surgical and radiation protection glove all-in-one!
- Lead-free, powder-free, latex-free

See glove comparison chart above.

#925923	FREEGUARD 1
#925924	FREEGUARD 2

Specify size: 6, 6½, 7, 7½, 8, 8½, or 9.



\*Not made with natural rubber latex.

Any order placed before 3 p.m. EST, if in stock, ships the same day!  
Prices, product availability, and specs subject to change.

Pub. #6935REV0223 Product #11192  
©2022 AliMed® Inc. All rights reserved.

