



## High-quality PPE storage *without* the costly price tag

Easily accessible PPE for both staff and patients is critical for compliance while reducing the spread of infection. But the cost of must-have PPE organizers can be high depending on your facility's needs. AliMed's complete line of **Infection Prevention Stations** and **PPE Organizers** offer all the durability and high-quality storage options you need without the hefty price tag.

### Patient Protection

For patient/public use in waiting areas, points of entry, lobbies, physician offices, and more



#### A./B. Infection Prevention Stations

- Holds earloop face masks, tissues or gloves, and hand sanitizer bottle or wipes canister
- Choose Counter/Wall model or with Stand

Download a **FREE** patient education sign at [AliMed.com/ic-signs](http://AliMed.com/ic-signs)

### Staff Protection

For staff use in isolation or procedure rooms, patient care or exam areas, or wherever PPE is required



#### C. Isolation Stations

- Stores gloves, gowns, and earloop face masks
- Mount on the wall or use with optional Door Hanger (F.), sold separately


#### D. PPE Organizers

- Stores gloves, gowns, and earloop and shield face masks
- Mount on the wall or use with optional Door Hanger (F.), sold separately

#### E./F. Bulk PPE Organizers and Door Hanger

- Stores gloves and earloop face masks with additional storage for bulky items
- Mount on the wall or use with optional Door Hanger, sold separately

**See back for specifications and ordering information**



INFECTION PREVENTION STATIONS and PPE ORGANIZERS		DIMENSIONS	# OF PPE BOXES EACH UNIT HOLDS				
ITEM #	DESCRIPTION	OVERALL (W x H x D)	GLOVES	GOWNS	SHIELD FACE MASKS	EARLOOP FACE MASKS	ADDITIONAL BINS
A. #960877	INFECTION PREVENTION COUNTER/WALL STATION	18 $\frac{1}{4}$ " x 12 $\frac{3}{8}$ " x 5"	1	—	—	2	—
B. #960828	INFECTION PREVENTION STATION w/STAND	18 $\frac{1}{4}$ " x 58 $\frac{1}{4}$ " x 18"	1	—	—	2	—
C. #960826	ISOLATION STATION	17 $\frac{1}{4}$ " x 25 $\frac{1}{2}$ " x 4 $\frac{1}{2}$ "	3	1	—	1	—
D. #960827	PPE ORGANIZER	23 $\frac{1}{4}$ " x 20 $\frac{1}{16}$ " x 4 $\frac{5}{8}$ "	2	1	1	1	—
E. #960825	BULK PPE ORGANIZER	22" x 19 $\frac{1}{8}$ " x 4 $\frac{1}{2}$ "	3	—	—	1	2
F. #960824	OPTIONAL DOOR HANGER FITS #960826, #960827, #960825						

Cold-rolled steel • Beige • Consumables not included • Full specifications at [AliMed.com/ic](http://AliMed.com/ic)

## AliMed® Glove Dispensers



**Wire Glove Dispensers** Wire construction allows easy identification of contents. Assorted sizes fit up to 4 of any size glove box. Keyholes for easy wall mounting. Available in **Horizontal** or **Vertical** models.

$\frac{3}{16}$ " thick • Coated 6-gauge steel • Mounting hardware included  
• Gloves not included • White

ITEM #	DESCRIPTION	W	H	D
<b>HORIZONTAL</b>				
#922657	SINGLE	11 $\frac{1}{4}$ "	4"	3 $\frac{1}{4}$ "
#924744	DOUBLE	11"	11"	3 $\frac{1}{2}$ "
#924746	TRIPLE	11"	14"	3 $\frac{1}{2}$ "
#924748	QUADRUPLE	11"	22"	3 $\frac{1}{2}$ "
<b>VERTICAL</b>				
#924745	DOUBLE	11"	8"	3 $\frac{3}{4}$ "
#924747	TRIPLE	16"	8"	3 $\frac{3}{4}$ "
#924749	QUADRUPLE	22"	8"	3 $\frac{3}{4}$ "



**Clear Plastic Glove Dispensers** Made of clear, durable PETG plastic for easy viewing. Fits glove boxes of any size with keyholes for easy mounting on any wall horizontally or vertically. **Single** model is spring-loaded. 2/cs.

$\frac{3}{32}$ " thick • PETG plastic • Mounting hardware included  
• Gloves not included • Clear

ITEM #	DESCRIPTION	W	H	D
#935935	SINGLE	5 $\frac{1}{2}$ "	10"	3 $\frac{3}{4}$ "
#935936	DOUBLE	10 $\frac{1}{2}$ "	10"	3 $\frac{3}{4}$ "
#935937	TRIPLE	15 $\frac{3}{4}$ "	10"	3 $\frac{3}{4}$ "
#935939	QUADRUPLE	22 $\frac{1}{4}$ "	10"	3 $\frac{3}{4}$ "

Any order placed before 3 p.m. EST, if in stock, ships the same day!  
Prices, product availability, and specs are subject to change.

Publication #7839 ©2020 AliMed® Inc. All rights reserved. • Printed in USA.