



Packaging Includes

- Door Hanger (no assembly required)
- Rubber stickers to protect door

Product Specifications

- 2"W x 12"H (5 cm x 30.5 cm)
- Weight: 1 lb. (0.45 kg)
- Hangs over standard door up to 2 $\frac{7}{16}$ " thick (6.2 cm)
- Beige

For Use With


AliMed Isolation Stations and PPE Organizers:

- #960826 AliMed Isolation Station
- #960827 AliMed PPE Organizer
- #960825 AliMed Bulk PPE Organizer

Care and Cleaning Instructions

- Clean with a non-abrasive cloth using a diluted solution of mild, non-abrasive detergent or standard hospital disinfectants.
- Dry with a soft, non-abrasive cloth.

Mounting Instructions



Visit [AliMed.com/ic](https://www.AliMed.com/ic) for our full line of PPE Organizers

Unit is intended to hold specified AliMed PPE dispensers. Not recommended for other uses.
Content provided is subject to change. Drawing dimensions are approximate.

Pub. #7859 Product #11567
©2020 AliMed® Inc. All rights reserved.