

# A safer way to transfer patients

- **Greater comfort and stability—even for your largest patients**
- **Economical—use for multiple transfers per patient stay**

### PPS Glide Air-Assisted Lateral Transfer System

Similar to an air hockey table, the innovative air-assisted technology of the PPS Glide reduces friction by circulating low pressure, high-volume air through hundreds of micro-perforations. This creates a cushion of air to help smoothly, comfortably, and safely lift and glide patients from surface to surface, drastically reducing pull force and the risk of staff injuries compared to traditional transfer devices.

Offered in both a reusable and disposable option, the PPS Glide is ideal for any hospital protocol.

Both options can remain with the patient for the duration of their stay for multiple transfers or as a repositioning aid for boosting, turning, or proning.

**See reverse  
for details** 



Roller Storage Tote sold separately (compatible with both models)



Storage Cart sold separately (compatible with both models)

## PPS Glide Single Patient Use



UP TO  
**1200**  
lbs.



**PPS Glide Single Patient Use** is completely disposable and best for areas where faster cleanup and cross-contamination are a greater concern, such as in the operating room or emergency department.

**Glide:** 1,200-lb. capacity • **MR-safe** • **Pump:** MR-unsafe • 1-year warranty

- #960926 Regular, 34"W x 78"L 10/cs
- #960927 Large, 39"W x 78"L 10/cs
- #960928 Bariatric, 50"W x 78"L 5/cs

### Accessories (sold separately)

#960929 120V Pump

Optional Roller Storage Tote and Storage Cart at [AliMed.com](http://AliMed.com)

## PPS Glide Reusable



UP TO  
**1000**  
lbs.

**PPS Glide Reusable** is best for use on patient floors during longer-term stays where devices are less likely to be soiled and can be easily machine-washed between patients, reducing facility waste.

**Machine-washable** • **MRI-compatible** • 1-year limited warranty • **Regular:** 500-lb. capacity • **Large:** 700-lb. capacity • **Bariatric:** 1,000-lb. capacity • **Pump:** specs at left

- #938605 Regular, 28"W x 78"L
- #938606 Large, 32"W x 78"L
- #938607 Bariatric, 46"W x 78"L

### Accessories (sold separately)

#960929 120V Pump

Optional Roller Storage Tote, Storage Cart, and Protective Sheets at [AliMed.com](http://AliMed.com)

## Your everyday needs. Our unique capabilities.

AliMed is your **go-to source** for **everyday, hard-to-find, or made-to-order** medical products designed to improve patient outcomes, spanning the continuum of care.



- ▶ Friendly, experienced, customer-first service
- ▶ U.S.-based and in-house to help
- ▶ Exceptional responsiveness



- ▶ Contract pricing—affiliated with major IDNs and GPOs
- ▶ Volume discounts
- ▶ Order by phone, email, fax, or EDI



- ▶ 98% in-stock rate—we have thousands of items ready to ship
- ▶ Same-day shipping on in-stock items\*
- ▶ No minimum requirements—ever!

\*Any order placed before 3 p.m. EST, if in stock, ships the same day!  
Prices, product availability, and specs are subject to change.

Pub. #7934REV0222 • Product #11592  
©2022 AliMed® Inc. All rights reserved.



[AliMed.com](http://AliMed.com)

[customerservice@AliMed.com](mailto:customerservice@AliMed.com)

800.225.2610