

# Premium pressure relief at an economical price



**Head and Neck**



**Upper Extremity**



**Lower Extremity**



**Body**

*See inside for product listing* 



# AliBlue™ Gel

## 100% Gel Positioners

AliBlue Gel Positioners feature the same premium pressure-relieving properties as other leading gels on the market but are more economical for a cost-effective positioning option. Made of 100% polyurethane-based viscoelastic gel. 18-month warranty.



100% viscoelastic gel



Impervious to fluids




Easy to clean



Shop our full line of AliBlue™ Gel, including position-specific sets

[AliMed.com/aliblu-gel-positioners](http://AliMed.com/aliblu-gel-positioners)



Adult, Standard

Adult, High

Pediatric

Neonatal

Adult, Standard

Adult, High

Pediatric

Neonatal

### A. Head Donuts

### B. Horseshoes



### C. Contoured Headrests



### D. Prone Headrests



### E. Contoured Armboard Pads



### F. Standard Armboard Pads



### G. Mayo Stand Covers



### H. Toboggan Liners



### I. Ulnar Forearm Protectors, Rectangle



### J. Ulnar Nerve Protectors, Oval



### K. Graduated Thoracic Positioners



### L. Axillary Rolls



### M. Chest Rolls



### N. Contoured Chest Rolls



### O. Dome Positioners



### P. Knee Crutch Pads



### Q. Heel Cup Positioners

HEAD AND NECK	EACH*
<b>A. HEAD DONUT</b>	
ADULT, STANDARD	#2970017924
ADULT, HIGH	#2970017986
PEDIATRIC	#2970017925
NEONATAL	#2970017988
<b>B. HORSESHOE</b>	
ADULT, STANDARD	#2970017938
ADULT, HIGH	#2970017993
PEDIATRIC	#2970017939
NEONATAL	#2970017999
<b>C. CONTOURED HEADREST</b>	#2970017937

MD • 65 California residents, see below.

HEAD AND NECK (CONTINUED)	EACH*
<b>D. PRONE HEADREST</b>	
MEDIUM	#2970017935
LARGE	#2970017936
<b>UPPER EXTREMITY</b>	<b>EACH*</b>
<b>E. CONTOURED ARMBOARD PAD</b>	
SMALL	#2970017944
MEDIUM	#2970017943
LARGE	#2970017962
<b>F. STANDARD ARMBOARD PAD</b>	
SHORT	#2970017947
LONG	#2970017945
WIDE LONG	#2970017946

MD • 65 California residents, see below.

UPPER EXTREMITY (CONTINUED)	EACH*
<b>G. MAYO STAND COVER</b>	#2970017953
<b>H. TOBOGGAN LINER</b>	#2970017981
<b>I. ULNAR FOREARM PROTECTOR, RECTANGLE</b>	#2970017956
<b>J. ULNAR NERVE PROTECTOR, OVAL</b>	#2970017955
<b>BODY</b>	<b>EACH*</b>
<b>K. GRADUATED THORACIC POSITIONER</b>	#2970017998
<b>L. AXILLARY ROLL</b>	
SMALL	#2970017964
LARGE	#2970017951
X-LARGE	#2970017961

MD • 65 California residents, see below.

BODY (CONTINUED)	EACH*
<b>M. CHEST ROLL</b>	
PEDIATRIC, 2"W x 6"L x 2"H	#2970017980
PEDIATRIC, 2"W x 8"L x 2"H	#2970017983
PEDIATRIC, 2"W x 10"L x 2"H	#2970017992
PEDIATRIC, 2"W x 12"L x 2"H	#2970017982
PEDIATRIC, 3"W x 6"L x 3"H	#2970017987
PEDIATRIC, 3"W x 8"L x 3"H	#2970017985
PEDIATRIC, 3"W x 10"L x 3"H	#2970017984
PEDIATRIC, 3"W x 12"L x 3"H	#2970017942
ADULT, 4"W x 12½"L x 3"H	#2970017941
ADULT, 4"W x 16½"L x 3"H	#2970017929
ADULT, 4"W x 20½"L x 3"H	#2970017940
ADULT, 6"W x 12"L x 5"H	#2970017957
ADULT, 6"W x 16"L x 5"H	#2970017928
ADULT, 6"W x 20"L x 5"H	#2970017927

MD • 65 California residents, see left page.

BODY (CONTINUED)	EACH*
<b>N. CONTOURED CHEST ROLL</b>	
SMALL	#2970017958
MEDIUM	#2970017959
LARGE	#2970017960
<b>LOWER EXTREMITY</b>	<b>EACH*</b>
<b>O. DOME POSITIONER</b>	#2970017926
<b>P. KNEE CRUTCH PAD</b>	
SMALL	#2970017949
MEDIUM	#2970017948
LARGE	#2970017952
<b>Q. HEEL CUP POSITIONER</b>	#2970017950

MD • 65 California residents, see left page.

\*All products sold individually unless otherwise noted.

**AliBlue™ Gel**

More styles and pressure mapping

See back cover

California Proposition 65 Statement  
65 ⚠️ WARNING: This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.

\*All products sold individually unless otherwise noted.

# Specialty AliBlue™ Gel Products

LOWER EXTREMITY	EACH*
<b>R.</b> ANKLE STIRRUP PAD, PAIR	
SMALL	#2970017930
LARGE	#2970017989
<b>S.</b> LITHOTOMY BOOT LINER, PAIR	#2970017996
<b>T.</b> FRACTURE FOOT PAD	#2970017997
OVERLAY	EACH*
<b>U.</b> VAC PAC UNIVERSAL COVER	#2970017954

MD • 65 California residents, see below.



**R. Ankle Stirrup Pads**



**S. Lithotomy Boot Liners**

\*All products sold individually unless otherwise noted.

Scan to shop the full line of AliBlue™ Gel



**T. Fracture Foot Pads**



**U. Vac Pac Universal Covers**



## Small puncture in your gel? We got it covered!

**Gel Repair Kits** extend the life of your gel positioners with easy-to-apply patches that quickly repair small cuts and punctures. 4"W x 4"L. Cut to size.


#2970021380 GEL REPAIR KIT 12/PACK

MD • 65 California residents, see below.

## AliBlue Gel vs. other leading gels

In a side-by-side comparison illustrated by the pressure mapping images at right, AliBlue Gel performed as well as or better than the leading premium and economical gels in relieving pressure on the head, shoulders, pelvis, and heels.

Test subject: Male, 6'0", 185 lbs.  
Position: Supine



Any order placed before 3 p.m. EST, if in stock, ships the same day!  
Prices, product availability, and specs subject to change.

California Proposition 65 Statement  
65 ⚠️ **WARNING:** This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.



Pub. #7939REV0124 Product #11594  
©2024 AliMed® Inc. All rights reserved.

**AliMed.com**  
Expanded product offering online

**800.225.2610**  
Knowledgeable customer care team