



# Lighter for easier staff handling and fewer injuries



**Head Donuts and Horseshoes**  
Available in Adult and Pediatric




**Body Positioners**  
Chest Rolls shown



**Upper Extremities**  
Standard and Contoured Armboards shown



**Lower Extremities**  
Heel Pads shown



**Positioning Sets**  
Basic Positioning Set shown

■ **Lightweight gel-topped foam**



■ **Impervious to fluids**



■ **Bacteriostatic surface**



## AliLite Positioners

AliLite Precontoured Positioners reduce pressure by distributing weight over a larger surface area. Up to 50% lighter than comparable gel-only products, AliLite gel-topped foam positioners are easier to handle, impervious to fluids, naturally bacteriostatic, and feature a stable, nonskid base. AliLite gel is made of dry viscoelastic polymers and includes a 24-month warranty against manufacturing defects.

### For all your patient positioning needs:

- Head and Neck
- Body
- Upper Extremity
- Lower Extremity
- Overlays
- Positioning Sets

See reverse for a complete listing and ordering information.

# AliLite™ Positioners

DESCRIPTION	SIZE	ITEM #
<b>HEAD AND NECK</b>		
DONUT	ADULT	#931555
DONUT	PEDIATRIC	#931556
HORSESHOE	ADULT	#931553
HORSESHOE	PEDIATRIC	#931554
PRONE HEADREST	ADULT	#931696
PRONE HEADREST	PEDIATRIC	#931697
SUPINE HEADREST	ADULT	#931874
SUPINE HEADREST	PEDIATRIC	#937884
<b>BODY</b>		
CHEST ROLL	6"W x 12"L x 4"H	#931559
CHEST ROLL	6"W x 16"L x 4"H	#931558
CHEST ROLL	6"W x 20"L x 4"H	#931557
CHEST ROLL	6"W x 16"L x 6"H	#937881
CHEST ROLL	6"W x 20"L x 6"H	#937882
DOME POSITIONER		#931871
LATERAL POSITIONER		#931875
PRONE POSITIONER	FULL BODY	#937887
PRONE POSITIONER	TORSO	#937886

DESCRIPTION	SIZE	ITEM #
<b>UPPER EXTREMITY</b>		
ARMBOARD PAD	SHORT	#931551
ARMBOARD PAD	LONG	#931552
CONTOURED ARMBOARD PAD	SHORT	#937885
CONTOURED ARMBOARD PAD	LONG	#931560
<b>LOWER EXTREMITY</b>		
HEEL PADS, PAIR		#931873
KNEE SUPPORT	S	#937889
KNEE SUPPORT	L	#937888
LATERAL LEG POSITIONER		#937883
<b>SURFACE OVERLAY</b>		
SURGICAL SURFACE	SHORT	#931872
SURGICAL SURFACE	MEDIUM	#937880
<b>PATIENT POSITIONING SET</b>		
PATIENT POSITIONING SYSTEM		#932990
LATERAL POSITIONING SET		#937922
PRONE POSITIONING SET		#937921
SUPINE POSITIONING SET		#937925
CARDIOVASCULAR POSITIONING SET		#937924



## Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

### The AliMed Advantage

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items\*
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

\*Any order placed before 3 p.m. EST, if in stock, ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #7595 Product #11420

©2018 AliMed® Inc. All rights reserved. • Printed in USA.



Email: [customerservice@AliMed.com](mailto:customerservice@AliMed.com)

Call: 800.225.2610

Fax: 800.437.2966