

## Better pressure distribution

Plus, more comfort and durability  
than traditional foam



*T-Foam™ retains its compressed shape long enough for capillary vessels to nourish skin tissue*

### Best for patients who:

- ✓ Sit for long periods
- ✓ Are prone to pressure injuries
- ✓ Need postural support
- ✓ Have fragile skin



Choosing the right cushion – read our blog:  
[AliMed.com/choosing-a-cushion](https://www.AliMed.com/choosing-a-cushion)

■ **Compressed shape**—allows capillaries to nourish skin


■ **Open-cell structure** helps reduce pressure points

■ **Softens in reaction to body heat**, molds to user's shape

### T-Foam Wheelchair Cushions

Unlike standard foam wheelchair cushions, “slow recovery” T-Foam (memory foam) Cushions retain their compressed shape long enough to allow the user's capillary vessels to nourish tissues on the weight-bearing surface, improving blood flow and reducing the risk of pressure injuries. With its closely interconnected cells, T-Foam's open-cell structure also provides better impact protection and pressure distribution than standard foam and permits air to circulate freely, preventing heat buildup and perspiration.

Available in three firmness ratings—**Soft, Medium, or Hard**—for optimal pressure distribution based on the weight of the user. Choose cushions with basic T-Foam or combined with a Solid Seat Insert™ or high-density foam base to avoid bottoming out.

See reverse for  
ordering information 


# T-Foam™ combinations for superior performance



## Basic T-Foam Cushions

All-day comfort and pressure relief  
Black stretch knit cover • Meets CAL 117


T-Foam Cushions		Soft	Medium	Hard
<b>2"H</b>	18"W x 16"D	#163S	#163M	#163H
	16"W x 16"D	—	#163M-JR	—
<b>3"H</b>	18"W x 16"D	#164S	#164M	#164H
	16"W x 16"D	#164S-JR	#164M-JR	—
<b>4"H</b>	20"W x 16"D	—	#2125	#2126
	18"W x 16"D	#172S	#172M	#172H



## T-Foam Cushions with Solid Seat Insert™

Stable base, choose your firmness and thickness  
Add 1½"H of foam plus insert • Black stretch knit cover with nonslip bottom • Meets CAL 117


- #1222 18"W x 16"D x 1"H, Soft
- #1224 18"W x 16"D x 1"H, Medium
- #1225 16"W x 16"D x 1"H, Medium
- #1228 18"W x 16"D x 2"H, Soft
- #1229 16"W x 16"D x 2"H, Soft
- #1231 16"W x 16"D x 2"H, Medium



## Deluxe T-Foam Combo™ Cushions

Soft comfort without bottoming out  
3"H overall • Navy water-resistant cover with nonslip bottom  
• Meets CAL 117

- #1175 18"W x 16"D
- #1176 16"W x 16"D



## T-Foam Combo™ Solid Seat Insert™ Cushions

Maximum stability without bottoming out  
3½"H overall • Black water-resistant cover with nonslip bottom  
• Meets CAL 117

- #1200 18"W x 16"D
- #1201 16"W x 16"D

## Custom Cushions to YOUR specifications

We can customize any cushion in a range of sizes, with foam combinations, and more. Call us to learn more.

Visit [AliMed.com/premium-foam-cushions](http://AliMed.com/premium-foam-cushions) for our full line of wheelchair cushions.

Any order placed before 3 p.m. EST, if in stock, ships the same day!  
Prices, product availability, and specs are subject to change.

Pub. #7237REV0519 Product #11214  
©2019 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com / customerservice@AliMed.com / 800.225.2610