

Optimize cervical visualization for sedated patients



Watch a video demonstration!
AliMed.com/cvh

- Reduce risk of incidental trauma
- Use during intraoperative X-rays
- Single use for infection control

Cervical Visualization Harnesses

Retract shoulders for better visualization of C5–6, C6–7, and other vertebrae for intraoperative X-rays. Harness can be pulled to either side by hook-and-loop strap at the chest to expose pelvis for bone graft. The body portion of the harness is held by the patient's weight to help eliminate incidental trauma caused by having to strap the patient down.

Adjustable hook-and-loop straps extend from the back of the body, up and over the shoulders, and then crisscross in front before attaching to the sides of the table just above the ankles. May also be used as a safety harness for CT and MRI procedures. One size fits all.

Pad: 18½"W x 44"L • Straps: 66"L • Disposable • Nonskid foam • Nonsterile • MR-safe

Ordering Information

#920877

Why AliMed?



Customer-First Service

- ▶ **Friendly, personal, experienced**
We're ready to assist you with product questions, suggestions, or special requests.
- ▶ **U.S.-based and on site**
You'll speak with real people who have real answers—with no wait!
- ▶ **Fast response, fast resolution**
We'll make every effort to solve your issue—first call, first time.



Fulfilling Your Every Need

- ▶ **We have what you need now...**
Curated by expert product managers, we have a vast in-stock portfolio of everyday products.
- ▶ **Or we can find it...**
We can often source hard-to-find items through our many partners.
- ▶ **Or even customize it**
With on-site capabilities, we can customize just about anything.



In Stock and Out the Door

- ▶ **98% in-stock rate**
We have thousands of in-stock items ready to ship when you're ready to order.
- ▶ **Same-day shipping**
In-stock orders placed by 3 p.m. EST ship same day, every day!
- ▶ **No minimum order**
Need 1 or 1,000? No problem and no minimum requirements—ever!



*Call us today
to see how
we can help you.
800.225.2610*

Any order placed before 3 p.m. EST, if in stock, ships the same day!

Pub. #7519REV1219 Product #11385

©2019 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com / customerservice@AliMed.com / 800.225.2610