



WHAT IS YOUR IDEAL

**BACTERIA**

**KILL TIME?**



# **PDI SANI-CLOTH DISPOSABLE WIPES**

EPA-registered, Hospital-grade germicidal solution for disinfecting hard, nonporous surfaces

Wipe Sizes: Large: 6" W x 6 3/4"L; X-Large: 8"W x 14"L



## **Kill Time: 2 minutes SUPER SANI-CLOTH**

- High-alcohol formula
  - Tested effective against 22 microorganisms
- |         |                            |
|---------|----------------------------|
| 92 4232 | Large Canister, 12/cs      |
| 92 7775 | Extra-Large Canister, 6/cs |

## **Kill Time: 3 minutes SANI-CLOTH AF3**

- Alcohol-free
  - Tested effective against 44 microorganisms
- |         |                                     |
|---------|-------------------------------------|
| 93 7935 | AF3 Large Canister Case, 12/cs      |
| 93 7936 | AF3 Extra-Large Canister Case, 6/cs |

## **Kill Time: 4 minutes SANI-CLOTH BLEACH**

- 1:10 Sodium Hypochlorite
  - Tested effective against 50 microorganisms
- |         |   |
|---------|---|
| 93 7933 | Bleach, Large Canister Case, 12/cs      |
| 93 7934 | Bleach, Extra-Large Canister Case, 6/cs |

## **Kill Time: Under 5 minutes SANI-CLOTH PLUS**

- Low-alcohol formula
  - Tested effective against 13 microorganisms
- |         |                            |
|---------|----------------------------|
| 92 4230 | Large Canister, 12/cs      |
| 92 4231 | Extra-Large Canister, 6/cs |

## **Kill Time: 10 minutes or less SANI-CLOTH HB**

- Alcohol-free
  - Tested effective against over 100 microorganisms
- |         |                            |
|---------|----------------------------|
| 92 4756 | Large Canister, 12/cs      |
| 92 4757 | Extra-Large Canister, 6/cs |