

Notifies with a personalized message



- Soothing, familiar voice reduces patient anxiety
- Record in any language—perfect for non-English speaking patients
- Easily record sounds

IQ VOICE Alarm

With personalized voice messaging capabilities, the IQ VOICE Alarm plays familial instructions or calming sounds in lieu of startling alarms to encourage patients to remain seated or in bed, reducing patient anxiety and falls.

Feature-rich and easy-to-use, the alarm gives caregivers the option to record and play a 20-second message before reverting to one of four tones—helping to reduce alarm fatigue in facilities with multiple monitoring systems.

Once activated, the alarm's protected settings and "Always On" technology eliminate the need for staff reconfiguration, ensuring the alarm is functioning properly.

With easy mounting to any bed or chair, the IQ VOICE Alarm can be used as a pull-cord or sensor alarm, which is compatible with AliMed sensors, including pressure sensor pads, roll-control belts, seatbelts, or cushions to create a complete fall monitoring system.

To order, see reverse side.

IQ VOICE Alarm



Voice Alarm Specifications:

Functionality	Sensor or Pull-Cord	Pull-Cord included along with storage compartment; for use with a sensor, see below
Alert Indicator	Yes	Located at top of alarm
Low Battery Indicator	Yes	Flashes red when batteries need to be changed
Reset Button	Yes	Dual reset buttons located on both sides of unit; must be pressed simultaneously
Nurse-Call Integration	Yes	Compatible with most nurse-call systems
Number of Tones	4	4 unique alarm tones
Silent Setting	Yes	For use with nurse-call
Volume	High/Low	95dB/90dB
Size	Average	3½"W x 5½"H x 1½"D
Power Source	4 AA Batteries	Included, average 6 months between battery changes*

*Results may vary, based on alarm use.

Compatible with AliMed sensors:

Pressure Sensor Pads

- 30-day bed and chair pads
- 6-month bed and chair pads
- 12-month bed and chair pads

Belts

- Roll-control belts for beds
- Seatbelts for chairs

Cushions

- BreakAway lap cushion
- Sit-Straight™ wedge cushion
- Sensor cushion

Ordering Information: #2032 IQ VOICE Alarm

Place your order before 3 p.m. EST and, if in stock, it ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #5098REV0419 Product #11143

©2019 AliMed® Inc. All rights reserved. • Printed in USA



Email: customerservice@AliMed.com

Call: 800.225.2610

Fax: 800.437.2966