


Style: Coat

Full front/back protection in a one-piece style with fully adjustable 6"W Back-Saver Belt for lower lumbar support, helping to relieve back strain and fatigue by redistributing weight from shoulders to waist.



1. Apron Weight Options: Coats offer full 0.5 mm Pb-equivalent front protection and 0.25 mm Pb-equivalent back protection.

Select Material	Weight*	Material	Procedure Length	Disposal Requirements
FeatherLite Lead Free	7.2 lbs.	100% antimony	Long	Environmentally friendly—no special disposal
Ultralight Lead Free	11.9 lbs.	100% antimony	Long	
Lightweight Lead	14.2 lbs.	Lead-based alloy	Short to medium	Requires special disposal or recycling
Standard Lead	16 lbs.	100% lead	Short	

*Based on Medium Front Protection Apron.

2. Sizing Options:


How to Measure for Coat Aprons

Chest (C): Measure circumference around widest part of chest.

Waist (W): Measure circumference of waist.

Hips (H): Measure circumference around widest part of hips.

Length (L): Measure length of apron from the top of shoulder to the hem (suggest 1" above the knee). Adjust according to personal preference.



Sizing Chart for Coat Aprons

Size (Unisex)	Chest Size	Waist Size	Hip Size	Front/Rear Length
S	33" - 36"	34" - 37"	35" - 38"	36"/32"
M	37" - 40"	38" - 41"	39" - 42"	38"/34"
L	41" - 44"	42" - 45"	43" - 46"	40"/36"
XL	45" - 48"	46" - 49"	47" - 50"	42"/38"
Tailored-Fit (Record measurements) →				

Tailored-Fit Measurements

Record your measurements in inches if requesting Tailored-Fit option, and select "Tailored-Fit," below left.

Chest:

Waist:

Hips:

Length:

3. Color Options:

Choose from 400+ colors/patterns! [Click here to see color swatch choices.](#)

Color Swatch #:	Color Name:
-----------------	-------------

4. Thyroid Shield

Attached Thyroid Shield (same weight and color/print as apron)

5. Monogramming

Apron Monogramming

Lettering Style:

Block

ALL CAPS

Script

Line 1 (16 characters max. including spaces and punctuation):

Line 2 (16 characters max. including spaces and punctuation):

Thyroid Shield Monogramming

One Line Only (10 characters max. including spaces and punctuation):

6. Number of aprons you are requesting with the same specifications:

(For Tailored-Fit, please submit a separate form for each order.)

7. Additional Requests:

Please enter any specific special order selections including embroidery symbols, logos, and monogram thread colors:

8. Contact Information

Contact Person:

Sales Rep (if applicable):

Email Address:

Phone:

Facility Name:

Billing Address:

Address 2:

City:

State:

Zip:

(Shipping location, if different than billing)

Facility Name:

Shipping Address:

Address 2:

City:

State:

Zip:

(If applicable)

Account:

PO:

Contact me to process this order via credit card

Pub #7793-E

Perfect Fit Aprons and Thyroid Shields are not returnable. Sales are not final until approved by AliMed. Visit AliMed.com/billing-and-payments for sales tax details. We will email you a quote within 3-4 business days.