


1. Style Options:

Select Standard or Flex.



Standard: Adjustable, crisscross shoulder straps distribute weight evenly, providing extra comfort by lessening stress on neck, shoulders, and back.



Flex: Elasticized stretch back for a closer fit to the body and greater weight distribution across the neck, shoulders, and back.

2. Apron Weight Options:

All apron styles offer full 0.5 mm Pb-equivalent front protection.

Select Material	Weight*	Material	Procedure Length	Disposal Requirements
FeatherLite Lead Free	6.5 lbs.	100% antimony	Long	Environmentally friendly—no special disposal
Ultralight Lead Free	7.2 lbs.	100% antimony	Long	
Lightweight Lead	8.8 lbs.	Lead-based alloy	Short to medium	Requires special disposal or recycling
Standard Lead	11.3 lbs.	100% lead	Short	

*Based on Medium Front Protection Apron.

3. Sizing Options:

How to Measure for Front Protection Aprons


When measuring, stand tall with arms down and feet together.

Chest (C): Measure circumference around widest part of chest.

Waist (W): Measure circumference of waist.

Hips (H): Measure circumference around widest part of hips.

Length (L): Measure length of apron from the sternal notch (center of collarbones) to the hem (suggest 1" above the knee). Adjust according to personal preference.



Sizing Chart for Front Protection Aprons

Select Size (Unisex)	Apron Size (W x L)	Fits Dress Size	Fits Suit Size	Max Chest/Hip	Fits Max. Height
S	20" x 32"	4-6	34" - 36"	36"	5'6"
M	22" x 34"	8-10	38" - 40"	40"	5'10"
L	24" x 36"	12-14	42" - 44"	44"	6'0"
XL	26" x 38"	16-18	46" - 48"	48"	6'2"
XXL	28" x 40"	20-22	50" - 52"	52"	6'4"
Tailored-Fit (Record measurements) →					

Tailored-Fit Measurements

Record your measurements in inches if requesting Tailored-Fit option, and select "Tailored-Fit," below left.

Chest:

Waist:

Hips:

Length:

4. Color Options:

Choose from 400+ colors/patterns! [Click here to see color swatch choices.](#)

Color Swatch #:	Color Name:
-----------------	-------------

5. Thyroid Shield

Attached Thyroid Shield (same weight and color/print as apron)

6. Monogramming

Apron Monogramming

Lettering Style:

Block

ALL CAPS

Script

Line 1 (16 characters max. including spaces and punctuation):

Line 2 (16 characters max. including spaces and punctuation):

Thyroid Shield Monogramming

One Line Only (10 characters max. including spaces and punctuation):

7. Number of aprons you are requesting with the same specifications:

(For Tailored-Fit, please submit a separate form for each order.)

8. Additional Requests:

Please enter any specific special order selections including embroidery symbols, logos, and monogram thread colors:

9. Contact Information

Contact Person:

Sales Rep (if applicable):

Email Address:

Phone:

Facility Name:

Billing Address:

Address 2:

City:

State:

Zip:

(Shipping location, if different than billing)

Facility Name:

Shipping Address:

Address 2:

City:

State:

Zip:

(If applicable)

Account:

PO:

Contact me to process this order via credit card

Pub #7793-C

Perfect Fit Aprons and Thyroid Shields are not returnable. Sales are not final until approved by AliMed. Visit AliMed.com/billing-and-payments for sales tax details. We will email you a quote within 3-4 business days.