


## Style: Reverse Wraparound

Pairs full front/back protection with a stretch-back design. Elasticized back offers a snug fit that hugs the body and shifts weight off the shoulders and to the waist. Hook-and-loop shoulder and waist fastener straps allow for adjustment for a closer fit. Shrug your shoulders prior to securing straps for a “shrug-and-cinch” custom fit.



## 1. Apron Weight Options:

Coats offer full 0.5 mm Pb-equivalent front protection and 0.25 mm Pb-equivalent back protection.

Select Material	Weight*	Material	Procedure Length	Disposal Requirements
FeatherLite Lead Free	8.79 lbs.	100% antimony	Long	Environmentally friendly—no special disposal
Ultralight Lead Free	9.57 lbs.	100% antimony	Long	
Lightweight Lead	12.01 lbs.	Lead-based alloy	Short to medium	Requires special disposal or recycling
Standard Lead	13.47 lbs.	100% lead	Short	

\*Based on Medium Front Protection Apron.

## 2. Sizing Options:


### How to Measure for Reverse Wraparound Aprons

**Chest (C):** Measure circumference around widest part of chest.

**Waist (W):** Measure circumference of waist.

**Hips (H):** Measure circumference around widest part of hips.

**Length (L):** Measure length of apron from the top of shoulder to the hem (suggest 1" above the knee). Adjust according to personal preference.



### Sizing Chart for Reverse Wraparound Aprons

Size (Unisex)	Chest Size	Waist Size	Hip Size	Front/Rear Length
S	33" - 36"	34" - 37"	35" - 38"	36"/32"
M	37" - 40"	38" - 41"	39" - 42"	38"/34"
L	41" - 44"	42" - 45"	43" - 46"	40"/36"
XL	45" - 48"	46" - 49"	47" - 50"	42"/38"
Tailored-Fit (Record measurements) →				

### Tailored-Fit Measurements

Record your measurements in inches if requesting Tailored-Fit option, and select “Tailored-Fit,” below left.

Chest:

Waist:

Hips:

Length:

## 3. Color Options:

Choose from 400+ colors/patterns! [Click here to see color swatch choices.](#)

Color Swatch #:	Color Name:
-----------------	-------------

## 4. Thyroid Shield

Attached Thyroid Shield (same weight and color/print as apron)

## 5. Monogramming

Apron Monogramming

Lettering Style:

Block

ALL CAPS

*Script*

Line 1 (16 characters max. including spaces and punctuation):

Line 2 (16 characters max. including spaces and punctuation):

Thyroid Shield Monogramming

One Line Only (10 characters max. including spaces and punctuation):

## 6. Number of aprons you are requesting with the same specifications:

(For Tailored-Fit, please submit a separate form for each order.)

## 7. Additional Requests:

Please enter any specific special order selections including embroidery symbols, logos, and monogram thread colors:

## 8. Contact Information

Contact Person:

Sales Rep (if applicable):

Email Address:

Phone:

Facility Name:

Billing Address:

Address 2:

City:

State:

Zip:

*(Shipping location, if different than billing)*

Facility Name:

Shipping Address:

Address 2:

City:

State:

Zip:

*(If applicable)*

Account:

PO:

Contact me to process this order via credit card

Pub #7793-F

*Perfect Fit Aprons and Thyroid Shields are not returnable. Sales are not final until approved by AliMed. Visit [AliMed.com/billing-and-payments](http://AliMed.com/billing-and-payments) for sales tax details. We will email you a quote within 3-4 business days.*