

## Safely disinfect your AliBlue™ Gel

Approved Cleaning Products from AliMed



### DisCide® Ultra

Kills pathogens in one minute



### CaviCide® Spray and CaviWipes™ Towelettes

Kills organisms in three minutes



### PureGreen24™

24-hour residual protection, nontoxic



### PDI® Sani-Cloth® Disposable Wipes

Three formulas for every disinfecting need



### Protex® Spray and Protex® Ultra Wipes

Alcohol-, bleach-, and solvent-free

**SHOP NOW!**

#### Other Recommended Generic Cleaners

- Bleach
- Quaternary Ammonium (i.e., ammonium chloride)
- Isopropyl Alcohol
- Chlorohexidine Gluconate
- Povidone-Iodine (This will discolor the membrane slightly but not damage it)

#### Do Not Use

- Hydrogen Peroxide

*Please note: Always follow your facility's cleaning protocol before and after each use.*



## Small puncture in your gel? We got it covered!

**AliBlue Gel Repair Kits** extend the life of your gel positioners with easy-to-apply patches that quickly repair small cuts and punctures. Cut to size. 4"W x 4"L. 12 patches per pack.

**SHOP NOW!**