

**No setup.
No batteries to replace.
No worries.**



- **Easy functionality with no setup**
- **Two-year battery life**
- **Economical**

Worry-Free[®] Sensor Alarm

With effortless functionality, the Worry-Free Sensor Alarm is virtually foolproof with no setup to worry about—simply pull the plastic starter tab and it's ready for immediate use providing fast, dependable fall monitoring.

Once activated, the alarm's "Always On" technology eliminates the need for staff reconfiguration, ensuring the alarm is functioning properly. Plus, you'll never have to replace a battery again—the permanent battery is designed to last up to two years under typical usage and the low-battery indicator is triggered when five days of life remain.

Equipped with a secure clip for easy mounting, the Worry-Free Sensor Alarm is a perfect local alarm for use with wheelchairs, chairs, or beds not connected to nurse-call systems, and is compatible with AliMed sensors, including pressure and ultra-thin sensor pads, roll-control belts, seatbelts, or cushions to create a complete fall monitoring system.

To order, see reverse side.

Worry-Free® Sensor Alarm



Sensor Alarm Specifications:

Functionality	Sensor	For use with a sensor (See below)
Low Battery Indicator	Yes	Turns red when 5 days of battery life remain
Reset Button	Yes	Recessed reset button
Volume	High/Low	100dB/85dB
Size	Compact	2 7/8" W x 4 1/8" L x 1 1/2" D
Power Source	Permanent Battery	2 year battery life under typical usage*
Mounting Hardware	Yes	Clip and strap
Warranty	Yes	30 days

*Results may vary, based on alarm use.

Compatible with AliMed sensors:

Pressure Sensor Pads

- 30-day bed and chair pads
- 6-month bed and chair pads
- 12-month bed and chair pads

Ultra-Thin™, Touch-Sensitive Sensor Pads

- 30-day pads
- 6-month pads

Belts

- Roll-control belts for beds
- Seatbelts for chairs

Cushions

- BreakAway Lap cushion
- Sit-Straight™ wedge cushion
- Sensor cushion

Ordering Information: #713307 Worry-Free® Sensor Alarm

Place your order before 3 p.m. and, if in stock, it ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #7143REV0918 Product #11147

©2018 AliMed® inc. All rights reserved. • Printed in USA



AliMed.com/fm

Call: 800.225.2610

Fax: 800.437.2966

