

Our versatile hook-and-loop table straps Imaging or Surgical

Call for
made-to-order
sizes!



Soft and strong.
Trim to any length.
Four widths available!

- Radiolucent for surgical and imaging environments
- Versatile—available in multiple widths and lengths
- Single-patient use for infection control

Universal Table Straps

Extraordinarily strong and highly versatile for safely securing just about any patient during any procedure. The red hook adheres to either side of the soft blue loop to accommodate any application. Use the **One-Closure** strap for tables without side rails—easily shorten by overlapping or trimming the blue loop. The **Two-Closure** model is ideal for tables with rails, eliminating the need to wrap the strap under the table—simply cinch up the ends for smaller patients. Use the **Three-Closure** model for procedures requiring quick and easy patient access.


Pre-configured to match your table and save you time!

A. Simplest

One-piece • One closure
Encircles patient and table




One closure




B. Strongest

One-piece • Two closures
Loops through side rails




Two closures




C. Easy Rapid Patient Access

Two-piece • Three closures
Loops through side rails



Three closures



See reverse for ordering information →


Strong closure for maximum stability

- Soft against skin
- Versatile
- Single-patient use

Pre-configured to match your table and save you time!


A. Simplest

One-piece • One closure
Encircles patient and table




B. Strongest

One-piece • Two closures
Loops through side rails



C. Easy Rapid Patient Access

Two-piece • Three closures
Loops through side rails



A. SIMPLEST

ITEM #	WIDTH x LENGTH	QUANTITY
#2970011607	2" x 76"	ROLL
#2970011611	2" x 96"	ROLL
#2970013589	2" x 108"	ROLL
#2970011608	3" x 76"	ROLL
#2970011612	3" x 96"	ROLL
#2970013590	3" x 108"	ROLL
#2970011609	4" x 76"	ROLL
#2970011613	4" x 96"	ROLL
#2970013591	4" x 108"	ROLL
#2970011610	6" x 76"	ROLL
#2970011615	6" x 96"	ROLL
#2970013592	6" x 108"	ROLL

MD • 65 California residents, see below.

B. STRONGEST

ITEM #	WIDTH x LENGTH	QUANTITY
#2970013580	2" x 76"	ROLL
#2970013581	2" x 96"	ROLL
#2970013582	2" x 108"	ROLL
#2970011673	3" x 65"	ROLL
#2970013583	3" x 96"	ROLL
#2970013584	3" x 108"	ROLL
#2970011674	4" x 65"	ROLL
#2970013585	4" x 96"	ROLL
#2970013586	4" x 108"	ROLL
#2970011675	6" x 65"	ROLL
#2970013587	6" x 96"	ROLL
#2970013588	6" x 108"	ROLL

MD • 65 California residents, see below.

C. EASY RAPID PATIENT ACCESS

ITEM #	WIDTH x LENGTH	QUANTITY
#2970011616	2" x 47"-76"	ROLL
#2970011617	3" x 47"-76"	ROLL
#2970013593	3" x 76"-108"	ROLL
#2970011618	4" x 47"-76"	ROLL
#2970013594	4" x 76"-108"	ROLL
#2970011619	6" x 47"-76"	ROLL
#2970013595	6" x 76"-108"	ROLL

MD • 65 California residents, see below.

Made-to-order sizes available... **Call 800.225.2610**

Any order placed before 3 p.m. EST, if in stock, ships the same day!
Prices, product availability, and specs subject to change.

California Proposition 65 Statement
65 **WARNING:** This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.



Pub. #7494REV0124

©2024 AliMed® Inc. All rights reserved.