

Accurate results, early detection, improved patient outcomes



Note: Filaments shown in color for visual example only

- Repeatability ensures consistent, accurate results

- Disposable for single-patient use and infection control

Disposable Monofilament

Monofilaments are the most effective testing tool for early detection of sensory loss or neuropathy in diabetic patients. Choosing the right one for consistent, reliable results is critical.

The **AliMed Disposable Monofilament** provides superior testing consistency and accuracy using a proprietary filament that is resistant to environmental fluctuations—unlike nylon filaments, which can be affected by changes in humidity and temperature that can compromise results.


The calibrated 10-gram 5.07 monofilament is attached to a paperboard handle with 5 site testing protocol map on the back.

Perforated blister pack for protective storage.

Ordering Information

#32745 Disposable Monofilament, 20/pk

For volume pricing,
email orders@alimed.com.

See reverse for
more sensory testing
and foot care products. 

Stay a step ahead with assessment tools and protective foot care products from AliMed



- **Tests static and moving 2-point discrimination**

#32640 Touch-Test® Two-Point Discriminators




- **Evaluates current, diminishing, or returning sensory levels**

#63232 Touch-Test® Sensory Evaluators




- **Helps relieve pressure on the foot**

#63135 AliMed® Diabetic Walkers



- **Maximum protection without bottoming out**

#6227 AliMed® Duo/Laminate D-Soles™



- **Heat-moldable for full contact with plantar surface**

#66376 AliMed® FREEDOM® Trilaminare Diabetic Insoles



- **Stabilizing heel and moderate arch support**

#64201 AliMed® FREEDOM® Full Contact Insoles

Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

The AliMed Advantage:

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items*
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

*Place your order before 3 p.m. and, if in stock, it ships the same day!

Prices, product availability, and specs are subject to change.

Pub. #7542REV0418 Product #11395

©2018 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com

Call: 800.225.2610

Fax: 800.437.2966