

Disposable solutions for maximizing infection control

AliStrap® Soft

p. 2

Universal Table Straps

p. 5

AliStrap®

p. 2

AliStrap® Precut

p. 7

**Universal Table Straps** p. 5**Disposable Surgical Straps** p. 4**Universal Econo-Straps** p. 4**AliMed.com**

Expanded product offering online


**800.225.2610**

Knowledgeable customer care team



AliStrap®—the most versatile back-to-back hook-and-loop strap


Casting




Patient Safety




Equipment Cover Security




Body Positioning




Fragile Skin Applications




Pediatric Line Positioning



Catheter Security



Armboard Positioning




The Right Strap. Right Length. Every Time!

AliStrap and **AliStrap Soft** are strong non-adhesive, self-adhering, back-to-back hook-and-loop straps with multiple applications whether you need an everyday strap for securing just about anything or one tailored for more fragile, gentler positioning jobs—**just unroll, cut, wrap, and press!**

The original Our best-selling AliStrap features infinite versatility for any general use application yet is strong enough for even your most difficult strapping jobs, from patient immobilization to heavy-duty equipment security.

A softer solution AliStrap Soft with its fully padded, skin-safe, soft edges is preferred for your most sensitive skin applications. Featuring an internally padded loop with a built-in foam layer nestled between the hook and loop, AliStrap Soft provides extra protection and comfort for the patients who need it most, from pediatric to geriatric.



No sticky residue With no bacteria-trapping sticky tape residue to leave behind, AliStrap can be easily repositioned again and again.

Both styles are disposable for single use and available in convenient cut-to-length rolls in a variety of widths to best suit your desired application.




Precut lengths available—individually wrapped and ready to use! p. 7

AliStrap®—best for general use

AliStrap® With cut-to-fit versatility and a strong, self-adhering back-to-back hook-and-loop closure, AliStrap is our best-selling and most versatile positioning strap that is ideal for securing, immobilizing, or positioning patients, equipment, cords, and more. Low-profile, medical-grade hooks and non-abrasive soft edges will not cut into or abrade skin, instead yielding and flexing against soft tissue to reduce subcutaneous shear. Economical and disposable, these straps are perfect for single use, preventing cross-contamination and offering easy cleanup. Simply cut any length from the continuous 30'L roll, wrap, and press together to form a tight, strong bond—even holds when wet. Available in five widths. Blue.

30'L ROLL	
95-905	¾"W NARROW
95-900	1½"W REGULAR
95-902	3"W WIDE
95-904	5"W X-WIDE
95-908	10"W MAXX

Call for volume pricing or custom sizes.



AliStrap® Soft—preferred for fragile skin

AliStrap® Soft is a softer, smoother alternative to abrasive straps, combining the strength and security of our best-selling AliStrap with fully padded, skin-safe edges that are gentler against sensitive or fragile skin. With a non-adhesive, self-adhering, back-to-back hook-and-loop closure, AliStrap Soft uses low-profile, smoother medical-grade hooks that are less abrasive than traditional hooks yet strong enough to stabilize patients during surgery, imaging, and other procedures. Padded loop center with built-in foam center and flexible, yielding edges reduces skin shear and won't damage sensitive skin. Simply cut any length from the continuous 30'L roll, wrap, and press together to form a tight, strong bond—even holds when wet. Available in two widths. Beige.

30'L ROLL	
930382	1½"W PADDED
930163	2"W PADDED

Call for volume pricing or custom sizes.



1. Unroll




2. Cut to length



3. Wrap and press!



Single-Patient-Use Table Straps



■ Padded foam for pressure relief

Disposable Surgical Straps offer a clean strap that's always ready to use. Strong black web strap with hook-and-loop closure for maximum adjustment. 3/8"-thick hypoallergenic foam pad adheres to strap with hook-and-loop, offering protection between patient and strap. For use on tables with side rails. Not made with natural rubber latex.

Strap: 2"W x 53"L • Pad: 3"W x 25"L x 3/8" thick

937545 DISPOSABLE SURGICAL STRAP



■ Ideal for CT and MRI!

■ Skin-safe padding

Universal Econo-Straps combine six feet of easy-catch loop on one side, doughy-soft padded sponge on the other, and 10" of strong hook material at the end—delivering maximum skin protection, ultimate adjustability, and excellent closure strength. Wrap around any body part from finger to foot. Cut for the best fit. No rough edges to irritate skin. Use for imaging, surgery, or critical care.

72"L • Radiolucent and artifact-free • Metal-free • MR-safe

7536 3"W
7539 4"W

Custom sizes available...CALL!



■ Cut to any length

AliStrap® is a continuous roll of hook-and-loop strapping that cuts to any length you need. Simply wrap around patient and table or surface, then press together to form a tight, secure bond. Adhesive-free strap eliminates the risk of bacteria-trapping residue being left on patients' skin or table surface. Flexible yielding edges are gentle on skin and reduce shear. Not made with natural rubber latex.

30'L roll


95-902 3"W
95-904 5"W
95-908 10"W



■ Secure, comfortable, economical

Covidien Devon™ Disposable Body/Knee Straps Strong fabric construction for secure and comfortable patient positioning. For use on operating room tables and gurneys. Secures with hook-and-loop closure. 12/cs. 3"W x 60"L • Plastic D-ring

938133 DISPOSABLE BODY/KNEE STRAP, 12/cs



■ Soft, plush feel

Halyard O.R. Table/Stretchers Straps

Completely adjustable, single-use strap securely holds patient. Simply feed the strap through the metal ring and press the hook-and-loop fastener to the strap when the preferred tension is achieved. For use on tables and stretchers. Not made with natural rubber latex. 12/bx, 2 bx/cs.

3"W x 60"L • Plush polyester with durable Sontara® layer • Metal D-ring

938134 OR TABLE STRETCHER STRAP, 12/bx, 2 bx/cs




Universal Table Straps

Strong and versatile for maximum stability




Universal Table Straps are extraordinarily strong, highly versatile, and pre-configured to match your table, allowing you to quickly and safely secure patients during any procedure. The red hook side adheres to either side of the soft blue loop for a strong hold that is safe against skin. Disposable for fast, easy cleanup and to help reduce cross-contamination risk. Easily shorten the **One-Closure** strap by overlapping or trimming the blue loop. With **Two-Closures**, secure smaller patients by cinching up the side rail ends. **Three-Closure** strap allows for quick and easy patient access. All are available in multiple widths and lengths for patients of any size.

Radiolucent • Metal-free • MR-safe




Trim to any length—four widths available!



Strong and durable—yet soft against skin


3 configurations to match your table

A. Simplest
One-piece • One closure
Encircles patient and table




A. SIMPLEST			
ITEM #	W x L	ITEM #	W x L
91-277	2" x 76"	91-279	4" x 76"
91-281	2" x 96"	91-283	4" x 96"
932579	2" x 108"	932581	4" x 108"
91-278	3" x 76"	91-280	6" x 76"
91-282	3" x 96"	91-284	6" x 96"
932580	3" x 108"	932582	6" x 108"

B. Strongest
One-piece • Two closures
Loops through side rails



B. STRONGEST			
ITEM #	W x L	ITEM #	W x L
932570	2" x 76"	91-598	4" x 65"
932571	2" x 96"	932575	4" x 96"
932572	2" x 108"	932576	4" x 108"
91-597	3" x 65"	91-599	6" x 65"
932573	3" x 96"	932577	6" x 96"
932574	3" x 108"	932578	6" x 108"

C. Easy Rapid Patient Access
Two-piece • Three closures
Loops through side rails



C. EASY RAPID PATIENT ACCESS	
ITEM #	W x L
91-285	2" x 47"-76"
91-286	3" x 47"-76"
932584	3" x 76"-108"
91-287	4" x 47"-76"
932585	4" x 76"-108"
91-288	6" x 47"-76"
932586	6" x 76"-108"

Optimize cervical visualization for sedated patients

- ✓ Reduce risk of incidental trauma
- ✓ Use during intraoperative X-rays
- ✓ Single use for infection control



Cervical Visualization Harnesses retract shoulders for better visualization of C-5-6, C-6-7, and other vertebrae for intraoperative X-rays. Harness can be pulled to either side by hook-and-loop strap at the chest to expose pelvis for bone graft. After X-rays are completed, harness can be cut at shoulders to release traction. The body portion of the harness is held by the patient's weight to help eliminate incidental trauma caused by having to strap the patient down. Adjustable hook-and-loop straps extend from the back of the body, up and over the shoulders, and then crisscross in front before attaching to the sides of the table just above the ankles. One size fits all.

Pad: 18½"W x 44"L • Straps: 66"L • Disposable
• Nonskid foam • Nonsterile • MR-safe

920877 CERVICAL VISUALIZATION HARNESS




Watch a video demonstration!
AliMed.com/cvh

Patient's weight secures body portion of harness



Single-Patient Use Armboard Straps




Soft, ultra-smooth nylon strap is extremely comfortable against skin

■ Ultra-smooth

Soft-Touch Armboard Straps Ultra-smooth woven straps are soft for maximum comfort and fragile skin protection in an economical, single-use strap. Minimal stretch ensures a firm hold while providing some conformity.

3"W or 1½"W x 32"L • Woven fabric • Hook-and-loop closure • Not made with natural rubber latex • White

3"W		
938162	12/cs	
93816204	48/cs	
93816212	144/cs	
1½"W		
938164	12/cs	
93816404	48/cs	
93816412	144/cs	




Resists stretch for secure hold with a soft, hook-sensitive exterior and padded interior

■ Stretch resistant

Soft-Secure Armboard Straps Padded straps feature a stretch-resistant core, for a secure hold. Soft, hook-sensitive outer surface laminated to a padded foam interior is gentle against skin.

3"W or 1½"W x 32"L • Hook-sensitive outer/foam-padded inner • Hook-and-loop closure • Not made with natural rubber latex • White/Blue

3"W		
938158	12/cs	
93815804	48/cs	
93815812	144/cs	
1½"W		
938160	12/cs	
93816004	48/cs	
93816012	144/cs	



■ Strong closure

Covidien Devon™ Single-Use Armboard Straps Ideal for use on O.R. tables and gurneys. Strong construction ensures secure and comfortable patient positioning. 48/cs.

1½"W x 32"L • Hook-and-loop closure • Not made with natural rubber latex • White/Blue

938132 SINGLE-USE ARMBOARD STRAP, 48/cs



Shop full line of Armboard Pads

AliMed.com/armboards



AliStrap® Precut Positioning Straps

Save O.R. and Imaging prep time

No cutting, no measuring, no hassles

Prepping your patients for O.R. or Imaging procedures can be a lengthy process, so having to waste valuable time searching for scissors or measuring tape to cut just the right size strap can be an added frustration.

Convenient precut lengths

AliMed's best-selling AliStrap is now available in convenient precut lengths designed to speed up your prep time while minimizing disruption to your sterile environment.

All the features of AliStrap

With a strong, self-adhering, back-to-back hook-and-loop closure and skin-safe, soft edges, these Precut Positioning Straps are safe enough to use against patient's skin yet strong enough to secure accessories—all without leaving behind sticky, bacteria-trapping residue. Perfect for single-use applications.

- ✓ No cutting or measuring required
- ✓ Individually wrapped and ready to use
- ✓ Skin-safe, soft edges, strong closure

Great for plastic surgery



More AliStrap options pp. 2, 3



See our full line of Positioning Products online



AliMed.com/positioners



Secure patients, equipment, cords, and more



Fully padded edges are gentle on fragile skin

■ No cutting, no measuring, no hassles



AliStrap® Precut Positioning Straps are available in original AliStrap or fully padded AliStrap Soft for fragile skin. Straps are individually wrapped for greater infection control and ready to use right out of the box to quickly and easily secure arms to armboards, IV lines and catheters or hot/cold therapy devices to extremities, or even accessories to your O.R. table. Sold in boxes of 12 or cases of 120.

Radiolucent • Metal-free • MR-safe

	ITEM #	ITEM #
ALISTRAP PRECUT	12/BX	120/CS
1½"W x 30"L	926725	939163
1½"W x 38"L	939182	939183
3"W x 30"L	926726	939164
ALISTRAP SOFT PRECUT	12/BX	120/CS
1½"W x 30"L	926727	939165
1½"W x 38"L	939184	939185
2"W x 30"L	926728	939166
2"W x 38"L	939206	939207


Custom sizes and volume discounts available...CALL!

Custom sizes and volume discounts available...CALL!



Custom to *YOUR* specifications

- Length
- Width
- Materials
- Closure



Surgical and Imaging Straps




Simple ordering and fast turnaround
@ [AliMed.com/custom](https://www.AliMed.com/custom)



Positioning solutions that improve patient outcomes

For 50 years, AliMed has been an industry leader in developing and delivering innovative gel and foam positioners engineered to redistribute pressure, maximize stability, reduce shear, and improve imaging visualization.

For our full line of Patient Positioners, visit [AliMed.com/positioners](https://www.AliMed.com/positioners)



Single-Use Foam Positioners



Gel Positioners



Radiolucent Positioners



AliMed.com



800.225.2610