

Premium pressure relief at an economical price



Head and Neck



Upper Extremity



Lower Extremity



Body

See inside for full product listing ➤



100% Gel Positioners

Our new AliBlue Gel Positioners feature the same premium pressure-relieving properties as other leading gels on the market but are more economical for a cost-effective positioning option. Made of 100% polyurethane-based viscoelastic gel. 18-month warranty.



100% viscoelastic gel



Impervious to fluids



Easy to clean



Adult



Pediatric



Adult



Pediatric

A. Head Donuts

B. Horseshoes

C. Contoured Headrests




D. Chest Rolls



E. Contoured Chest Rolls



F. Ankle Stirrup Pads



G. Knee Crutch Pads




H. Heel Cup Positioners




I. Dome Positioners



J. Contoured Armboard Pads



K. Standard Armboard Pads




L. Mayo Stand Covers




M. Ulnar Nerve Protectors, Oval



N. Ulnar Forearm Protectors, Rectangle



O. Surface Overlays



P. Vac Pac Universal Covers

ALIBLUE POSITIONER, GEL HEAD AND NECK

A. HEAD POSITIONER DONUT		
ADULT, 8"OD x 3"ID x 1 1/4"H		#960880
PEDIATRIC, 5 1/2"OD x 2 1/4"ID x 1 1/4"H		#960881
B. HORSESHOE DONUT		
ADULT, 8"OD x 3"ID x 1 1/4"H		#960894
PEDIATRIC, 5 1/2"OD x 2 1/4"ID x 1 1/4"H		#960895
C. CONTOURED HEADREST		
7" x 7" x 3 1/4" (1 1/2" CONTOURED)		#960893

ALIBLUE POSITIONER, GEL BODY

D. CHEST ROLL		
PEDIATRIC, 12" x 3" x 3"		#960898
12 1/2" x 4 3/8" x 3 3/4"		#960897
16 1/2" x 4 3/8" x 3 3/4"		#960885
20 1/2" x 4 3/8" x 3 3/4"		#960896
12" x 6" x 5"		#960915
16" x 6" x 5"		#960884
20" x 6" x 5"		#960883
E. CONTOURED CHEST ROLL		
SMALL, 12" x 6 1/4" x 5"		#960916
MEDIUM, 16" x 6 1/4" x 5"		#960917
LARGE, 20" x 6 1/4" x 5"		#960918

ALIBLUE POSITIONER, GEL LOWER EXTREMITY


F. ANKLE STIRRUP PADS		
SMALL, 12 1/2" x 1 1/4" x 3/8"		#960886
G. KNEE CRUTCH PAD		
SMALL, 11" x 9" x 3/8"		#960905
MEDIUM, 17" x 11" x 1/4"		#960904
LARGE, 22" x 20" x 1/2"		#960910
H. HEEL CUP POSITIONER, 6" x 3 1/4" x 1 1/4"		
		#960907
I. DOME POSITIONER, 14" x 6" x 3"		
		#960882

ALIBLUE POSITIONER, GEL UPPER EXTREMITY

J. CONTOURED ARMBOARD PAD		
SMALL, 16" x 5 1/4" x 2 1/4"		#960900
MEDIUM, 20" x 5 3/4" x 2 1/4"		#960899
LARGE, 24" x 5 3/4" x 2 1/4"		#960920
K. STANDARD ARMBOARD PAD		
SHORT, 15" x 4" x 1/4"		#960903
LONG, 24" x 4" x 3/8"		#960901
WIDE LONG, 24" x 6 1/2" x 3/8"		#960902
L. MAYO STAND COVER, 21" x 18 1/2" x 5/8"		
		#960911
M. ULNAR NERVE PROTECTOR		
OVAL, 15 1/8" x 6" x 1/2"		#960913
N. ULNAR FOREARM PROTECTOR		
RECTANGLE, 18" x 13" x 1/4"		#960914


ALIBLUE POSITIONER, GEL OVERLAYS

O. SURFACE OVERLAY		
SMALL, 20" x 20" x 1/2"		#960889
MEDIUM, 46" x 20" x 1/2"		#960879
FULL LENGTH, 72" x 20" x 1/2"		#960887
TORSO SECTION		
w/PERINEAL CUTOUT, 39 1/2" x 20" x 1/2"		#960888
FOOT SECTION		
w/PERINEAL CUTOUT, 20" x 20" x 1/2"		#960921
HEAD SECTION, 9" x 20" x 1/2"		
		#960890
P. VAC PAC UNIVERSAL COVER, 37" x 28" x 1/2"		
		#960912



AliBlue™ Gel vs. other leading gels

In a side-by-side comparison illustrated by the pressure mapping images below, AliBlue Gel performed as well as or better than the leading premium and economical gels in relieving pressure on the head, shoulders, pelvis, and heels. Test subject: Male, 6'0", 185 lbs.



Complete your reusable patient positioning setup

FULLY SUBMERSIBLE—THE ONLY HIGH LEVEL DISINFECTABLE STRAP!


VeriClean™ O.R. Table Straps with Mid-Panel The first surgical table strap with ultrasonically welded seams—no punctures, no holes, no stitching—allowing for complete and repeated submersion in high-level central processing disinfectant and reducing the risk of cross-contamination.

2"W thermoplastic urethane belt • Mid-panel: 5"W x 18"L • Cam-style metal buckle • Not made with natural rubber latex • 3-year warranty

#937065 MED., 96"L #937067 LG., 120"L
#937066 MED./LG., 108"L #937068 XL, 132"L



Ultrasonically welded seams for better infection control



To learn more and watch the video, visit AliMed.com/vericlean



Shop the full line at AliMed.com/surgical-straps-ecatalog



Any order placed before 3 p.m. EST, if in stock, ships the same day!
Prices, product availability, and specs are subject to change.

Pub. #7939REV0222 Product #11594
©2022 AliMed® Inc. All rights reserved.