

All the gear for full radiation protection

**LIMITED TIME—
SAVE OVER \$250!***

Pass on the savings to your NEW residents! Each Starter Kit contains all the personal radiation protection your new residents need to stay safe. Plus, no sharing means no cross-contamination!

Apron — Ultralight lead-free weight reliever apron with shrug-and-cinch design evenly distributes weight for less neck and back stress. Full 0.5 mm Pb-equivalent front protection available in four sizes.

Thyroid Shield — Lead-free with no binding to prevent chafing. Full 0.5 mm Pb-equivalent protection. One size. Color matches apron.

Glasses — Ultralight wraparound design for added side protection. 0.75 mm Pb-equivalent protection. Fits medium to large faces. Black.

SPECIAL PRICE!

\$525

~~REG. PRICE
\$786.50~~



ORDERING AND APRON SIZING CHART

KIT #	SIZE	APRON (W x L)	DRESS SIZE	SUIT SIZE	MAX. CHEST/HIP	MAX. HT.
#937202STK	S	20" x 32"	4-6	34"-36"	36"	5'6"
#937203STK	M	22" x 34"	8-10	38"-40"	40"	5'10"
#937204STK	L	24" x 36"	12-14	42"-44"	44"	6'0"
#937205STK	XL	26" x 38"	16-18	46"-48"	48"	6'2"

**Terms: Offer valid thru Dec. 31, 2022. Kit # includes Apron, Thyroid Shield, and Glasses. No additional discounts or substitutions. Choose contract price or Starter Kit promotion, whichever is lower. Sizes set by color, not interchangeable for this promotion. Monogramming not available. Kit ships same day when ordered by 3 p.m. EST.*

65 California residents, see below.



Gloves sold separately, [click here for details](#)

Order at [AliMed.com/starter-kit](https://www.ali-med.com/starter-kit)

California Proposition 65 Statement
65 **WARNING:** This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.

©2022 AliMed® Inc. All rights reserved. • Pub. #7943REV0222-A
 Prices, product availability, and specs are subject to change.