

Conforms to the patient's body for secure, stable positioning

✓ **Vacuum suction locks small beads into place**

✓ **Lightweight and flexible for easy positioning**

✓ **Stain-resistant and easy to clean**



Bean Bag with Cutout



Medium Universal Bean Bag



Small Universal Bean Bag

Bean Bag Positioners use vacuum suction to conform small, loosely packed beads to the patient's body, creating a solid foundation for safe, stable positioning. Lightweight and easy to position with superior flexibility, precise shaping, and maximum conformity. Stain-resistant, easy to clean, and nonporous. Available in two styles: **Universal Bean Bags** offer versatile positioning; **Bean Bags with Cutout** offer a space for extremities to lie in a more natural position. Select styles available **with Gel** for protection against shear. All styles require **Bean Bag Valve** (sold separately) and in-wall vacuum unit (not included) to apply suction—see "How to Use" section at right. Not made with natural rubber latex.



Available with gel top for additional shear protection

Small: 10"W x 20"L • **Medium:** 20"W x 20"L • **Large:** 30"W x 40"L • **Full Length:** 30"W x 76"L • **X-Large:** 36"W x 40"L • **Bariatric:** 46"W x 76"L

UNIVERSAL BEAN BAG			
	SIZE	WEIGHT	QTY
#2970014084	SMALL	2 lbs	EACH
#2970014085	MEDIUM	4 lbs	EACH

(MD) • 65 California residents, see below.

FOR ALL BEAN BAGS		
#2970013600	BEAN BAG VALVE	EACH
(MD) • 65	California residents, see below.	

BEAN BAG w/CUTOUT			
	SIZE	WEIGHT	QTY
#2970014086	LARGE	6 lbs	EACH
#2970012695	LARGE w/GEL	17 lbs	EACH
#2970014935	FULL LENGTH	10 lbs	EACH
#2970014087	X-LARGE BEST SELLER!	6.5 lbs	EACH
#2970013008	X-LARGE w/GEL	23 lbs	EACH
#2970014088	BARIATRIC	8 lbs	EACH
#2970012696	BARIATRIC w/GEL	31 lbs	EACH

(MD) • 65 California residents, see below.

How to Use



Bean Bag Valve

Form the Bean Bag around the patient in the desired position. Then, connect the Bean Bag Valve to your in-wall vacuum unit and apply suction to lock the Bean Bag in place. Unseal the vacuum valve to let air in and allow the beads move freely again for adjustment or to prepare for the next application.

California Proposition 65 Statement

65 ⚠ **WARNING:** This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.

Pub. #8140

©2023 AliMed® Inc. All rights reserved.