



AliMed®

Helping people...help people™ ...for 38 years

Effective and economical

BASIC PULL-CORD ALARM

One of the most economical ways to signal staff that a patient is about to wander. Monitor a patient who is seated or lying down. Easy to mount on beds and chairs with its metal clip or the included hook-and-loop fastener. Clip the activation cord to the patient's clothing. If the patient starts to rise, the pin connector pulls out of the unit, activating the alarm.

#78125 Basic Pull-Cord Alarm

ACCESSORIES (SOLD SEPARATELY)

Repl. Mounting Holster w/Strap, 2/pk #75424

- Mounts easily to beds and chairs
- Screws and nylon strap included

Repl. Pull-Cords, 6/pk #75382

- Cords easily adjust from 18" to 36"

Universal Alarm Mount Clip #77936

- Clips on bed, wheelchair, or IV stand
- Holds almost all alarms
- Spring-loaded clamp

Anti-Tamper Clips, 6/pk #74961

- Help prevent patient from releasing alarm clip
- Lock to deter removal
- Two-step release reduces tampering

Repl. 9V Batteries, 8/pk #71088



KEY BENEFITS

- ▶ Most economical wander-prevention alarm
- ▶ Use to monitor patients in beds or chairs
- ▶ Mount unit with metal clip or hook-and-loop fastener
- ▶ Cord adjusts from 18" to 36"
- ▶ Reinsert pin connector or use on/off switch to deactivate
- ▶ Two volume settings



297 High Street, Dedham, MA 02026 USA
www.AliMed.com
E-mail: info@alimed.com

Order Toll Free 800.225.2610
24-Hour Fax 800.437.2966
General Info 781.329.2900