

## Optimize cervical visualization for sedated patients



**Patient's weight secures  
body portion of harness**

- **Reduce risk of incidental trauma**
- **Use during intraoperative X-rays**
- **Single use for infection control**

### Cervical Visualization Harnesses

Retract shoulders for better visualization of C5–6, C6–7, and other vertebrae for intraoperative X-rays. Harness can be pulled to either side by hook-and-loop strap at the chest to expose pelvis for bone graft. The body portion of the harness is held by the patient's weight to help eliminate incidental trauma caused by having to strap the patient down.

Adjustable hook-and-loop straps extend from the back of the body, up and over the shoulders, and then crisscross in front before attaching to the sides of the table just above the ankles. One size fits all.

*Pad: 18½"W x 44"L • Straps: 66"L • Disposable • Nonskid foam • Nonsterile • MR-safe*

### Ordering Information

#92 0877

## Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

### The AliMed Advantage:

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items\*
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

### Operating Room/Medical Imaging/Critical Care product lines:

Anesthesia/Respiratory	Interventional Radiology	PACS Accessories	Surgical Instruments
Bed Positioning and Safety	Lighting	Positioners	Surgical Positioning
Clinical Furnishings	MRI	Radiation Protection	Surgical Table Accessories
Diagnostics	Medical Carts	Safe Patient Handling	Ultrasound
Ergonomics	MediChoice®	Storage and Organization	Wound Care
Imaging Accessories	Operating Room Safety	Straps	
Infection Control	Orthopedics	Support Surfaces	

### Rehabilitation/Orthopedics/Skilled Nursing product lines:

Activities of Daily Living	Clinical Lighting	Lower Extremities	Speech and Cognition
Bathing and Toileting	Contracture Management	Materials and Tools	Spinal Bracing
Beds and Accessories	Ergonomics	Mobility Aids	Therapeutic Exercise
Carts and Storage	Fall Management	Modalities	Upper Extremities
Clinical Assessment	Fracture Bracing	Pediatric Orthoses	Wheelchair Cushions
Clinical Diagnostics	Infection Control	Safe Patient Handling	Wheelchairs and Positioners
Clinical Furniture	Lighting	Skin Protection	Wound Care and Edema

\*Any order placed before 3 p.m. EST, if in stock, ships the same day!

Pub. #7519 Product #11385

©2017 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com

Call: 888.625.4633

Fax: 800.437.2966