

Sounds when the patient rises

- ✓ Early warning of patient movement
- ✓ Bed and chair sensor pads are cushioned to increase patient comfort
- ✓ Strong adhesive strips keep pad in place for fewer false alarms

Standard Bed and Chair Sensor Pad Monitor Systems

Give staff early warning of patient movement, helping prevent falls from beds and chairs. As soon as the patient moves off the bed or chair pad, the monitor triggers, giving staff time to respond and keep the patient safe.

AliMed Sensor Pads are water-resistant, and thinner for discreet comfort. Chair Sensor Pad is 10" x 15", Bed Sensor Pad is 10" x 30".

Customize pad placement for specific patient needs. 6-month guarantee. (Also available in 30-day guarantee for bed and chair sensor pads, and 18 month for chair sensor pads. Systems available in 30-day, 6-month, and 18-month guarantees.)

Strong adhesive strips keep the large sensor pads securely in place, allowing patients to shift positions without setting off the monitor. Fewer distractions with fewer false alarms.

Systems are available with AliMed alarms, including our TR2 Alarms featuring "Always On" technology that resists patient tampering.




6-Month Systems, Replacement Pads, and Accessories



TR2 Alarm with:

#77070 Bed Sensor Pad System
#77069 Chair Sensor Pad System



#79366 Bed Sensor Pads
#79367 Bed Sensor Pads, 10/cs
#77061 Chair Sensor Pads
#77063 Chair Sensor Pads, 10/cs



Partners in the Continuum of Care

AliMed is a leading developer, manufacturer, and distributor of high-quality patient care, rehabilitation, orthopedic, operating room, diagnostic imaging, and ergonomic products. For over 40 years, AliMed has served the hospital, long-term care, and rehabilitation markets with the broadest selections of specialty products.

Through established partnerships with more than 500 suppliers, a reputation for exceptional customer service, and a comprehensive range of unique products, AliMed distributes more than 70,000 different products to healthcare facilities and businesses worldwide.

AliMed provides the following:

- One of the largest portfolios of products with over 70,000 line items
- Service to all medical markets with an extensive array of products including those that are difficult to find
- Competitive pricing through contracts
- Field sales and dedicated customer service support
- Same-day shipping for in-stock items*
- Reduced acquisition costs through EDI ordering
- E-commerce and specialized catalogs to easily find and order products through

* if ordered before 3PM EST.

Offering a wide selection of products:**

Rehabilitation/Orthopedics/Skilled Nursing product lines:

Activities of Daily Living	Materials and Tools	Skin Protection	Wheelchairs, Cushions, and Positioners
Fall Management and Patient Safety	Mobility Aids	Speech and Cognition	Workstation Ergonomics
Contracture Management	Orthoses and Orthotics	Splints and Braces	Wound Care
	Patient Transfer	Therapeutic Exercise	

Operating Room/Diagnostic Imaging/Critical Care product lines:

Clinical Furnishings	Operating Room Safety	Patient Handling/Safety	Straps
Diagnostics	Organization and Storage	Positioners	Surgical Positioning
Imaging Accessories	PACS Accessories	Radiation Protection	Surgical Table Accessories
Infection Control			

****Visit AliMed.com for our complete product selection!**

297 High Street, Dedham, MA 02026 USA • AliMed.com

Call: 800.225.2610

Fax: 800.437.2966

Email: info@AliMed.com