



AliMed®

"Helping people help people"

Sounds only where it's needed

ALIMED® WIRELESS BED AND CHAIR SYSTEMS (WITH BOTH LOCAL AND REMOTE ALARMS)

AliMed's Wireless Bed and Chair Sensor Pad Systems utilize new thin-profile, high-reliability sensor pads connected to the local bed/chair-side Alarm/Transmitter Unit (#2130). When a patient attempts to rise, the alarm sounds and a signal is transmitted to the Remote Receiver/Alarm Unit (#2131) located up to 100' away. The receiver sounds an alarm and flashes a light. Transmitter and receiver can be coded to pairs to allow multi-unit installation.

#2144 Wireless Chair Sensor Pad Alarm System

#2143 Wireless Bed Sensor Pad Alarm System

ACCESSORIES (SOLD SEPARATELY)

Repl. Chair Sensor Pads	#79210
Repl. Bed Sensor Pads	#79366
Repl. Alarm/Transmitter Unit	#2130
Repl. Receiver/Alarm Unit	#2131
Repl. C Batteries, 4/pk	#5299



KEY BENEFITS

- ▶ Provides early warning of patient's attempts to get out of bed or chair
- ▶ Alarm/Transmitter Unit (#2130) can alarm, chime, or remain silent while transmitting to Remote Receiver/Alarm Unit (#2131) up to 100' away
- ▶ Bright red locator light on Alarm/Transmitter Unit (#2130)
- ▶ Multiple units can be used together yet controlled separately—ideal for hospitals, long-term care facilities, and nursing homes
- ▶ One-year warranty



297 High Street, Dedham, MA 02026 USA
www.AliMed.com
E-mail: info@alimed.com

Order Toll Free 800.225.2610
24-Hour Fax 800.437.2966
General Info 781.329.2900