



AliMed®

Helping people...help people™...for over 40 years

Protect the diabetic and arthritic foot

Why diabetic insoles?

Avoid ulceration of the diabetic foot—brought on by an unyielding shoe insole, inadequate arch support, or a poorly fitted shoe.

To keep feet properly supported and protected, it's important to use a pair of insoles that offer total contact, support, and comfort.

Insoles like the **Duo/Laminate D-Sole™**, **FREEDOM® Trilaminat 1 Insole**, and **FREEDOM® Full-Contact Insole** offer maximum cushioning and superior pressure distribution for the entire plantar surface by utilizing Plastazote®, which has gained a reputation for beneficial use with the neuropathic foot and in managing the arthritic foot. An incorporated heel cup and longitudinal arch combine for stabilization, support, and a comfortable fit. Heat-mold for a customized fit. For additional relief to high-pressure areas, add accommodations such as scaphoid or met pads. Protect feet and save time by keeping our diabetic insoles inside each pair of shoes.

Sizing and ordering information on back.

Duo/Laminate D-Soles™

- ▶ Plastazote layer offers total contact and maximum pressure distribution



- ▶ Urethane base prevents bottoming out
- ▶ Tapered forefoot reduces toe box cramping

FREEDOM® Trilaminat 1 Insoles

- ▶ Perfect for the insensitve/hypersensitive foot
- ▶ Carefully selected materials ensure shape retention
- ▶ Heat moldable
- ▶ Medicare approved—HCPCS Code A5512



FREEDOM® Full-Contact Insoles

- ▶ Feature a deeply cupped, stabilizing heel and moderate longitudinal arch
- ▶ Resilient and conforming materials protect and accommodate insensitive feet
- ▶ Molds with heat and holds to give total contact for maximum pressure distribution




297 High Street, Dedham, MA 02026 USA
www.AliMed.com
email: info@AliMed.com

Order Toll Free 800.225.2610
24-Hour Fax 800.437.2966
General Info 781.329.2900

Duo/Laminate D-Soles™




Urethane thickness gradually tapers in the forefoot.



Size	Wmns	Mens
A	5-8	4-6
B	8½-10	6½-8
C	10½-11	8½-10
D	11½+	10½+


#KWW6227 ¼" Thick
 #KWW6226 ⅛" Thick
 Specify size.

FREEDOM® Trilaminate 1 Insoles



⅛" fine-celled ProCell™ EVA

⅛" supportive Poron®



⅜" ProCell™ EVA base

Size	Wmns	Mens
XS	5-7	—
S	8-9	6-8
M	10-11	8-9
L	—	10-11
XL	—	12-13
2XL	—	14-15

#KWW66376* Trilaminate 1 Insoles
 #KWW6100 Casting Pad
 #KWW70389 Convection Oven
 *Specify size.

FREEDOM® Full-Contact Insoles



Size	Wmns	Mens	pr	6 pr/cs
XS	5-7	—	#KWW64201	#KWW64211
S	8-9	6-7	#KWW64202	#KWW64212
M	10-11	8-9	#KWW64203	#KWW64213
L	—	10-11	#KWW64204	#KWW64214
XL	—	12-13	#KWW64205	#KWW64215
2XL	—	14-15	#KWW64206	#KWW64216



297 High Street, Dedham, MA 02026 USA
 www.AliMed.com
 email: info@AliMed.com

Order Toll Free 800.225.2610
 24-Hour Fax 800.437.2966
 General Info 781.329.2900