



## Gel Cushions



### T-Gel™ Cushion



**T-Gel Cushion** Excellent for applications requiring minimal thickness. T-Gel polymer gel “slides” laterally, effectively reducing harmful shear forces and skin-damaging friction. Low-profile design allows for use on any sitting surface.

- #1321 16"x16" smooth
- #1319 18"x16" smooth
- #1530 16"x16" checkerboard
- #1529 18"x16" checkerboard

### T-Gel Plus™ Cushion:



**T-Gel Plus Cushion** combines the benefits of T-Gel with those of slow-recovery T-Foam for even better pressure distribution and shear reduction. 1" thick medium-resistance T-Foam layer magnifies the positive effects of the gel at only half the weight of all gel cushions.

- #1292 18"x16" smooth
- #1531 18"x16" checkerboard

## T-Gel™




### KEY BENEFITS:

- ▶ Gel has similar consistency to human tissue—excellent pressure reduction function
- ▶ Best way to reduce risk of skin breakdown
- ▶ Delivers optimal shear reduction


### Lightweight Gel Composite Cushion

Raised circles on top layer are for lightweight shear and pressure protection. Two inch Open-Cell Urethane Foam with 3/4" lightweight gel polymer top layer. Washable nylon cover with a nonslip base.

- #10950 16" x 16" x 2 3/4"
- #10951 18" x 16" x 2 3/4"



## Aqueous Gel and Foam combination cushions



**Elite Cushions** for moderate to high-risk users. DuraGel™ inserts for relief at critical pressure points. Moderate slope from front to back with medial/lateral thigh positioners and lateral pelvic supports for optimal posture. AT model is an anti-thrust cushion with a pre-ischial bar for greater pelvic stability. Viscoelastic high-density for long term use. Polymer DuraGel inserts will not leak, stay in place, and are non-toxic. Includes black, fluid-proof stretch cover. Meets CAL 117. Two-year warranty.

- #1410 Elite
- #1411 Elite AT

Specify 16" x 16"; 18" x 16"; 18" x 18"; 20" x 18"; 20" x 20"; 22" x 18"; 24" x 18"; 22" x 20"; or 24" x 20".



**Econo-Gel Wheelchair Cushion** Two-part gel design reduces pressure and increases comfort. High-density foam with aqueous gel in middle. Safety ties. Wipe-clean vinyl cover. Optional polyester outer cover. Anti-slip bottom.

- Vinyl Only  
#75 1155 18" x 16" x 2"
- Vinyl w/Polyester Cover  
#8784 18" x 16" x 2"



**Gel-Foam Cushion** Bottom chamber contains resilient foam to prevent bottoming out. Aqueous gel-filled top chamber evenly distributes weight for pressure relief. Sealed with incontinence-proof vinyl. Washable polyester cloth cover ties in place. Meets CAL 117B. Replacement guaranteed for one full year. 2 1/2" H.

- #75 1010 16" x 16"



297 High Street, Dedham, MA 02026 USA  
www.AliMed.com  
E-mail: info@alimed.com

Order Toll Free 800.225.2610  
24-Hour Fax 800.437.2966  
General Info 781.329.2900