

Built for your demands

- **Organized, portable, and secure storage**
- **Durable construction**
- **Multiple styles and locking options available**

AliMed Medical Carts offer convenient, organized, and secure access to medications, supplies, and other treatment needs—anytime, anywhere—helping to meet the everyday demands of your department while maximizing staff efficiency for better patient outcomes. Choose from **Emergency, Anesthesia, Isolation, Treatment, or Procedure Carts**.



Standard Features on all AliMed Carts:

- Double-wall construction
- Lock styles for every security level
- Durable powder-coat paint
- Ball-bearing drawer slides with full-extension drawers
- Solid or two-tone frames with choice of color
- All-swivel casters



1 CART STYLES

Start by choosing your **Series** type—made from lightweight aluminum or durable steel.



STANDARD

Durable and affordable
Economically priced and built to last—guaranteed. Durable 12-gauge steel construction is designed to handle any situation for everyday storage you can count on.



SELECT

30% lighter—loaded with features

Aluminum construction is lighter than our Standard Series for effortless handling. Larger work surface plus movable mounting panels to load up on accessories.



WIDE

Over 40"W for maximum storage

Lightweight aluminum construction featuring an extra-wide frame and 30" of vertical drawer space for ample storage.



MINI TOWER

Standard functionality, compact size

All the functionality and features of our Standard Series in a compact, fully mobile size so you can easily maneuver through congested areas.



MINI CART

Space-saver model

Economical Mini Carts feature a space-saving, narrow design that easily fits under counters or at bedsides. Ideal for when mobility needs are limited.

See reverse for:

Lock Styles, Colors, Accessories, and Customization details

2 LOCK STYLES

Five locking options offer the security level that's right for your department. All carts (except Emergency models) ship with key override option unless otherwise specified.



Breakaway Locks
Low security, best for quick access



Key Locks
Low security, best for quick access



Push-Button Locks
Medium security, best for multi-person use



Electronic Locks
Medium to high security, best for monitored access



Proximity Locks
High security, best for monitored access

3 CART COLORS

Choose from five drawer colors available in matching solid or two-tone beige frame for easy recognition or departmentalization.


Solid Frame



Two-Tone Frame



Standard Colors



4 ACCESSORIES

Accessory Packages

SAVE when you order one of our pre-bundled **Emergency**, **Anesthesia**, or **Procedure** packages featuring our most popular, must-have accessories.



Individual Accessories

Choose just the accessories you need—50+ options available!



Custom to YOUR specifications

Have a unique need? Design your own cart with your choice of drawer configurations, expanded color choices, locking options, extended warranty, and accessories.

12 Custom Color Options



Simple ordering and fast turnaround @ AliMed.com/custom

Visit AliMed.com/mc for details.

Prices, product availability, and specs are subject to change.

Pub. #7680 Product #11513
©2019 AliMed® Inc. All rights reserved. • Printed in USA.



AliMed.com

Call: 800.225.2610

Fax: 800.437.2966