

Minimize cross-contamination risk with Hygienic radiation protection

- **Protect your thyroid while guarding against infection**
- **Full 0.5 mm Pb-equivalent protection**



Ultralight Washable

Washable up to 170 times in temps exceeding 140° F*



Ultralight Disposable

Single use eliminates cleaning and disinfection



Antimicrobial

Resists bacterial growth—ideal for sharing



Fluidproof

Impervious to liquids and antimicrobial

Total protection

The thyroid is one of the most radiation-sensitive organs in the head and neck. Wearing Hygienic Thyroid Shields during procedures helps to reduce harmful radiation doses by as much as 45%¹ while guarding against dangerous infection-causing bacteria.

Made from a proprietary lead-free alloy featuring full 0.5 mm Pb-equivalent protection, **AliMed Hygienic Thyroid Shields** are extremely lightweight, free of top binding to help prevent chafing, and available in your choice of infection protection—from **Ultralight Washable** or **Disposable** options to **Antimicrobial** or **Fluidproof** models.

¹American Thyroid Association, Policy Statement on Thyroid Shielding During Diagnostic Medical and Dental Radiology, 2013.

See reverse for full product details and order information. ➤

Visit AliMed.com/radiation-protection

See our full line of Radiation Protection Aprons, Eyewear, and Gloves.



*Data on file.

Hygienic Thyroid Shields



Ultralight Washable Thyroid Shield

- Ripstop nylon cover
- Machine-washable

#938186 Without Embroidery
#938187[†] With Embroidery

Specify color, see chart below.

[†]Not available for same-day shipping.



Ultralight Disposable Thyroid Shield

- Spunbound polypropylene cover
- Light Blue only

#937929

Monogramming not available.



Antimicrobial Thyroid Shield

- Linen-textured fabric
- Light Blue only

#934833

Monogramming not available.



Fluidproof Thyroid Shield

- Polyurethane fabric
- Antimicrobial

#939123

Monogramming not available.

Specify color, see chart below.

Visit AliMed.com/radiation-protection for more colors.



Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

The AliMed Advantage

- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts
- Field sales and dedicated customer service support
- Same-day shipping for in-stock items*
- Efficient ordering through EDI

*Any order placed before 3 p.m. EST, if in stock, ships the same day!
Prices, product availability, and specs are subject to change.

Pub. #7232REV0618 Product #11211
©2018 AliMed® Inc. All rights reserved. • Printed in USA.