

The only high-level, fully disinfectable strap



Strap with Mid-Panel shown

- Eliminates places where bacteria can hide
- No stitching reduces the risk of cross-contamination
- Fits all surgical tables with side rails

VeriClean™ O.R. Table Straps



Properly cleaning and disinfecting **conventional table straps** is challenging. Most rubber laminate O.R. straps are stitched together, creating punctures, perforations, and thread fibers that can be a breeding ground for bacteria. It is nearly impossible to disinfect those buried and hidden areas.

VeriClean is the first ultrasonically welded surgical strap. No stitching. No puncture holes. No hidden seams. Completely submersible in approved disinfectants*—every conceivable surface can be disinfected, leaving nowhere for bacteria to hide.

One less host for bacteria. In an independent lab study, VeriClean exceeded the FDA-recommended endpoint used to define high-level disinfection.¹ Conventional O.R. table straps fell short of achieving high-level disinfection.²

**Follow manufacturer's disinfection instructions.*

See reverse side for additional product details and ordering information.



¹ Data on file.

² High-level disinfection requires a 6.0 Log Reduction in microorganisms. Class II Special Controls Guidance Document: Medical Washers and Medical Washer-Disinfectors; Final Guidance for the Medical Device industry and FDA Review Staff.

VeriClean™ O.R. Table Straps

- **Welded construction eliminates stitching and reduces cross-contamination risk**
- **Multiple lengths available to accommodate patients of all sizes**
- **Fits all surgical tables with side rails**
- **3-year warranty**



Reduces the risk of cross-contamination



Withstands high-level disinfection

Order your VeriClean Strap

STRAP WITH MID-PANEL, 5"W x 18"L mid-panel

- #93 7065 Medium, 96"L
- #93 7066 Medium/Large, 108"L
- #93 7067 Large, 120"L
- #93 7068 X-Large, 132"L



Strap with Mid-Panel

BARIATRIC STRAP, 11"W X 30"L mid-panel

- #93 7678 140"L



Bariatric Strap

SINGLE STRAP, 2"W

- #93 7121 Medium, 96"L
- #93 7122 Medium/Large, 118"L
- #93 7123 X-Large, 132"L



Single Strap

All straps: 2"W thermoplastic urethane belt • Cam-style metal buckles • Not made with natural rubber latex

Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

The AliMed Advantage:

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items*
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

*Any order placed before 3 p.m. EST, if in stock, ships the same day!

Pub. #7199REV1117 Product #11193

©2017 AliMed® Inc. All rights reserved. • Printed in USA.