

The *secure* solution for Trendelenburg Positioning

Simple. Secure. Reusable.



■ **Cradles patient to maximize pressure redistribution**

■ **Unobstructed chest and IV access**

SecureFit[™] TPS

The secure, easy-to-use solution that dramatically reduces setup time, resists slipping, and evenly redistributes pressure throughout the patient's body.

Engineered to minimize slipping

The SecureFit TPS cradles the body's natural curvature, reducing peak focal pressure over the shoulder. No ancillary shoulder braces are required.

Maximum pressure redistribution

The **TPS Gel Pad** conforms to each patient's individual anatomy and redistributes pressure to provide support while minimizing slipping.

Simple, reusable, and easy to clean

With only a few key components and smooth, easy-to-clean surfaces, AliMed SecureFit TPS reduces setup time, which helps improve efficiency and lower costs.

IV access simplified

Access IV lines quickly and easily through convenient openings between arm support and gel pad.

To order, see reverse side.

SecureFit™ TPS

TPS Arm Supports and Toboggans

Integrated, easy-to-clean vinyl-covered foam pads work with toboggans to protect ulnar nerve while providing easy IV access to the arms.



Toboggan

TPS Gel Pad

Viscoelastic gel conforms to patient's anatomy and redistributes pressure for maximum support with minimal slipping.



- **Simple, secure, reusable**
- **Unobstructed chest and IV access**
- **Set up in less than a minute**
- **No shoulder braces required**

Specifications:

TPS Foundation: 20"W x 17"L x 6"H • **TPS Gel Pad:** 22½"W x 39"L x 1½" thick (1" gel surface with ½" Polyfoam base)
• Perineal cutout: 9"W x 13"L • **TPS Stabilizing Wedge Pads:** 1" Memory Foam over Polyfoam • **Cover:** Black anti-shear elastic • Antimicrobial • Fluid-proof fabric • Fluid-resistant sewn seams • CAL 117-compliant • **Stabilizing Posts:** 9¼"L x 5/8" diameter • Stainless Steel • **Toboggan:** 10"L x 8"H x 10¼"D • **Toboggan Arm Supports:** 25"L x 4½"H Polyfoam pads • Patent Pending

Ordering Information:

#93 8040 SecureFit TPS

Replacement Items:

- #93 8042 TPS Gel Pad
- #93 8043 TPS Stabilizing Wedge Pads, pair
- #93 8041 TPS Foundation
- #93 8046 Right Stabilizing Wedge
- #93 8047 Left Stabilizing Wedge
- #93 8044 Stabilizing Posts, pair
- #93 8045 Toboggans, pair

TPS Stabilizing Wedge Pads

Memory foam pads minimize the incidence of poor blood flow to the shoulder.

TPS Foundation

Strong, durable, nonskid foundation is contoured to cradle patient's natural curvature, and locks onto the surgical table with clamps for a secure hold.



Pub. #7253 Product #11218

©2015 AliMed® inc. All rights reserved. • Printed in U.S.A.