



 **AliMed[®]**

Your global source for healthcare products

Table of Contents

POSITIONERS	3
STRAPS	9
SUPPORT SURFACES	15
SURGICAL TABLE ACCESSORIES	16
SAFE PATIENT HANDLING	18
RADIATION PROTECTION	20

POSITIONERS

Nonslip Overlay

Nonslip Surface Overlay or Disc Keep positioners and patients in place with this innovative, autoclavable material—use with positioners to ensure patient security at any angle during procedures, reducing the need for tape or other adhesives. Reusable and antimicrobial, the dual-sided silicone is ultrathin and flexible to conform to any table's support surface. Surface Overlay is easy to clean and remove. Non-adhesive material is safe against skin and won't leave residue on surfaces that could collect bacteria. Trim to any size and works even when wet. Both the Overlay and Disc stabilize gel positioners, coated foam, and vinyl-covered positioners from shifting.



Overlay: 20"W x 46½"L (cut to size) • 1-mm thick • Pre-cut Disc: 9" diam. • 1.5-mm thick • Autoclavable

[938108](#) Overlay
[938109](#) Disc

Bean Bag Positioners

Bean Bag Positioners use vacuum suction to sculpt a solid foundational support. Unseal the vacuum valve to let air in, and the beads move freely again. Lightweight and easy to position. AliMed's urethane material works together with the beads for the right flexibility, precise shaping, and maximum conformity. Stain-resistant, easy to clean, and nonporous. Available in two styles: **Bean Bags with Cutout** and **Universal** (without cutout). Bean Bags with Cutout offer a space for extremities to lie in a more natural position. Universal Bean Bags can be used in a variety of ways for quick and easy positioning. Not made with natural rubber latex.



Universal		w/Cutout	
935231	Small, 10"W x 20"L	935233	Lg. w/Shoulder Cutout, 30"W x 40"L
935232	Med., 20"W x 20"L	935234	X-Lg. w/Shoulder Cutout, 36"W x 40"L
		937565	Full Length w/Shoulder Cutout, 30"W x 76"L
932644	Valve Repair	935235	Bariatric w/Shoulder Cutout, 46"W x 46"L

SecureFit TPS (Trendelenburg Positioning System)

SecureFit® TPS provides a stable, easy-to-use solution that dramatically reduces setup time, resists slipping, and evenly redistributes pressure. The contoured design cradles the scapula region, avoiding the brachial plexus, reducing peak focal pressure over the shoulder. SecureFit TPS fits all standard-width tables. Best suited for Robotic, Laparoscopic, and Open procedures where Trendelenburg is used. Table clamps not included.



TPS Foundation: 20"W x 17"L x 6"H • Polycarbonate construction • TPS Gel Pad: 22½"W x 39"L x 1½" thick (1" gel surface with ½" Polyfoam base) • Perineal cutout: 9"W x 13"L • TPS Stabilizing Wedge Pads: 1" Memory Foam over Polyfoam • Cover: Black antishair elastic • Antimicrobial • Fluid-proof fabric • Fluid-resistant sewn seams • CAL 117-compliant • Stabilizing Posts: 9¼"L x ¾" diameter • Stainless Steel • Toboggan: 10"L x 8"H x 10 ¼"D • Toboggan Arm Supports: 25"L x 4½"H Polyfoam pads • Patent Pending

[938040](#) SecureFit TPS

A-Line Support

A-Line Support improves access to the radial artery by positioning the wrist just right—every time. It allows skin to be stretched taut with artery presented optimally and thumb abducted—all of which facilitate cannulation. Deluxe model includes a skin-friendly, 1"W padded strap to minimize movement, reduce ulceration risk and secure the hand and wrist to the positioner. Use your own straps with our Standard version.



3"W x 10 ½"L • Closed-cell foam

[95-984](#) Standard A-Line Support, 30/cs

[95-985](#) Deluxe A-Line Support, 12/cs

Arthroscopic Leg Holders – Single-Patient-Use

Single-Patient-Use Well Leg Holders cradle the nonaffected leg to help eliminate discomfort and circulatory impairment that can result when the leg hangs over the edge of the table.



Single-Patient-Use Operative Leg Holders Positioning aid for operative and well leg during Arthroscopy. Knee holder secures and stabilizes operative leg in the Arthroscopic Leg Holder. 5/cs.

Well Leg: 14"W x 17"L with 6½"W x 2" cutout and 90° angle • Polyfoam • White •
Operative: 3"W x 11"OD x 7"ID • White

[931959](#) Well Leg Holder, Uncovered, each

[931961](#) Operative Leg Holder, 5/cs

Arthroscopic Leg Holders – Vinyl-Coated

Vinyl-Coated Arthroscopic Leg Holders are fluid-proof, seamless positioners that are economical and easy to clean. **Well-Leg Holder** protects against nerve damage to non-operative leg. Available in two sizes. **Operative Holder** has locator grooves help keep it centered in clamp. Not made with natural latex rubber.



Vinyl-coated polyethylene • **Well-Leg Short:** 8"W x 4"H x 4"D • **Well-Leg Long:** 8"W x 8"H x 8"D •
Long offers additional coverage along thigh and lower leg • **Operative Holder:** 3"W x 11"OD x 7"ID

[931938](#) Vinyl-Coated Well Leg Holder, Long

[931939](#) Vinyl-Coated Well Leg Holder, Short

[931940](#) Vinyl-Coated Operative Leg Holder



Arthroscopic Leg Holders – Vinyl-Covered

Upholstered Well Leg Holders cradle the nonaffected leg to help eliminate discomfort and circulatory impairment that can result when the leg hangs over the edge of the table. Upholstered version is Polyfoam covered with black conductive vinyl cover. Simply wipe clean and store for reuse. Reversible—either side can sit on top edge of table.



Polyfoam • 14"W x 17"L with 6½"W x 2" cutout and 90° angle

[931960](#) Well Leg Holder

Feldon Hand Prepper

Feldon Hand Prepper™ designed by Paul Feldon, M.D. This very firm, concave block slips under the arm of the prone or supine patient. 5" elevation on a stable base.

6"W x 5"H x 6"D • Water-resistant seams • Black conductive vinyl cover

[95-054](#) Feldon Hand Prepper



Foot/Hand Preppers

Foot/Hand Preppers Large base to ensure stability. Concave shape holds arm or leg securely. Slanted side supports shoulder. Bacteriostatic black conductive cover with water-resistant seams.

12"W x 15"L x 7"H

[95-055](#) Foot/Hand Prepper



Adjustable Head Donut

Adjustable Donuts Customize the size of center cutout to meet a range of positioning needs. Size and shape vary in diameter and surface area. Select from a solid surface or remove inserts to widen cutout in center of donut. 36/cs.

9"W x 10"L x 1"H (Includes 3 cutouts: 2 ¼"W x 2 ¼"L; 3 ¼"W x 4"L; 3 ¾" x 5"L) • Polyfoam

[95-549](#) Adjustable Donut



Foam Donut

Reusable Head Donuts offload pressure from occipital protrusion to protect the patient from pressure-related issues, such as pressure ulcers and possible blindness. Protect the head, neck, and ears. Recommended for supine, lateral, and lithotomy. Also use to protect knees in prone. Conductive Black Vinyl cover is durable, and easily wipes clean.

S: 1 ½"H x 6" diam. L: 1 ½"H x 9" diam.

[930484](#) Foam Donut, Small

[95-681](#)

Foam Donut, Large



Positioning Wedges

Vinyl-Covered Positioners Designed specifically for radiology and the O.R., our Radiolucent Blue Vinyl and Conductive Black Vinyl covered foam positioners are easy to clean with water-resistant seams.


Flame resistant CAL 117-complaint • Antibacterial vinyl cover

[927798](#) 45 Degree Wedge, Blue Vinyl, 7"W x 30"L x 7"H

[931711](#) 45 Degree Wedge, Black Vinyl, 7"W x 24"L x 7"H

[931716](#) 27 Degrees Wedge, Black Vinyl, 10"W x 10"L x 6"H

[930774](#) High Perineal Wedge, Black Vinyl, 18"W x 20"L x 9"H



AliLite Gel Positioners

AliLite™ Precontoured Gel Positioners reduce pressure by distributing weight over a larger surface area. More than 50% lighter than comparable gel-only products, AliLite™ gel-foam positioners are easier to handle, impervious to fluids, naturally bacteriostatic, and feature a stable, nonskid base. AliLite™ gel is made of dry viscoelastic polymers.



Body Positioners

- [931872](#) Table Pad, 3/4 Length
- [937880](#) Table Pad, Full Length
- [937881](#) High Chest Roll, 6"W x 16"L x 6"H
- [937882](#) High Chest Roll, 6"W x 20"L x 6"H
- [931559](#) Chest Roll, 6"W x 12"L x 4"H
- [931558](#) Chest Roll, 6"W x 16"L x 4"H
- [931557](#) Chest Roll, 6"W x 20"L x 4"H
- [931871](#) Dome Positioner

Upper Extremity

- [937885](#) Contoured Armboard Pad, Short
- [931560](#) Contoured Armboard Pad, Long
- [931551](#) Armboard Pad, Short
- [931552](#) Armboard Pad, Long

Lower Extremity

- [937883](#) Lateral Leg
- [931875](#) Lateral Positioner
- [937886](#) Prone Positioner, Torso
- [937887](#) Prone Positioner, Full Body
- [932990](#) Patient Prone Positioning System
- [937888](#) Knee Supports, Large
- [937889](#) Knee Supports, Small
- [931873](#) Heel Pads, pair

Head

- [931553](#) Horseshoe, Adult
- [931554](#) Horseshoe, Pediatric
- [931555](#) Donut, Adult
- [931556](#) Donut, Pediatric
- [931696](#) Prone Headrest, Adult
- [931697](#) Prone Headrest, Pediatric
- [931874](#) Supine Headrest, Adult
- [937884](#) Supine Headrest, Pediatric

Positioning Sets

- [937921](#) Prone Positioning Set
- [937922](#) Lateral Positioning Set
- [937924](#) Cardiovascular Positioning Set
- [937925](#) Supine Positioning Set



AliSoft Head Positioners

AliSoft Head Positioners offload pressure on the head, neck, and spine during prone procedures. 2" soft Polyfoam top layer for extra skin protection and higher-density bottom layer for stability. Includes two cutout slots for ventilation tubing. Vacuum-sealed to save storage space. 12/cs.

10 3/4"W x 9"L • Polyfoam

[925933](#) Adult, 7"H

[925934](#) Pediatric, 6"H



Disposable Foam Heel Pads

Disposable Foam Heel Pads offer pressure reduction to ankle, heel, and foot. Secure with hook-and-loop attachment. Universal sizing. Available in standard packaging or in a compressed, vacuum-sealed version to reduce storage space. 12 pr/bx.

8 1/2"W x 8 1/2"L x 1 1/4"H • Individually wrapped • Convoluted foam

[937136](#) Foam Heel Pads

[937136EU](#) Foam Heel Pads, Vacuum-sealed



Disposable Arm Cradles

Disposable Arm Cradles are ideal for comfortable positioning of the full length of the patient's arms and hands. Can be used in either supine or prone positions. Cutouts and slices cradle the arm's natural bent position to help prevent tissue breakdown and nerve damage. Available in standard packaging or in a compressed, vacuum-sealed version. 12 pairs - ships in 2 cases of 6 pairs (total of 24 pieces).

5"W x 24"L x 3"H • Standard polyfoam

[930377](#) Arm Cradle

[930377EU](#) Arm Cradles, Vacuum-sealed



Jackson Table Positioning Set

Jackson Table Positioning Sets include all the necessary positioning components in one convenient kit, saving O.R. setup and cleanup time. Made of economical Polyfoam, these single-patient-use sets help reduce the risk of cross-contamination between procedures while ensuring patients are properly and consistently positioned during any surgical procedure using specialty O.R. tables. Vacuum-sealed to reduce storage space. Not made with natural rubber latex. 6/cs.


Kit contains: (A) arm cradles (pair), (B) chest pillow, (C) hip and thigh pad covers (4).

[925936](#) Jackson Table Positioning Set



Wilson Frame Positioning Set

Wilson Frame Positioning Sets include all the necessary positioning components in one convenient kit, saving O.R. setup and cleanup time. Made of economical Polyfoam, these single-patient-use sets help reduce the risk of cross-contamination between procedures while ensuring patients are properly and consistently positioned during any surgical procedure using specialty O.R. tables. Vacuum-sealed to reduce storage space. Not made with natural rubber latex. 6/cs.



Components: (A) drape protector, (B) pad sleeve covers, (C) bar protector, (D) AliSoft head positioner, (E) arm cradles (pair). Available in three sets:

[925938](#) Basic (A-C Only)

[925937](#) Economy (A-D Only)

[925939](#) Standard (All)

Andrews Table Positioning Set

Andrews Table Positioning Sets include all the necessary positioning components in one convenient kit, saving O.R. setup and cleanup time. Made of economical Polyfoam, these single-patient-use sets help reduce the risk of cross-contamination between procedures while ensuring patients are properly and consistently positioned during any surgical procedure using specialty O.R. tables. Vacuum-sealed to reduce storage space. Not made with natural rubber latex. 6/cs.



Components: (A) AliSoft head positioner, (B) knee pads, (C) seat cover, (D) chest cover, (E) thigh strap, (F) arm cradles (pair), (G) heel protectors (pair).

[925940](#) Andrews Table Positioning Set

STRAPS

AliStrap®

AliStrap and AliStrap Soft are all-purpose positioning straps with a self-adhering, back-to-back hook-and-loop closure that won't leave behind sticky, bacteria-trapping residue—making them ideal for securing, immobilizing, or positioning patients, equipment, and more. Low-profile, medical-grade hooks and non-abrasive soft edges will not cut into or abrade skin, instead yielding and flexing against soft tissue to reduce subcutaneous shear. AliStrap Soft's padded-loop center provides extra protection for at-risk skin. Economical and disposable, these straps are perfect for single use, preventing cross-contamination and offering easy cleanup. Simply cut any length from the continuous 30'L roll, wrap, and press together to form a tight, strong bond—even holds when wet.



AliStrap, 30'L Roll

- [95-905](#) ¾"W Narrow
- [95-900](#) 1 ½"W Regular
- [95-902](#) 3"W Wide
- [95-904](#) 5"W X-Wide
- [95-908](#) 10"W Maxx
- [95-906](#) Combo Pack, One each of ¾", 1 ½", 3"W, and 5"W

AliStrap Soft, 30'L Roll

- [930382](#) 1 ½"W Padded
- [930163](#) 2"W Padded



AliStrap® Precut Positioning Straps

AliStrap® Precut Positioning Straps Our best-selling AliStrap and fully padded AliStrap Soft for fragile skin are now available in convenient precut lengths for Imaging and Surgery, saving valuable prep time searching for scissors or measuring tape while minimizing disruption to your sterile environment. Self-adhering, back-to-back hook-and-loop closure allows you to quickly and easily secure arms to armboards, IV lines and catheters or hot/cold therapy devices to extremities, or even accessories to your O.R. table. Individually wrapped for greater infection control and ready to use right out of the box. Sold in boxes of 12 or cases of 120.



AliStrap	12/bx	10 boxes/cs	AliStrap Soft	12/bx	10 boxes/cs
1 ½"W x 30"L	926725	939163	1 ½"W x 30"L	926727	939165
1 ½"W x 38"L	939182	939183	1 ½"W x 38"L	939184	939183
3"W x 30"L	926726	939164	2"W x 30"L	926728	939166
			2"W x 38"L	939206	939207

VeriClean O.R. Table Strap

VeriClean O.R. Table Strap is the first 100% completely disinfectable strap for use in the operating room. Typical straps are sewn together, allowing for bacteria to permeate. This makes it extremely difficult to completely disinfect the strap without compromising its structural integrity. Because infection control is critical, this strap is entirely cloth- and thread-free.



Thermoplastic urethane belt • Welded seams • Cam-style metal buckles • Not made with natural rubber latex • 3-yr. warranty • Mid-panel measures 5"W x 18". Follow manufacturer's disinfection instructions.

- [937065](#) Medium, 96"L
- [937066](#) Medium/Large, 108"L
- [937067](#) Large, 120"L
- [937068](#) X-Large, 132"L

O.R. Table Strap, Deluxe Padded

Deluxe Padded Straps have a ½"-thick T-Foam™ padding (also called memory foam) between mid-panel and patient for maximum pressure reduction. T-Foam molds to the body and offers more uniform contact, helping to reduce pressure. Simply wipe clean this neoprene-coated vinyl strap with soap and water, then disinfect.



Straps: 2"W • Mid-panel: 5"W x 18"L or 6"W x 24"L • Neoprene-coated vinyl • Sewn seams • Cam-style metal buckles • Not made with natural rubber latex

- [933002](#) 96"L, 5" Mid-Panel
- [933004](#) 108"L, 6" Mid-Panel
- [933005](#) 120"L, 6" Mid-Panel
- [933006](#) 132"L, 6" Mid-Panel

O.R. Table Straps, Wipeable

Wipeable Straps feature a smooth surface that reduces opportunity for bacterial contamination. Simply wipe clean with soap and water, then disinfect. Cam-style metal buckles. Not made with natural rubber latex.



Strap: 2"W neoprene-coated vinyl • Mid-panel: 5"W x 18"L • Sewn seams

- [931088](#) 96"L
- [933008](#) 108"L
- [933009](#) 120"L
- [933010](#) 132"L

O.R. Table Straps, Basic

Basic O.R. Table Straps offer the same great quality and workmanship as the original Amsco® strap, but for less money! Heavy-duty, reinforced conductive vinyl mid-panel provides stability and redistributes pressure evenly across the abdomen. 2" W nylon strap and cam-style buckles with teeth for secure locking. Long enough to go around the patient and table. Sewn seams. Total length is 96".



Mid-panel is 5"W x 18"L • Not made from natural rubber latex

[95-082](#) Basic O.R. Table Strap

Disposable Surgical Strap

Disposable Surgical Straps offer a clean strap that's always ready to use. Strong black web strap with hook-and-loop closure for maximum adjustment. 3/8" thick hypoallergenic foam pad adheres to strap with hook-and-loop, offering protection between patient and strap. For use on tables with side rails. Not made with natural rubber latex.



Strap: 2"W x 53"L • Pad: 3"W x 25"L x 3/8" thick

[937545](#) Disposable Surgical Strap

AliGel™ Surgical Table Strap

AliGel Surgical Table Strap cushions and restrains the patient. Adjusts to fit all operating tables. One-piece, urethane Strap attaches to the surgical table side rails and secures back on itself with hook- and-loop closures. Thick gel mid-panel provides soft contact with patient. Not made with natural rubber latex.



4 1/2"W x 56"L • Gel mid-panel: 38"L x 3/8" thick

[931332](#) AliGel™ Surgical Table Strap

Nonmetallic O.R. Table Straps

Nonmetallic O.R. Table Straps eliminate protruding buckles. Heavy-duty conductive vinyl. Secure around patient with hook-and-loop -closures. Adjust to fit all surgical tables with side rails. Compression pads are ideal for reducing imaging artifacts on gas-distended patients. Not made with natural rubber latex.



Metal-free • Standard: 5"W x 70"L • Mid-panel: 18"L • Bariatric: 10"W x 120"L • Mid-panel: 24"L

[925470](#) Bariatric, 10"W [95-577](#) Standard, 5"W



Black Conductive Positioning Strap

The Black Conductive Positioning Strap One-piece black conductive vinyl strap that threads around side rails with strong hook-and-loop closures. Metal-free. 75"L. Available in two widths: 3"W or 5"W.



- [95-080](#) 3"W x 75"L
- [95-081](#) 5"W x 75"L

VeriClean™ Single Straps

VeriClean™ Single Straps are the first O.R. table straps with ultrasonically welded seams allowing for fully submersible, high-level disinfection. Free of punctures, holes, or stitching, allows for easier cleaning and fewer potential sources of contamination. Wipe down or send to Central Processing. Can be secured around O.R. table side rails or with strap hooks. Fits all surgical side rails.



2"W • Thermoplastic urethane belt • One cam-style metal buckle • Welded seams • Not made with natural rubber latex • 3-year warranty

- [937121](#) Medium, 2"W x 96"L
- [937122](#) Medium/Large, 2"W x 118"L
- [937123](#) X-Large, 2"W x 132"L

Economy Strap

The Economy Strap fits all surgical tables at a fraction of the cost of straps from table manufacturers. Features a 2"W neoprene-coated nylon strap and single cam-style buckle with teeth for secure locking. Choose from three lengths: 96"L, 118"L or 132"L. Not made from natural rubber latex.



- [933011](#) 96"L
- [933012](#) 118"L
- [933013](#) 132"L

Bariatric Super Strap

Bariatric Super Strap features four 55"L conductive, sewn neoprene-coated vinyl straps wrap around surgical table side rail. The 1" thick pad cushions the patient for proper pressure-management protection. This extra-strong, extra-wide strap provides a secure padded fit, designed to accommodate your largest patients. Not made with natural rubber latex. Sewn seams.



Mid panel: 12"W x 30"L x 1" thick • Usable strap: 80" diam. • Cam-style metal buckles • 1,000-lb. capacity

- [932447](#) Bariatric Super Strap

Bariatric VeriClean™ O.R. Table Strap

Bariatric VeriClean™ O.R. Table Strap offers the same infection-fighting features as standard VeriClean, but in a larger size to accommodate bariatric patients. VeriClean™ Bariatric is stitch-free, cloth-free, and made of durable urethane with welded seams. It's the only bariatric O.R. table strap that can be fully disinfected. The welded seams eliminate a potential hiding spot for bacteria. VeriClean™ can be completely submerged in approved disinfectants, making it ideal for cases that require more than a surface wipe-down. Fits all surgical tables with side rails or strap hooks. Features sewn seams, four strong cam-style buckles, and 2" straps. Not made with natural rubber latex.



Mid panel: 12"W x 30"L x 1" thick • Usable strap: 80" diam. • Cam-style metal buckles • 1,000-lb. capacity

[937678](#) Bariatric VeriClean™ O.R. Table Strap

Universal Straps

Universal Straps are unique extraordinarily strong straps that can accommodate just about any patient strapping job! The entire strap length is double-sided loop material, so the hook can securely close on any area on either side. Shorten the strap to accommodate the smaller patient by overlapping the one-piece circumferential strap. Loosen the overlap for the larger patient. To permanently shorten a strap, simply scissor-trim the blue loop material to the appropriate length. The versatile strap of choice!



One Piece Circumferential

91-277	2"W x 76"L
91-281	2"W x 96"L
932579	2"W x 108"L
91-278	3"W x 76"L
91-282	3"W x 96"L
932580	3"W x 108"L
91-279	4"W x 76"L
91-283	4"W x 96"L
932581	4"W x 108"L
91-280	6"W x 76"L
91-284	6"W x 96"L
932582	6"W x 108"L

One-Piece, Side Rail Hold

932570	2"W x 76"L
932571	2"W x 96"L
932572	2"W x 108"L
91-597	3"W x 65"L
932573	3"W x 96"L
932574	3"W x 108"L
91-598	4"W x 65"L
932575	4"W x 96"L
932576	4"W x 108"L
91-599	6"W x 65"L
932577	6"W x 96"L
932578	6"W x 108"L

Two-Piece, Side Rail Hold

91-285	2"W x 47"-76"L
932583	2"W x 76"-108"L
91-286	3"W x 47"-76"L
932584	3"W x 76"-108"
91-287	4"W x 47"-76"L
932585	4"W x 76"-108"L
91-288	6"W x 47"-76"L
932586	6"W x 76"-108"

Soft-Stretch Armboard Straps

Soft-Stretch Armboard Straps Single-use, smooth padded straps conform to any anatomy due to their moderate stretch core. Soft, hook-sensitive outer surface laminated to a gentle-padded foam interior. Infinite hook-and-loop attachment for secure fit on any board/arm.

3"W or 1 ½"W x 32"L • Hook-sensitive outer/foam-padded inner • Hook-and-loop closure • Not made with natural rubber latex • Beige

	3" Wide	1 ½" Wide
12/cs	938168	938166
48/cs	93816804	93816604
144/cs	93816812	93816612



Soft-Touch Armboard Straps

Soft-Touch Armboard Straps These ultra-smooth woven straps offer maximum comfort and fragile skin protection in an economical, single-use strap. Minimal stretch ensures a firm hold while providing some conformity.

3"W or 1½"W x 32"L • Woven fabric • Hook-and-loop closure • Not made with natural rubber latex • White

	3" Wide	1 ½" Wide
12/cs	938162	938164
48/cs	93816204	93816404
144/cs	93816212	93816412



Soft-Secure Armboard Straps

Soft-Secure Armboard Straps These soft, secure padded straps are resistant to stretch, offering a gentle, yet secure hold. Soft, hook-sensitive outer surface laminated to a gentle padded foam interior. Infinite hook-and-loop attachment for secure fit on any board/arm. White/Blue.

3"W or 1 ½"W x 32"L • Hook-sensitive outer/foam padded inner • Hook-and-loop closure • Not made with natural rubber latex

	3" Wide	1 ½" Wide
12/cs	938158	938160
48/cs	93815804	93816004
144/cs	93815812	93816012



Reusable Armboard Straps

Reusable Armboard Straps Hold patient's arms securely to armboard during surgery. Broad surface area reduces pressure points on arm – less chance of skin breakdown! Constructed of black conductive vinyl with hook-and-loop closure – comfortable against skin. One size fits most. Cleans easily with disinfectant. Not made with natural rubber latex. Sold in pairs.



4"W • Available in 30" and 36" lengths

[931902](#) 30"L [937546](#) 36"L

SUPPORT SURFACES

Gel Surgical Surface

Gel Anti-Shear Surgical Surfaces. Exposed layer of Akton® polymer gel (3/8" thick) is heat-sealed over a high-density Polyfoam pad (2" thick), providing constant support beneath patients without bottoming out. Anti-shear, elastic cover with fluid-proof, welded seams encases foam layer. Not made with natural rubber latex.



[937867](#) AMSCO C-Max Series [937865](#) Skytron 3600 Series
[937864](#) Skytron 3500 [937866](#) Skytron 6700 Series

Stretcher Pads

Basic Stretcher Pads are the most economical choice with 3" or 4" Polyfoam for patient comfort over a short period time. Include fluid-resistant seams for wipe-down ease and an anti-static conductive vinyl cover. Cover is stain- and tear-resistant and designed for long-lasting durability and easy maintenance.



Deluxe Stretcher Pads are the best pressure management option for when patients are on stretchers for longer periods of time. Deluxe pads have a T-Foam™ top layer with Polyfoam that lends itself to longer patient use, providing stability, support, and maximum pressure distribution.

Bacteriostatic and flame-retardant. (NFPA 16CFR part 1632) Class A. Meet Boston Fire Code. Not made with natural rubber latex. Include antimicrobial properties.

Model: Midmark/OMI/Pedigo 5822002

[910394PB](#) 3" Thick, Basic [910394DBW](#) 3" Thick, Deluxe
[910395PB](#) 4" Thick, Basic [910395DBW](#) 4" Thick, Deluxe

SURGICAL TABLE ACCESSORIES

Freedom® 360 Socket

FREEDOM® 360 Socket adjustability makes it easier and quicker to reposition a patient. Keep other settings constant—adjust location on the table rail, adjust the angle of the accessory, or adjust the height of the accessory. EZ-On connection attaches any-where along the rail, even between standoffs and accessories, allowing quicker setups and breakdowns. Store the clamp on the accessory and never have to look for the clamp. Quickly identify the desired position with clear markings. Orange safety line alerts staff when socket is not fully secured.



5 ¼"W x 2 ½"H x 3 ½"D • Anodized aluminum • 2 lbs. • Accepts round posts up to 11/16" diameter • Accepts flat blades up to 1"W x ¼" thick • Fits all O.R. tables with standard U.S. side rails ¾"W x 1 ½"H • Rotates 360° in 5° increments • 3-year limited warranty

[937932](#) FREEDOM 360 Socket

Freedom® Clark Socket

Freedom® Clark Socket This super-lightweight socket is up to 50% lighter than stainless steel sockets. Accepts blades AND posts. EZ-On connection attaches any-where along the rail, even between standoffs and accessories, allowing quicker setups and breakdowns. Ergonomic handle enables stronger torque and is easy to see amidst drapes and equipment. Rotational increments in 5° offer 72 possible positions for maximum placement options. Orange safety line alerts staff when socket is not fully secured.



3 ½"H x 2 ½" diam. • Anodized aluminum construction • 0.8 lb. • Accepts round posts up to 11/16" diameter • Accepts flat blades up to 1"W x ¼" thick • Fits all O.R. tables with standard U.S. side rails (¾"W x 1 ½"H) • Rotates 360° in 5° increments • 3-year limited warranty

[934725](#) FREEDOM Clark Socket

Clark Socket Plus

The **Clark Socket Plus** family has all the features of the Original Clark Socket, plus it holds FLAT blades too! By accommodating posts AND blades, you reduce socket inventory and save money—all while using a product you know and trust. Durable stainless steel ensures a long-lasting socket life. Our EZ-On model can be placed anywhere along the rail, allowing for quicker setups and breakdowns.



2 ¾"W x ¾"H x 1 ¾"D • Stainless steel • 2 lbs. • Accepts round posts up to 11/16" diameter • Accepts flat blades up to 1"W x ¼" thick • Fits all O.R. tables with standard U.S. side rails (¾" x 1 ½") • Rotates 360° in 9° increments • 3-year limited warranty • EZ-On Model: Connects anywhere along the rail • Slide-On Model: Slides on between standoffs

[937130](#) Clark Socket Plus, Slide-On
[937344](#) Clark Socket Plus, EZ-On
[937344EU](#) Clark Socket Plus, EZ-On (fits EU side rails)

Original Clark Socket

Original Clark Socket A tried-and-true design that facilities have relied upon for years—and still do today. Securely accommodates round posts. Solid, stainless steel construction ensures long life. Teeth won't chip or wear down—ever.

2 3/4"W x 7/8"H x 1 3/4"D • Stainless steel • 2 lbs. • Accepts round posts up to 11/16" diameter • Fits all O.R. tables with U.S. side rails (3/8"W x 1 1/8"H) • Rotates 360° in 12° increments • Slides on between standoffs • 3-year limited warranty

[936673](#) Clark Socket



Lightweight Sockets

Lightweight Sockets are incredibly light compared to stainless steel yet still offer a high level of durability. They accommodate a variety of posts—choose Flat, Round, or Combo. Ergonomic handle enables stronger torque, and high-visibility orange knob makes it easy to find in a sea of drapes and equipment. Accessories mount vertically only. 2/cs.

Flat Blade Socket: Accepts 1/2"W x 1"L blades • **Round Post Socket:** Accepts 5/8" diam. post • **Combo Socket:** Accepts both 1"W x 1/4" thick blade or 5/8" diam. post • All sockets: Anodized aluminum construction • 5.6 oz. • Fits O.R. tables with standard U.S. side rails 3/8"W x 1 1/8"H • 3-year limited warranty

[934706](#) Flat Post Socket, 2/cs

[934707](#) Round Post Socket, 2/cs

[934726](#) Combo Socket, 2/cs



Arthroscopic Leg Holders

Arthroscopic Leg Holders easily and effectively immobilize leg for Arthroscopic procedures. Vertical stress post adjusts 7" to 11"H. Lateral stabilizing blade accommodates a 4 1/2" to 8 1/2" leg diameter and pivots 360°. Saves setup time by being able to leave holder sized for an average leg. Single-Use Pads (sold separately) provide cushion and comfort, protecting skin and tissue from trauma. Round post, compatible with a Clark Socket (sold separately).

7"H x 2" OD Stress Post • 7"H x 11"L x 3 1/2" radius Stabilizing Blade • 10"H x 5/8" OD stainless steel mounting post • 4 lbs.

[936699](#) Arthroscopic Leg and Knee Holder

[930371](#) Disposable Pads



Toboggans & Polyfoam Liners

Toboggan Arm/Leg Guards are easy-to-use, easy-to-clean, crystal clear radiolucent shields that keep arms and legs on the table, help prevent interference with IV lines, and help protect the ulnar nerve against pressure, among other uses. Just slip the Toboggan under the pad and the patient’s weight holds it in place.



Convoluted Polyfoam Liners Single-patient-use liners offer protection without the risk of cross-contamination. Fits all Toboggan styles and sizes. Can be cut to size. 12 pr/cs.

Vivak® thermoplastic copolyester • Small: 10"L x 6"H x 10¼"D: M: 15¼"W x 8"H x 10¼"D
 • Large, 15¼"L x 10½"H x 10¼"D • Polyfoam Liners: 8"W x 20"L x 2" thick



- [931177](#) Toboggans, Small, Pair
- [916341](#) Toboggans, Medium, Single
- [91-634](#) Toboggans, Medium, Pair
- [931176](#) Toboggans, Large, Pair
- [930159](#) Convoluted Polyfoam Liners

SAFE PATIENT HANDLING

Patient Transfer Devices

Patient Roller Transfer Devices Easily transfer patient from one surface to another with very little effort and with as few as one caregiver. Patient can be transferred between slightly uneven surfaces (up to 3" height difference) without lifting. Long-lasting, trouble-free aluminum rollers with steel ball bearings and heavy-duty cover. Replacement black conductive vinyl covers fit most rollers on the market.



Short: 14"W x 25"L • 5 ½ lbs • Medium: 14"W x 30"L • 6 ¾ lbs • Long: 14"W x 67"L,
 • 13 lbs • Clean with soap and water or disinfectant • Aluminum rollers • Conductive vinyl
 • 440-lb. capacity • Black

	Roller	Roller/Rack Kit
Small	9-254	920889
Medium	9-729	920890
Large	9-728	920891

Soft Rollboards

Soft Rollboards with conveyor-belt mechanics allow for a smooth transfer of your patients with a gentle push—no lifting required to reduce strain on the neck and shoulders—while the supportive Soft Core provides comfort for delicate patients. All board frames are treated with microbe-inhibiting, silver-based Biomaster Antimicrobial Technology. For total infection control, Gold Rollboards also feature Ultra-Fresh™ ionic-silver antimicrobial protection on fabric covers to prevent the growth of odor- and stain-causing microbes. Each Gold model includes a three-year warranty.



Choice of Short or foldable Long model for easy storage or transfer in Fowler's position. Wide option also available in each length. All boards include handles for easy transport. Optional Disposable Covers, Cover Dispenser System, and Mounting Hangers sold separately.

Short Narrow: 15"W x 48"L • 3 lbs. • Short Wide: 19"W x 34"L • 6 lbs. • Long Narrow: 15"W x 68"L • Long Wide: 19"W x 68"L • 6.3 lbs. • Lightweight and radiolucent • MR-safe • Antistatic • Wipeable • Blue: 1-year warranty • Gold: 3-year warranty

	Short Narrow	Short Wide	Long Narrow	Long Wide
Blue	938656	938655	938657	938654
Gold	938652	938651	938653	938650

Patient Shifters

Anti-Stat Patient Shifters Help prevent nurse injuries with no-lift patient transfer devices. Shifters slide easily onto a bed or stretcher pad. Slippery coating protects the patient from sticking to the mover and shields against static buildup. Thin, warp-resistant, smooth polyethylene board has built-in handles. Semi-rigid board bridges uneven gaps. Antimicrobial version uses Zinc to reduce cross-contamination risk.

Conductive Shifters eliminate surprise “zaps” that can be painful and dangerous for patients and jeopardize electrical equipment. We guarantee that static charges won't build up on this shifter—no more painful shocks! Carbon-filled plastic actually conducts electricity--not a surface treatment.



250-lb. capacity • 72"L x 3/16" thick • 250-lb. capacity • **Narrow:** 18"W • 7 3/4 lbs. **Standard:** 22"W • 9 lbs. • **Rectangle:** 23 3/4 "W x 45 3/4 "L • 7 lbs.

9-719	Antistat Shifter, Standard 22"W, White
9-076	Antistat Shifter, Standard 22"W, Aqua
9-077	Antistat Shifter, Narrow 18"W, Aqua
9-182	Conductive Shifter, Standard 22"W, Black
938288	Conductive Shifter, Narrow 18"W, Black
938289	Conductive Shifter, Rectangle 18"W, Black

Bariatric Patient Shifters

Bariatric Shifters Best option for lateral transfers of patients up to 600 lbs. Reduce caregiver back strain by minimizing friction on transfer surfaces. Elongated hand holes for secure grip. Flexible plastic is great for uneven surfaces. We recommend using one caregiver per 100 lbs.

38"W x 72"L x 3/16"thick • 600-lb. capacity • Purple

[935944](#) Bariatric Shifter



Half Shifters

Half Shifters decrease transfer risk and help keep staff and patients safe. Transfer patients in any position—upright, lying, or those with epidural catheters or spinal anesthesia lines. Six built-in elongated hand holes provide a secure grip and also function as hanging holes. Bridges small gaps for uneven transfers. Warp-resistant for a smooth sliding surface. Store a Half Shifter in every room for staff safety and convenience. Not made with natural rubber latex.

23"W x 31"L x 0.08" thick • MR-safe • 1.7 lbs. • 250-lb. capacity • Radiolucent • Stain-resistant • Clean with hospital-approved disinfectant • White

[931827](#) Half Shifter



RADIATION PROTECTION

Disposable Thyroid Shields & Headbands

Ultralight Disposable Thyroid Shields or Headbands Durable enough to use multiple times, and economical enough to throw away. Full 0.5-mm Pb-equivalent lead-free radiation protection. Thyroid shields are recommended for staff who prefer a hygienic option for protection and are lightweight and comfortable. Head Bands have an elastic back to fit snugly and help protect vital brain tissue. 10/bx.

Spunbound polypropylene cover • Light Blue

[937929](#) Disposable Thyroid Shields

[937929](#) Disposable Headbands



Washable Thyroid Shields

Ultralight Washable Thyroid Shields can be washed in temperatures in excess of 140°F, up to 170 times. A clean, hygienic approach that minimizes risk of cross-contamination* Made of a proprietary, lead-free alloy for full 0.5 mm Pb-equivalent radiation protection. Extremely lightweight and comfortable.

Ripstop nylon cover • Machine-washable *Data on file.

[938186](#) Washable Thyroid Shield

Diamond pattern. Specify color: Burgundy, Charcoal Grey, Black, Royal Blue, Navy, Red

